


# Celebrity® DX Sport 4-wheel

Model number:	SC4401DXSP
Weight capacity:	350 lbs.
Maximum speed:	Up to 9.5 mph
Ground clearance:	3.25" (at motor)
Turning radius:	59.75"
Overall length:	47.5"
Overall width:	25.125"
Wheels:	
Front:	3.5" x 10" Pneumatic
Rear:	4" x 13" Pneumatic
Anti-tips:	Rear
Suspension:	Front & rear
Transport weight	
Total w/o batteries:	189 lbs.
Heaviest piece when disassembled:	63 lbs. (front section)
Drive train:	Sealed transaxle, 24 volt DC motor
Dual braking system:	Electronic, regenerative and electromechanical (electronic, regenerative disc brakes)
Std. Seating:	CRS Deluxe seat with headrest and limited recline
Battery charger:	Off-board, 5 amp (standard)
Per-charge range:	Up to 10.81 miles
Battery requirements:	(2) 12V deep cycle (AGM or gel cell type recommended)
Battery size:	40 AH: 32.5 lbs. <sup>1</sup> each
Warranty:	3-year limited warranty

Colors:



<sup>1</sup> Battery weight may vary based on manufacturer

## FEATURES

- State-of-the-art rear suspension
- Larger rear wheels
- Top speed of 9.5 mph
- Maximum weight capacity of 350 lbs.
- Ultra-simple, auto-latching lockup for quick disassembly/reassembly
- Easily accessible tie down points (for transport of unoccupied scooter)
- Tight turning radius of 59.75"
- Front basket and rearview mirror standard
- Includes standard angle-adjustable dual bulb curb light
- Front frame-mounted seat post for optimal stability
- High visibility reflectors
- Microprocessor-based controller
- Tiller rotation lock permits tiller to be used as a handle during transport

## OPTIONS

- Solid tires
- Forearm crutch holder
- Cane/crutch holder
- Wishbone crutch holder
- Oxygen tank holder
- Cup holder
- Rear basket
- Safety flag
- Walker holder
- Saddlebag
- Quad cane holder
- Jumbo front basket
- Weather cover

## Pride® Lifts & Ramps

- Backpacker® Plus
- Full Platform
- Silver Boom 300
- Backpacker® MV



These lifts are compatible with the Celebrity® DX Sport.

# CELEBRITY® DX SPORT

THE ULTIMATE BLEND  
OF STYLE AND PERFORMANCE

## CELEBRITY® DX SPORT

The **Celebrity® DX Sport** delivers easy operation, high performance, and sleek style. Loaded with standard features like an increased top speed, larger rear wheels and upgraded electronics, the Celebrity DX Sport is equipped to excite.

