

The Jazzy® 614 HD features patented Mid-Wheel 6® technology and Active-Trac® ATX Suspension for an outstanding level of performance, maneuverability and tight turning radius, all in one heavy-duty package. All-around performance indoors and out over most surfaces, makes the Jazzy 614 HD an excellent heavy-duty power chair choice.

#### **ACTIVE-TRAC® ATX SUSPENSION**





**6" CASTERS** 



**EXTENDED SEATING SIZES** 

# NORTH AMERICA'S FAVORITE POWER CHAIR

Learn more online.









gojazzy.com

## Specifications

| Weight Capacity         | 450 lbs.   |
|-------------------------|--|
| Maximum Speed           | Up to 4.25 mph   |
| <b>Ground Clearance</b> | 2.25"  |
| Turning Radius          | 22.75"   |
| Overall Length          | 35.5"  |
| Overall Width           | 25.25"   |
| Drive Wheels            | 14" solid  |
| Front Wheels            | 6" solid casters   |
| Rear Wheels             | 6" solid casters   |
| Suspension              | Active-Trac® (ATX)   |
| Drivetrain              | Two-motor, in-line,<br>Mid-Wheel 6®  |
| Braking System          | Regenerative and electro-mechanical  |
| Standard Electronics    | 50A, Dynamic Shark controller  |
| Battery Charger         | Off-board, 5 amp   |
| Per-Charge Range        | Up to 11 miles   |
| Battery Requirements    | (2) 12 volt, NF-22   |
| <b>Battery Weight</b>   | 37 lbs. each   |
| Base Weight             | 122 lbs.   |
| Standard Seat Weight    | 52 lbs.  |
| Standard Seat           | 20" x 18-20", 22" x 20"-22"  |
| Max. Seat Dimensions    | 24" x 24"  |
| Warranty                | Lifetime limited on frame;<br>13-month limited on batteries;<br>13-month limited on electronics;<br>13-month limited on drive motors |

### Colors



Candy Apple Red



Blue

#### Features and Benefits

- · Depth adjustable, limited recline high back seat available
- · Universal, fixed, in-line joystick mount
- · Flip-up, angle-adjustable foot platform
- Lap belt (50", 70", 80", 90")

## Seating Options\*

- · Deluxe contour, limited recline high back available sizes: up to 24"x24"
- Deluxe contour, limited recline, high back with solid seat pan, available sizes: 16"x16"-20"x22"
- · Synergy® Seat

## Options

- · Elevating leg rests
- · Heel loops
- Residual limb support
- · Swing- away footrests
- · Swing-away joystick

#### Accessories

- · Cane/crutch holder
- Cell phone holder
- · Cup holder
- Lap belt (50", 70", 80", 90")
- · Oxygen tank holder
- · Rear basket
- Saddle bag

- · Walker holder
- · Weather cover
- XLR USB charger

## Recommended Lifts

- Backpacker® Series
- · Zeus Series
- Outlander Full Platform
- Commander 450



<sup>\*</sup>Additional charges may apply.