

End-to-End Power Chair Solutions



	JAZZY AIR® 2	JAZZY® SELECT 6 2.0	JAZZY® PASSPORT	GO CHAIR®	JAZZY® ELITE ES PORTABLE	JAZZY® ELITE 14	JAZZY® ELITE HD	JAZZY® SELECT 6	JAZZY® 614 HD	
Weight Capacity	300 lbs.	300 lbs.	250 lbs.	300 lbs.	300 lbs.	300 lbs.	450 lbs. ¹¹	300 lbs.	450 lbs.	
Turning Radius¹	24.25"	23.76"	24.25"	25.4"	25.5"	24"	24"	21.875"	22.75"	
Width¹	23.19"	23.27"	23.5"	22.6"	22.75"	24.25"	24.25"	23.5"	25.25"	
Length¹	41.93"	43.27"	40"	28.5" ²	34.75" ²	35.5"	35.5"	35.375" ²	35.5" ²	
Maximum Speed³	Up to 4.0 mph (elevated and non-elevated position)	Up to 4.3 mph	Up to 3.6 mph	Up to 3.7 mph	Up to 4 mph	Up to 4 mph	Up to 4 mph	Up to 4 mph	Up to 4.25 mph	
Ground Clearance¹	1"	2.5" at footplate, 1.77" at motor bracket	3"	2"	1.5"	3"	3"	1.75"	2.25"	
Front Wheels	6" flat free	6" flat free	8" solid casters	6" solid casters	3" solid anti-tips	3" solid anti-tips	3" solid anti-tips	6" solid casters	6" solid casters	
Drive Wheels	10" flat free	10" flat free	12" Pneumatic drive	9" solid	9" solid	14"	14" solid	10" solid	14"	
Rear Wheels	6" flat free	6" flat free	3" solid anti-tips	3" solid anti-tips	6" solid casters					
Drivetrain	Two-motor, in-line, Mid-Wheel 6® drive	Two-motor, in-line, Mid-Wheel 6® drive	Two-motor, rear-wheel drive	Two-motor, in-line, rear-wheel drive	Two-motor, in-line, front-wheel drive	Two-motor, front-wheel drive	Two-motor, front-wheel drive	Two-motor, in-line, mid-wheel drive	Two-motor, in-line, Mid-Wheel 6®	
Braking System	"Intelligent Braking" electronic regenerative, disc park brake	"Intelligent Braking" electronic regenerative, disc park brake	Regenerative and electro-mechanical	Regenerative and electro-mechanical	Regenerative and electro-mechanical	Regenerative and electro-mechanical	Regenerative and electro-mechanical	Regenerative and electro-mechanical	Regenerative and electro-mechanical	
Suspension Type	Active-Trac® with Mid-Wheel 6® Technology	Active-Trac® with Mid-Wheel 6® Technology	Limited	Limited	Limited	Limited	Limited	Active-Trac® (ATX)	Active-Trac® (ATX)	
High-Back Seat (max. size)	20"W x 20"-22"D	20"W x 20"D	N/A	N/A	20"W x 20"D	20"W x 20"D	24"W x 24"D	20"W x 20"D	24"W x 24"D	
High-Back Seat-to-Floor Range	21.25" - 33"	21.25" - 23.25"	19"	N/A	21.5" - 23.5"	20" - 22"	20" - 22"	21.25" - 23.25"	20.5" - 23.5"	
Synergy® Seating	N/A	20"W x 20"D	N/A	N/A	N/A	20"W x 20"D	24"W x 24"D	20"W x 20"D	24"W x 24"D	
Synergy® Seat-to-Floor Range	N/A	16.5" - 18.5"	N/A	N/A	N/A	15.75" - 17.75"	15.75" - 17.75"	16.5" - 18.5"	17.75" - 19.75"	
Power Elevating Seat	Standard	Optional ⁴	N/A	N/A	N/A	N/A	N/A	Optional ⁴	N/A	
TRU-Balance®	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	
Specialty Controls	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A	
Standard Electronics	40A, PG VR2	40A, PG GC3 ⁴	50A, PG, VSI	40A, PG GC3	40A, PG GC3	50A, Dynamic Shark	60A, Dynamic Shark	40A, PG GC3 ⁴	50A, Dynamic Shark	
Battery Size (2 required)	12 volt, NF-22 Interceptor (included), 12 volt, U1 (optional), 40 Ah (optional)	12 volt, NF-22 Interceptor (included) or 12 volt, U1 (optional)	(1) 24V 18Ah lithium-ion (included) 25.2V 10.4Ah lithium-ion travel battery (optional)	18 Ah, 12 volt, deep-cycle (included)	12 volt, U1 (included)	12 volt, NF-22 Interceptor (included)	12 volt, NF-22 Interceptor (included)	12 volt, NF-22 Interceptor (included) or 12 volt, U1 (optional)	12 volt, NF-22 Interceptor (included)	
Standard Battery Charger	5A, off-board	3.5A, off-board	2A, off-board	2A, off-board	3.5A, off-board	5A, off-board	5A, off-board	3.5A, off-board	5A, off-board	
Per-Charge Range^{3,6}	Up to 16.8 miles (U1), Up to 19.2 miles (40 Ah)	Up to 19 miles ¹⁰	Up to 9.6 miles at 250 lbs. Up to 11 miles at 200 lbs. Up to 5.7 miles (10.4 AH travel battery at 250 lbs.) Up to 6.3 miles (10.4 AH travel battery at 200 lbs.)	Up to 13.2 miles at 200 lbs. Up to 8.7 miles at 300 lbs.	Up to 11 miles	Up to 9.8 miles	Up to 9.8 miles	Up to 15 miles	Up to 11 miles	
Battery Weight⁸	24.5 lbs. each (U1), 32.5 lbs. each (40Ah)	24.5 lbs. each	7 lbs.	14 lbs. each	25 lbs. each	37 lbs. each	37 lbs. each	25 lbs. each	37 lbs. each	
Weight of Heaviest Piece⁵	180 lbs. (with seat) ⁹	114 lbs.	54 lbs. (unit doesn't disassemble) ⁷	36 lbs.	48 lbs.	99 lbs.	99 lbs.	100 lbs.	122 lbs.	
Standard Seat Weight	N/A	50 lbs. (comfort, high-back)	N/A	30 lbs. (compact, foldable)	43 lbs. (comfort, high-back)	52 lbs. (high-back)	52 lbs. (high-back)	43 lbs. (comfort, high-back)	52 lbs. (high-back)	
Standard Seat Size	16" x 16" - 18", 18" x 18" - 20", 20" x 20" - 22"	16" x 16" - 18", 18" x 18" - 20", 20" x 18" - 20"	17.5" x 16.5"	18" x 17"	16" x 16" - 18", 18" x 18" - 20", 20" x 18" - 20"	18" x 18" - 20", 20" x 18" - 20", 22" x 20" - 22"	20" x 18" - 20", 22" x 20" - 22"	16" x 16" - 18", 18" x 18" - 20", 20" x 18" - 20"	20" x 18" - 20", 22" x 20" - 22"	
Warranty	5-year limited on frame; 24-month limited on electronics; 24-month on drive motors; 13-month on batteries	5-year limited on frame; 24-month limited on electronics; 24-month on drive motors; 13-month on batteries	5-year limited on frame; 13-month limited on electronics; 13-month on drive motors; 6-month on batteries	5-year limited on frame; 24-month limited on electronics; 24-month on drive motors; 13-month on batteries	5-year limited on frame; 24-month limited on electronics; 24-month on drive motors; 13-month on batteries	5-year limited on frame; 24-month limited on electronics; 24-month on drive motors; 13-month on batteries	5-year limited on frame; 24-month limited on electronics; 24-month on drive motors; 13-month on batteries	5-year limited on frame; 24-month limited on electronics; 24-month on drive motors; 13-month on batteries	5-year limited on frame; 24-month limited on electronics; 24-month on drive motors; 13-month on batteries	Lifetime limited on frame; 24-month limited on electronics; 24-month on drive motors; 13-month on batteries

1) Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.
 2) Without front riggings.
 3) Varies with user weight, terrain, battery type, battery condition and tire pressure.
 4) 40A, PG 6-key VR2 required with Power Elevating Seat Option.

5) Does not include seat or battery weights.
 6) Tested in accordance with ANSI/RESNA, WC Vol. 2, Section 4 and ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Testing conducted at maximum weight capacities. This specification can be subject to a variance of +10%, -5%.
 7) Seat weight included. Unit does not disassemble.

8) Battery weight may vary based on manufacturer.
 9) Seat weight included (because tools are required to remove seat). Battery weight not included.
 10) Tested with Interceptor U1 batteries.
 11) Maximum weight capacity is reduced to 400 lbs. (181 kg) on 18" deep seat configurations.

The information contained herein is correct at time of publication; we reserve the right to alter specifications without prior notice.

Beamsville, Ontario
 888-570-1113
 pridemobility.ca

