

LC-525 SERIES




FEATURES AND BENEFITS:

- Infinite-positions allows backrest and footrest to move independently from one another so that a variety of positions can be reached
- Standard large dual pockets
- Easy-to-use hand control
- Removable back
- Energy-efficient transformer, features self-diagnostic electronics
- Quiet and smooth lift system
- Wooden frame composed of engineered furniture grade laminate/hardwood for added strength and durability
- Heavy-duty lift actuator & scissor mechanism
- Integrated battery backup in case of power failure (batteries not included)
- Optional head pillow
- Optional lumbar pillow
- Standard with arm and headrest covers

AVAILABLE FABRICS

Cloud 9



BLACK CHERRY PACIFIC SPRUCE STONE WALNUT

Lexis Sta-Kleen®



CHESTNUT

(To be discontinued. Check for available stock. Not all fabric options may be available.)

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

MODEL	LC-525S	LC-525M	LC-525L
Size	Small	Medium	Large
Position	Infinite-Position, Zero Gravity, Trendelenburg Chaise Lounger 5'3" & Below	Infinite-Position, Zero Gravity, Trendelenburg Chaise Lounger 5'4" - 6'	Infinite-Position, Zero Gravity, Trendelenburg Chaise Lounger 5'9" - 6'2"
Weight Capacity	375 lbs.	375 lbs.	375 lbs.
Seat-to-Floor	19"	20"	21"
Seat Depth	18"	20"	22"
Seat Width	18"	20"	22"
Top of Back to Seat	25.5"	26.5"	27.5"
Back Style	Blown Fiber Split Back	Blown Fiber Split Back	Blown Fiber Split Back
Fully Padded Chaise	YES	YES	YES
Heat & Massage³	NO	NO	NO
Additional Options	NA	Footrest Extension (Std.)	Footrest Extension (Std.)
Seat Type	Sinuous spring seat	Sinuous spring seat	Sinuous spring seat
Distance From Wall	28"	29"	30"
Overall Width	32"	33"	35"
Chair Weight	123 lbs.	126 lbs.	128 lbs.
Hand Control	Easy Touch Infinite-Position with USB Port	Easy Touch Infinite-Position with USB Port	Easy Touch Infinite-Position with USB Port
Warranty	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated

Beamsville, Ontario • 888-570-1113 • pridemobility.com

