

Litestream[®] series



ACCESSORIES

REAR WHEELS




20", 22", 24", 25" and 26" high-performance Spinerogy[®] wheels feature double-wall rims with straight-pull spokes for reduced weight and optimal durability




22", 24", 25" & 26" Spinerogy[®] Light Extreme 12 spoke PBO Fiber wheels

AVAILABLE COLORS:



24", 25" and 26" Spinerogy[®] Sport Light Extreme 24 spoke PBO Fiber wheels

AVAILABLE COLORS:



20", 22", and 24" Composite mag wheels



24" lightweight, quick-release, one-arm drive system with aluminum handrim

CASTERS



4", 5" and 6" soft-roll casters smooth lumps and bumps in everyday use



5" aluminum casters with precision bearings and polyurethane tires



4" light-up casters



4" colored casters available in red, blue, green, yellow, orange or clear



4", 5", 6" & 8" black, smooth 3-spoke tires



6" pneumatic casters



8" pneumatic or flat-free polyurethane casters

HANDRIMS



20", 22", 24", 25" and 26" Projection Handrim (vertical or 45° oblique)



24" Spinerogy[®] FlexRims



22", 24", 25" and 26" Natural-Fit and 24" Natural-Fit LT Handrim (available in standard or super grip)

TIRES



24" & 26" Aggressive Knobby Tire



20", 22", 24" & 26" Everyday Gray Tire



20", 22", 24" & 25" Polyurethane Tire (24" also available in light gray)



22", 24", 25" & 26" 110 PSI Primo Tire



24" and 26" 110 PSI Primo Colored Tire



26" Kenda Mountain Bike Knobby Tires



24", 25" and 26" Schwalbe Marathon Plus Evolution



24", 25" and 26" Schwalbe Marathon Plus

(Aluminum, vinyl-coated & chrome handrims also available but not shown)

FORKS

WHEEL LOCKS



Frog Legs suspension fork



8" fork
Available as quick-release or high angle bolt-on



5" fork



Quick-release 60mm fork
(for use with 8" pneumatic or flat-free polyurethane casters only)



Pull-to-lock



Scissor lock



LTD/Traverse Wheel-lock extension



Push-to-lock



Pull-to-lock with extension



Push-to-lock with extension



Hill holder



Cable-linked, unilateral wheel locks (for use with one-arm drive system)

LEG RESTS



Elevating bi-directional leg rests with swing-away calf pads



60°, 70°, 70° tapered, 80° and 90° aluminum bi-directional swing-away leg rests



60° & 70° Hemi aluminum bi-directional swing-away leg rests



Flip-up aluminum footplates with heel loops



10-degree comfort adjust footplates in 4.25", 5" and 6" widths



Angle-adjustable footplates in 5.5", 6.5" and 7.5" widths



Nylon heel loops

BACK OPTIONS

TRANSIT SECUREMENT



Tension-adjustable upholstery



Standard nylon sling



Synergy® SHAPE back



The Transit Securement Package allows convenient securement of the Litestream® XF for safe transport. For additional information, please refer to the Transit Securement Package brochure.

ACCESSORIES AND PERSONALIZATION

ARMRESTS



Height-adjustable anti-tips



LTD Height-adjustable anti-tips



Composite side guards



Aluminum side guards



Combat, cane & crutch holder



Quantum® backpack

(Lap belt is also available but not shown)



12" tubular, swing-away, desk-length armrests (shown with optional gel arm covers)



Height-adjustable, flip-back, dual post armrests (desk or full length)



LTD Dual-post flip-back height adjustable arms



Height-adjustable drop-in armrests (desk or full length)

Beamsville, Ontario, **888-570-1113**
quantumrehab.com



QUANTUM REHAB®