

# Litestream<sup>®</sup> XF

The Litestream<sup>®</sup> XF Series of manual wheelchairs offer a unique, innovative, rigidizing double X folding frame that delivers the versatility and performance of a rigid frame while folding to a width of only 9" for easy transport. The Litestream XF Series accepts a broad selection of accessories to meet a wide variety of needs and personal preferences while remaining an excellent overall value. The Litestream XF is also available in an HD configuration.



Litestream<sup>®</sup> XF  
with Ultra-hemi configuration



Litestream<sup>®</sup> XF  
with HD configuration



Litestream<sup>®</sup> XF  
standard configuration



Litestream<sup>®</sup> XF LTD

# Litestream<sup>®</sup> series




# SPECIFICATIONS

	WEIGHT CAPACITY	SEAT-TO-FLOOR HEIGHTS	SEAT SIZES
LITESTREAM® XF	250 lbs.	15"-19.5" (front) 14.5"-19" (rear)	Width: 15"-21" Depth: 14"-20"
LITESTREAM® XF HD	350 lbs.	16"-19.5" (front) 14.5"-19" (rear)	Width: 18"-22" Depth: 18"-20"
LITESTREAM® XF ULTRA HEMI	250 lbs.	12.5"-16.5" (front) 12.5"-16.5" (rear)	Width: 15"-21" Depth: 14"-20"
LITESTREAM® XF LTD	250 lbs.	12.5"-18.5" (front) 12"-17" (rear)	Width: 15"-21" Depth: 14"-20"


## COLORS

### LITESTREAM® XF, HD AND ULTRA HEMI


#### METALLIC



#### NON-METALLIC



#### WATER TRANSFER PRINTING



### LITESTREAM® XF LTD

#### STANDARD



Black Abyss

#### OPTIONAL



Red Hot



Blue Sky

## OPTIONS/ACCESSORIES\* LITESTREAM® XF LTD

<b>REAR WHEELS:</b>	20", 22", 24" & 26" Composite mag 20", 22", 24" & 26" Spinergy® wire spoke
<b>TIRES:</b>	20", 22" & 24" Polyurethane tire 20", 22" & 24" Everyday gray tire
<b>HANDRIMS:</b>	20", 22", 24" & 26" Aluminum 20", 22", 24" & 26" Vinyl-coated 20", 22", 24" & 26" Projection 20", 22" & 24" Chrome 22", 24" & 26" Natural-Fit (Standard and Supergrip)
<b>AXLES:</b>	Quick-release, quad-release
<b>AXLE PLATES:</b>	Adjustable axle plate (standard) Amputee mounting (optional)
<b>CASTERS:</b>	4", 5", 6" & 8" Black polyurethane
<b>FORKS:</b>	5" or 8" bolt-on or quick-release forks (angle-adjustable) Frog Legs
<b>WHEEL LOCKS:</b>	Pull-to-lock Pull-to-lock with extensions Push-to-lock Push-to-lock with extensions
<b>ARMRESTS:</b>	Height-adjustable drop-in Dual-post, flip-back, height-adjustable arms Desk or full-length arm pads
<b>BACK CANES:</b>	8 degree with push handles Angle-adjustable, with push handles
<b>BACKS:</b>	Standard padded nylon sling with convenience pocket Tension-adjustable upholstery ADI backs
<b>FRONT RIGGINGS:</b>	Aluminum, bi-directional swing-away leg rests, 60° & 70° Hemi Flip-up aluminum footplates (with heel loops) Composite comfort adjust footplates Angle-adjustable footplates Nylon heel loops High pivot elevating leg rests with standard or offset calf pad
<b>ACCESSORIES:</b>	Height-adjustable, rear anti-tip wheels Lap belt Tool kit Combat cane/crutch holder Transit securement package

\*Not all options are available for all configurations.

# OPTIONS/ACCESSORIES\* LITESTREAM® XF, HD AND ULTRA HEMI

<b>REAR WHEELS:</b>	20", 22", 24", 25" & 26" Spinergy® wire spoke 22", 24" 25" & 26" Spinergy Light Extreme (available in red, yellow, blue, white, pink, green, orange and black) 24", 25" & 26" Spinergy Sport Light Extreme (available in red, yellow, blue, white, pink, green, orange and black) 20", 22" & 24" Composite mag
<b>TIRES:</b>	22", 24", 25" & 26" 110 PSI Primo tire 24" & 26" Aggressive knobby tire 20", 22", 24" & 26" Everyday gray tire 20", 22", 24" & 25" Polyurethane tire (24" also available in light gray) 26" Kenda mountain bike knobby tires (available in black or gray) 24" & 26" 110 PSI Primo colored tires (available in red, blue or yellow) 24", 25" and 26" Schwalbe Marathon Plus and Marathon Plus Evolution
<b>HANDRIMS:</b>	20", 22", 24", 25" & 26" Aluminum 22", 24", 25" & 26" Vinyl-coated 20", 22", 24", 25" and 26" Projection (vertical or 45° oblique) 22", 24", 25" & 26" Natural-Fit (available in standard or super grip) 24" Natural-Fit LT (available in standard or super grip) 20", 22", 24" Chrome
<b>AXLE PLATES:</b>	Standard (0°, 1.5°, 3° camber options) Amputee (0°, 1.5°, 3° camber options)
<b>CASTERS:</b>	4", 5", 6" & 8" composite 4", 5" & 6" soft roll 4" light-up casters 4" colored casters (available in red, blue, green, orange, yellow or clear) 8" pneumatic or flat-free polyurethane casters
<b>FORKS:</b>	5" or 8" quick release forks 5" or 8" bolt-on forks Frog Legs suspension forks 60mm forks (for use with 8" pneumatic casters only)
<b>WHEEL LOCKS:</b>	Pull-to-lock Pull-to-lock with extensions Push-to-lock Push-to-lock with extensions Scissor locks Hill holder
<b>ARMRESTS:</b>	12" tubular swing-away (desk length) Height-adjustable drop-in (desk or full length) Height-adjustable, flip-back, dual post armrests (desk or full length)
<b>BACK CANES:</b>	Straight, with or without push handles Angle-adjustable, with or without push handles 8 degree with push handles Straight reinforced
<b>BACKS:</b>	Standard nylon sling with convenience pocket Tension-adjustable upholstery ADI backs
<b>FRONT RIGGINGS:</b>	Aluminum, bi-directional swing-away leg rests, 60°, 70° and 90° Elevating leg rests with calf pads Comfort adjust footplates Angle-adjustable footplates Flip-up aluminum footrests Nylon heel loops
<b>ACCESSORIES:</b>	Height-adjustable, removable anti-tip wheels Composite side guards Lap belt Transit securement package Combat, cane and crutch holder

The information contained herein is correct at the time of publication. We reserve the right to change specifications without prior notice.