

Pride® Lifts

Mobility Made Easy.

END-TO-END



	Backpacker® AVP 2.0 Interior Lift	Backpacker® Plus Interior Lift	Commander 200 Interior Lift Hitch & Trunk	Outlander Full Platform Exterior Lift	Outlander DE Exterior Lift	Outlander XL Exterior Lift	Outlander TM Exterior Lift	Commander 450 Exterior Lift 2 & 3 Axis
BASE WEIGHT ⁵	193 lbs.	155 lbs.	Trunk Mount – 43 lbs. Hitch Mount – 49 lbs.	82 lbs.	84 lbs.	94 lbs.	79 lbs.	109 lbs.
MAX. WEIGHT CAPACITY	350 lbs.	350 lbs.	200 lbs.	350 lbs.	350 lbs.	350 lbs.	150 lbs.	450 lbs.
MAX. PRODUCT WIDTH	27"	27"	N/A	28"	27"	27"	24"	N/A
OVERALL PLATFORM LENGTH	42"-34" (41"-33" usable surface)	44" (43" usable surface)	N/A	43" (42" usable surface)	51" (42" usable surface)	56.75" (49" usable surface)	39" (37" max. product length)	N/A
MAX. TRAVEL DISTANCE FROM THE GROUND	37"	26"	70"	21"	21"	21"	21"	70"
MOTOR	Sealed 12 Volt DC	Sealed 12 Volt DC	12 Volt DC	Sealed 12 Volt DC	Sealed 12 Volt DC	Sealed 12 Volt DC	Sealed 12 Volt DC	12 Volt DC
WARRANTY	3-Year Transferable	3-Year Transferable	3-Year Transferable	3-Year Transferable	3-Year Transferable	3-Year Transferable	3-Year Transferable	3-Year Transferable
HITCH CLASS ^{4,6}	N/A	N/A	Class II, III & IV	Class II, III & IV	Class II, III & IV	Class II, III & IV	Class II, III & IV	N/A
ACCESSORY OPTIONS	<ul style="list-style-type: none"> Battery pack w/ charger Seat belt-style restraint used 	<ul style="list-style-type: none"> Battery pack w/ charger Seat belt-style restraint used 	Go Chair® docking kit	<ul style="list-style-type: none"> Battery pack w/ charger Swing-away adaptor⁶ Scooter hold-down arm Power chair hold-down bar Strap kit 	<ul style="list-style-type: none"> Battery pack w/ charger Swing-away adaptor⁶ Strap kit Scooter hold-down arm 	<ul style="list-style-type: none"> Battery pack w/ charger Swing-away adaptor⁶ Strap kit 	<ul style="list-style-type: none"> Battery pack w/ charger Swing-away adaptor⁶ 	N/A
VEHICLE PLATFORM	Full-size Van, Minivan, SUV & Crossover	Minivan, Full-size van ³	Sedan, Full-size Van, Minivan, Crossover & Hatch Back	Full-size Van, Minivan, Truck, Sedan, SUV and Crossover	Full-size van, Minivan, Truck, Sedan, SUV and Crossover	Full-size Van, Minivan, Truck, Sedan, SUV and Crossover	Full-size Van, Minivan, Truck, Sedan, SUV and Crossover	Full-size Van, Truck and SUV
COMPATIBLE PRODUCTS ^{1,2}	All Pride® Scooters & Pride®/Quantum® PWCs with contoured, Specialty & Synergy® seating And equivalent competitor products Except: Pursuit®, Pursuit® XL, Raptor	All Pride® Scooters & Pride®/Quantum® PWCs with contoured, Specialty & Synergy® seating And equivalent competitor products Except: Pursuit®, Pursuit® XL, Raptor	Go-Go® Travel Mobility Scooters, Jazzy® Elite ES Portable, Jazzy® Passport & Jazzy® Select Series Compatible with manual wheelchairs, most foldable power chairs and scooters.	All Pride® Scooters & Pride®/Quantum® PWCs with contoured, Specialty & Synergy® seating And equivalent competitor products Except: Pursuit®, Pursuit® XL, Maxima 3W & 4W, Raptor, Victory® and Celebrity® Scooter Lines	Victory® Twin, Victory® 10 3W & 4W, Victory® 10 LX w/CTS & Maxima 3W & 4W And equivalent competitor products Except: Pursuit®, Pursuit® XL, Raptor	Victory® 10 LX w/CTS, Pursuit® & Maxima 4W And equivalent competitor products Except: Pursuit® XL, Raptor	All Go-Go® Travel Mobility Products And equivalent competitor products Except: Go-Go® Folding Scooter	All Pride® Scooters, Power Chairs and some Quantum® units And equivalent competitor products Except: Pursuit® XL, Raptor
	INTERIOR LIFTS			EXTERIOR LIFTS				

1) Lifting weight does not include accessories or certain upgraded seating. Please contact Customer Service for greater detail.
2) Vehicle dependent.
3) May require vehicle modification.
4) Hitch weight carrying capacity varies by vehicle. Call for details.

5) May not include some accessories.
6) Class II, III and IV only.

To find a compatible lift for your vehicle and mobility product, use our new



pridemobility.com/lift-finder

Beamsville, Ontario • 888-570-1113
prideliftsandramps.com • pridemobility.com



The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.