



SPECIFICATIONS	Pursuit®	Pursuit® XL
FRONT TIRES	4" x 13" PNEUMATIC	13" PNEUMATIC
REAR TIRES	4" x 13" PNEUMATIC	13" PNEUMATIC
MAXIMUM SPEED²	UP TO 6 MPH HS: UP TO 9.5 MPH	UP TO 7.5 MPH
GROUND CLEARANCE^{1,3,4}	3.5"	5.25" at mid-frame
TURNING RADIUS⁴	68.5"	82.5"
LENGTH/WIDTH^{4,5}	53.5"/25"	56"/27.25"
SEAT-TO-GROUND HEIGHT RANGE	24.75"-26.25"	26.25"-27.75"
SEAT TYPE (STANDARD)	CRS DELUXE WITH LIMITED RECLINE (18" X 18")	HIGH-BACK RECLINER WITH SEAT SLIDE, SWIVEL, ADJUSTABLE ARMS AND HEIGHT-ADJUSTABLE HEADREST (18" X 18")
(3/4-WHEEL) WEIGHT^{3,4,6}	N/A/267 LBS.	259.9 LBS.
SEAT-TO-DECK HEIGHT RANGE	17.75"-18.75"	21"-22.5"
BATTERY WEIGHT⁸	53 LBS. EA.	70 LBS. EA. (G31) 60 LBS. EA. (G27) 53 LBS. EA. (G24)
BATTERY REQUIREMENTS⁷	(2) GROUP 24	(2) GROUP 31 ¹⁰ (110AH) (2) GROUP 27 (100AH) OR (2) GROUP 24 (75 AH)
PER CHARGE RANGE (MAX.)^{2,9}	UP TO 33.6 MILES HS: UP TO 22 MILES	G31: UP TO 31.83 MILES G27: UP TO 29.51 MILES G24: UP TO 22.13 MILES
BATTERY CHARGER	OFF-BOARD, 8A	OFF-BOARD, 8A
SUSPENSION	FRONT & REAR	FRONT & REAR
MAXIMUM WEIGHT CAPACITY	350 LBS.	400 LBS.
STANDARD BODY COLORS	CANDY APPLE RED, VIPER BLUE	CANDY APPLE RED, SILVER
WEIGHT OF HEAVIEST PIECE⁴	N/A UNIT DOES NOT DISASSEMBLE	205.5 LBS. (MAIN FRAME)

Footnotes

1. At midframe.
2. Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. This specification can be subject to a variance of (+ or -) 10%.
3. Does not account for anti-tip clearance.
4. Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3%.
5. Width of base (not including seat).
6. Includes standard seat weight.
7. Deep cycle (AGM or gel cell type recommended).
8. Battery weight may vary based on manufacturer.
9. Tested in accordance with ANSI/RESNA, WC Vol 2, Section 4 and ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
10. Option available through most providers. Not available through Pride Mobility®.

Dealer Information



PURSUIT®



PURSUIT® XL



**MODEL # S713ES CA (6 MPH)
S713ES HS (9.5 MPH)**



MODEL # SC714



FEATURES OF BOTH

- Black vinyl reclining high-back seat with sliders, swivel, adjustable armrests (height and width) and height-adjustable headrest
- Delta tiller with wrap-around handles, for left or right-handed operation, lets you rest your wrist comfortably while driving
- Easily accessible tie-down points (for transport of unoccupied scooter)
- Front and rear suspension
- Full lighting package includes directional signals and large lower headlight for pathway illumination
- Infinitely-adjustable tiller with conveniently located handle
- Large 13" pneumatic tires
- Powerful drivetrain features microprocessor-based controller; 140 amp (Pursuit), 200 amp (Pursuit XL)
- 24-volt, 4-pole DC motor for increased power
- Easy-to-read LED battery meter
- Easy-access tiller-mounted charger port
- Jumbo front basket (Pursuit only)

COLOR CHOICES

Pursuit



Candy Apple Red



Viper Blue

Pursuit XL



Candy Apple Red



Silver

ACCESSORIES

(See accessory trifold for more options)



*RAM® X-Grip® Cell Phone Holder **ACC145549**



Cup Holder **ACC120041**



50" Lap Belt **ACCBELT1002**



Rear Basket **ACCBSKT1010**



Rearview Mirror



Saddlebag
Small: **DWR1020D010**
Large: **DWR1020D016**



Safety Flag **ACCASMB2629**



Walker Holder **ACCASMB1018**



Weather Cover
(Large) **ACCCOVR1018**



XLR USB Charger
1pk: **ACC1704791**
12pk: **KIT1705402**