

# Engineered for Powerful Performance and Exceptional Style



The Pursuit® XL is an exciting blend of power and precision that features a 24-volt DC motor, full suspension, 13" low-profile tires and a hydraulic-sealed brake system. Standard luxury touches, such as an infinite adjustable wraparound delta tiller, high-back pillow top seat and LED lighting make it the perfect outdoor scooter.



Full Suspension



PRIDE® MOBILITY SCOOTERS  
LIVE YOUR BEST®

| [pridemobility.com](http://pridemobility.com) |





## Specifications 4-wheel Model # SC714

Weight Capacity	400 lbs.
Maximum Speed <sup>1</sup>	Up to 7.5 mph
Ground Clearance <sup>2</sup>	5.25" at mid-frame, 2.75" at motor
Turning Radius	82.5"
Overall Length <sup>2</sup>	56"
Overall Width <sup>2</sup>	27.25"
Seat-to-Deck Height	21"-22.5"
Seat-to-Ground Height	26.25"-27.75"
Front Tires	13" pneumatic
Rear Tires	13" pneumatic
Suspension	Front and rear
Range Per Charge <sup>1</sup>	Group 24 (70 AH): Up to 22 miles Group 27 (100 AH): Up to 30 miles Group 31 (110 AH): Up to 32 miles
Total Weight Without Batteries	254 lbs.
Heaviest Piece When Disassembled	N/A
Standard Seating	Type: Reclining high-back Weight: 45.5 lbs. Material: Black Vinyl Dimensions: 18" x 18"
Drivetrain	Rear-wheel drive, sealed transaxle, 24-volt DC motor
Brakes	Regenerative and electromechanical
Battery Requirements <sup>3</sup>	Size: Group 24 (75 AH) optional Group 27 (100 AH) optional Group 31 (110 AH) optional <sup>5</sup> Weight <sup>4</sup> : Group 24: 53 lbs. each Group 27: 60 lbs. each Group 31: 70 lbs. each
Battery Charger	8-amp, off-board
Colors	  Red Silver

## Features

- Delta tiller with wrap-around handles, for left or right-handed operation, lets you rest your wrist comfortably while driving
- Black vinyl reclining high-back seat with sliders, swivel, adjustable armrest (height and width) and height-adjustable headrest
- Easily accessible tie down points (for transport of unoccupied scooter)
- Front and rear suspension
- Front basket
- Full lighting package includes directional signals and lower headlight
- Hydraulic handbrake
- Infinitely-adjustable tiller angle with conveniently located handle
- Large 13" low-profile, pneumatic tires for excellent outdoor performance
- Powerful drivetrain; features 200 amp microprocessor-based controller
- Rearview mirror

## Accessories

- Cell phone holder
- Cup holder
- 50" Lap belt
- Rear basket
- Rear view mirror
- Saddlebag
- Safety flag
- Walker holder
- Weather cover
- XLR USB charger

<sup>1</sup> Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).

<sup>2</sup> Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).

<sup>3</sup> AGM or gel cell type required.

<sup>4</sup> Battery weight may vary based on manufacturer.

<sup>5</sup> Option available through most providers. Not available through Pride Mobility®.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for scooter only. Overall measurements will vary based on seating and accessory selections.