



**Q6 EDGE® 2.0**  
shown in *Lime Light* with  
12" iLevel® technology and  
Q-Logic 3 Advanced Drive  
Controls with Bluetooth®

# Q6 edge® 2.0

With features like standard 6 mph motors and advanced responsiveness, the Q6 Edge® 2.0 is engineered to meet the performance needs of the most active user. The Q6 Edge 2.0 accepts a complete range of seating and electronics options to deliver exceptional rehab capability at an excellent value, and is available with our iLevel® technology.

**QUANTUM®**

# Specifications

<b>DRIVE WHEELS</b>	Standard: 14" Solid Optional: 14" Pneumatic 14" Knobby
<b>CASTER WHEELS</b>	Front: 6" Solid (Optional: 8" Solid) Rear: 6" Solid (Optional: 8" Solid)
<b>ANTI-TIP WHEELS</b>	N/A
<b>MAXIMUM SPEED</b> <sup>1</sup>	Up to 6 mph
<b>GROUND CLEARANCE</b>	2.7" (battery tray)
<b>TURNING RADIUS</b> <sup>2</sup>	20.5" with 6" casters 23.25" with 8" casters
<b>OVERALL LENGTH</b> <sup>2</sup>	32.23" with 6" casters 40.48" with 8" casters
<b>BASE WIDTH</b>	24.4"
<b>SEATING SIZES</b>	
<b>Synergy®/TRU-Balance® 3 Static Seat</b>	<b>Width:</b> 10-22" <b>Depth:</b> 10-22"
<b>TRU-Balance® 3 Power Tilt</b>	<b>Width:</b> 12-22" <b>Depth:</b> 12-22"
<b>TRU-Balance® 3 Tilt &amp; iLevel®</b>	<b>Width:</b> 12-22" <b>Depth:</b> 12-22"
<b>SEAT-TO-FLOOR HEIGHTS</b>	
<b>Synergy®/TRU-Balance® 3 Static Seat</b> <sup>5</sup>	16.5-18.5"
<b>TRU-Balance® 3 Power Tilt</b> <sup>5</sup>	17.5-19.5"
<b>TRU-Balance® 3 Tilt &amp; iLevel®</b> <sup>5</sup>	17.5-19.5"
<b>TRU-Balance® Flex Seat</b> <sup>5</sup>	16.5"-19.5"
<b>MANUAL RECLINE</b> <sup>6</sup>	Yes
<b>BATTERY SIZE</b> <sup>3</sup>	NF-22
<b>BATTERY WEIGHT</b> <sup>7</sup>	38 lbs.
<b>AVAILABLE ELECTRONICS</b> <sup>8</sup>	120A Q-Logic 3 EX 75A Q-Logic 2 NE 75A Q-Logic 2 NE+
<b>BATTERY CHARGER</b>	8A Off-board
<b>MOTOR PACKAGES</b>	2-Pole Mid-Wheel 6®
<b>WEIGHT CAPACITY</b> <sup>4</sup>	300 lbs.
<b>BASE WEIGHT</b>	138.6 lbs. with 6" casters 148 lbs. with 8" casters

1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.  
 2) Without front riggings  
 3) Battery capacity may change when adding power positioning systems.  
 4) Weight capacity is dictated by the selection of seating.  
 5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.  
 6) Maximum dimensions of the Manual Recline vary by power base and seating system. Call for details.  
 7) Battery weight may vary +/- 3 lbs. based on manufacturer.  
 8) All Quantum Rehab® units can be ordered with specialty controls.  
 9) Call for availability

# Features

- Standard 6 mph motors offer advanced performance, precise control and responsiveness via innovative single-stage gearbox
- Bold aesthetics
- ATX Suspension for optimum ride and enhanced obstacle climbing
- Available with WC19 crash-test approved occupied transit kit option on tilt & recline and iLevel

# Options

- iLevel® option allows 12" of seat elevation while at walking speed (3.5 mph)
- Q-Logic 3 EX Drive Controls
- TRU-Balance® 3 Seating
- Trilogy vent tray<sup>9</sup>
- Compact vent tray<sup>9</sup>
- Stealth® Seating
- Contoured Seating
- i-Drive® Electronics
- 4" wide flat-free knobby extended traction tires<sup>9</sup>
- LED fender lights (standard with iLevel®)
- 12" x 12" seat sizes with iLevel (SOF only)
- 8" solid casters
- Accu-Trac Advanced Tracking Technology
- Stealth® Positioning Products

# Colors



# Hydrographics

