

# Q-LOGIC 2

DRIVE CONTROL SYSTEM



Connecting Lives



@QuantumRehab



youtube.com/QuantumRehab



facebook.com/QuantumRehab



QUANTUM®

(US) 866-800-2002 • (Canada) 888-570-1113 • [www.quantumrehab.com](http://www.quantumrehab.com)

QLGC-6-11-14

[www.quantumrehab.com](http://www.quantumrehab.com)

# Q-LOGIC 2 DRIVE CONTROL SYSTEM

The Q-Logic 2 Drive Control System by Quantum Rehab® is the most advanced power chair electronics package available. The Q-Logic Series creates an adaptable, expandable control system that serves even the most complex needs while delivering a host of exclusive features.

## FEATURES

- More programming and customization options for a smooth ride
- Larger 2.8" LCD screen
- Built-in ambient light sensor automatically adjusts screen brightness
- Controller armor available in color matching options
- Enhanced horn volume
- Bluetooth access to computer, tablet and phone functions


## Q-Logic 2 EX Enhanced Display



The optional Enhanced Display serves as a separate screen used with specialty drive controls to display system data. With the Enhanced Display's built-in environmental controls, Q-Logic 2 can operate a wide range of appliances that use infrared (IR) remote controls, such as televisions, and serves as a Bluetooth computer mouse emulator.



Drive Screen



Seating Screen

Q-Logic 2 Joystick

## Optional Modules



Q-Logic 2 EX Stand-Alone Joystick (SAJ) / EX Attendant Joystick



Q-Logic 2 EX Joystick for LED lights



Advanced Actuator Module (AAM)



Sip-N-Puff



Environmental Control Unit (ECU)

## Programming Devices



Hand-held Programmer



PC Programmer

## Q-Logic NE Series



Q-Logic NE



Q-Logic NE+

## Non-Expandable



NE and NE+ PG VR2 (4-key) PG VR2 (6-key) PG VR2 (9-key)

## Expandable



Q-Logic 2 EX Joystick Q-Logic 2 EX Joystick with Lights Q-Logic 2 EX Enhanced Display Q-Logic 2 EX Stand Alone Joystick / EX Attendant Joystick Q-Logic 2 Specialty Controls Interface Module

## Multi-Switch Systems



ASL 3-Switch Head Array Switch-It Head Array (3, 4 or 5 switch) 3-Switch Stealth Ultra Head Array Stealth i-Drive™ Tri-Array

## Proportional Drive Controls



ASL Micro Extremity Control Switch-It Proportional Joystick Stealth Mushroom Joystick ASL Rim Control Switch-It Touch Drive

## Single Switch



TASH Pillow Switch TASH Trigger Switch Buddy Button Egg Switch MicroLite Switch Switch-It Uniwobble Switch

\* Either the Q-Logic 2 Enhanced Display or Q-Logic Enhanced Display is still required for Sip N' Puff and single switch driving.



**Q-Logic 2 Specialty Controls Interface Module\***

The Q-Logic 2 Specialty Controls Interface Module offers an enhanced display at a lower cost and is an economical solution to provide specialty controls. On this module, specialty controls run standard through the Q-Logic 2 Drive Control System. Bluetooth® functionality also operates through Q-Logic 2. The module mounts to the TRU-Balance® 3 Power Positioning Systems accessory bar.

## Sip-N-Puff Options



Therafin Whisper-Lite Sip-N-Puff Hardware ASL Sip-N-Puff Head Array Switch-It Sip-N-Puff with Opti-Stop

## Additional Options



ASL Remote Stop Switch Switch-It Remote Stop Switch ASL Remote Attendant Control Joystick Handles