

PRIDE® **MOBILITY SCOOTERS**
LIVE YOUR BEST®

Get **Out**
There!



MOBILITY SCOOTER **CATALOG**

Beamsville, Ontario • 888-570-1113

pridemobility.com





Don't let a change in mobility affect your lifestyle. Live your best with Pride® Mobility Scooters and enjoy sleek, sporty style in a high-performance package.

With an array of models and colors to choose from, Pride Mobility Scooters make it easier for you to get out there and enjoy life.

What's Inside

• Go-Go® Travel Mobility	3-5	• Maxima	10
• Celebrity® Series	6	• Scooter Accessories	11
• Victory® Series	7-8	• Specifications	13-14
• Pursuit® Series	9		



Choosing the right Pride Mobility Scooter

When choosing a Pride Mobility Scooter, specific requirements and conditions can affect your process.

Indoor maneuverability

- Travel Mobility to Mid-Size Luxury scooters offer enhanced maneuverability

Stability

- 4-Wheel scooters provide increased stability

Indoor and outdoor versatility

- Mid-Size and Full-Size Luxury Scooters are great for both indoor and outdoor use

Outdoor capabilities

- Mid-Size and Full-Size Luxury Scooters offer excellent outdoor capabilities

Increased weight capacity

- Mid-Size and Full-Size Luxury Scooters provide increased weight capacity

Easy portability

- Travel Mobility Scooters are easy to disassemble, transport and store

Speed and performance

- Mid-Size and Full-Size Luxury Scooters offer increased speed and performance

GOOGO[®]

TRAVEL MOBILITY



Go-Go[®] Ultra X 3-wheel



Go-Go[®] Ultra X 4-wheel



Go-Go Elite Traveller[®]
4-wheel shown in Blue



Go-Go Elite Traveller[®]
3-wheel

3 & 4-WHEEL SCOOTER

GOOGO[®] ULTRA X

Features

- 260 lbs. weight capacity
- Feather-touch disassembly permits simple frame separation with only one hand into 5 super lightweight pieces for convenient transport
- Maximum speed up to 4 MPH
- Up to 6.9 miles per charge range (3W), 7.2 per charge range (4W)

Color Choices



Red

3 & 4-WHEEL SCOOTER

GOOGO[®] ELITE TRAVELLER[®]

Features

- 300 lbs. weight capacity
- Feather-touch disassembly permits simple frame separation with only one hand into 5 super lightweight pieces for convenient transport
- Two speed settings - high and low (indoor/outdoor)
- Maximum speed up to 4 MPH
- Up to 10.4 miles per charge range (3W), 9.7 per charge range (4W)
- Charger port in tiller

Color Choices

INTERCHANGEABLE
SHROUDS[®]



Red



Blue

The *First Name* in *Travel Mobility*®

3 & 4-WHEEL SCOOTER

GO-GO® LX WITH CTS SUSPENSION

Features

- 300 lbs. weight capacity
- Feather-touch disassembly permits simple frame separation with only one hand into 5 super lightweight pieces for convenient transport
- CTS Suspension (Comfort-Trac front and rear suspension)
- Maximum speed up to 4.4 MPH
- Up to 9.44 miles per charge range (3W), 10.5 miles per charge range (4W)
- Charger port in tiller
- Delta tiller with wraparound handles



Go-Go® LX w/CTS Suspension
3-wheel shown in Blue



Go-Go® LX w/CTS Suspension
4-wheel

Color Choices

INTERCHANGEABLE
SHROUDS*



Red



Blue

Optional



Orion
Yellow



Onyx
Black



Pearl
Pink



Sunburst
Orange

3 & 4-WHEEL SCOOTER

GO-GO® SPORT

Features

- 325 lbs. weight capacity
- Feather-touch disassembly permits simple frame separation with only one hand into 5 super lightweight pieces for convenient transport
- Maximum speed up to 4.7 MPH
- Up to 12.7 miles per charge range (3W), 12 miles per charge range (4W)
- Charger port in tiller
- Delta tiller with wraparound handles



Go-Go® Sport 4-wheel
shown in Blue



Go-Go® Sport 3-wheel

Color Choices

INTERCHANGEABLE
SHROUDS*



Red



Blue

pridemobility.com
pridescooters.com

Pride®
Mobility Products Co.
— LIVE YOUR BEST —

GO·GO®

TRAVEL MOBILITY

3-WHEEL SCOOTER

GO·GO®

FOLDING SCOOTER

Features

- 250 lbs. weight capacity
- Folds in three easy steps
- Maximum speed up to 4 MPH
- Up to 12.5 miles per charge range (lithium ion), 9.3 miles per charge range (sealed lead acid)
- Height adjustable tiller
- Lighting package
- Optional lithium ion battery
- Optional armrests

Color Choices



White

GO·CHAIR

Features

- 300 lbs. weight capacity
- One hand feather-touch disassembly
- 18" x 17", 5-point swivel seat with latch release
- Faster speed - up to 3.7 mph
- Compact turning radius 25.4"
- Dual swing-away, under seat storage bins
- Height & width adjustable armrests

Color Choices



Ruby Red



Sapphire Blue



Rose Quartz



Citrine Yellow



Amber Orange



Pearl White



Note: Transport of this item on commercial aircraft is subject to federal regulations and airline policy. Please consult your airline prior to your flight to ensure your Go-Go® Folding Scooter and its lithium-ion battery is permitted.



Feather-touch disassembly





Celebrity® LR 3-wheel



Celebrity® LR 4-wheel

3 & 4-WHEEL SCOOTER

CELEBRITY® LR

Features

- Deluxe, dual-adjustable delta tiller with wraparound handles
- Ultra-simple, auto-latching lockup for quick disassembly/reassembly
- Deluxe light package with turn signals
- Speed adjustable up to 6.5 mph
- 350 lbs. weight capacity
- Easily accessible tie-down points (for transport of unoccupied scooter)
- Front frame-mounted seat post for optimal stability
- Jumbo front basket and rearview mirror standard
- Larger battery capacity for optimal range

Color Choices



Candy
Apple Red



Onyx
Black

4-WHEEL SCOOTER

CELEBRITY® DX SPORT

Features

- State-of-the-art rear suspension
- Larger rear wheels
- Top speed of 9.5 mph
- Maximum weight capacity of 350 lbs.
- Ultra-simple, auto-latching lockup for quick disassembly/reassembly
- Easily accessible tie-down points (for transport of unoccupied scooter)
- Includes standard angle-adjustable dual bulb curb light
- Front basket and rearview mirror standard

Color Choices



Candy
Apple Red



Onyx
Black



VICTORY[®] SERIES

3 & 4-WHEEL SCOOTER



Victory[®] 10 3-wheel
shown in Blue



VICTORY[®] 10

Features

- 400 lbs. weight capacity
- Feather-touch disassembly permits simple frame separation with only one hand
- Maximum speed up to 5.25 MPH
- Up to 15.5 miles per charge range
- Delta tiller with wraparound handles
- Long lasting, high intensity LED headlight offers optimal pathway illumination
- Easy access tiller-mounted charger port
- Optional 18"x18" high-back seat
- Optional 22"x18" high-back seat

Color Choices



Candy
Apple Red



Viper
Blue



Dual front wheels
for added stability

VICTORY[®] Twin

Features

- 400 lbs. weight capacity
- Feather-touch disassembly permits simple frame separation with only one hand
- Maximum speed up to 5.25 MPH
- Up to 15.5 miles per charge range
- Auto-connecting front to rear cable
- Delta tiller with wraparound handles
- Pride's exclusive black, non-scuffing 10" tires
- Lightweight pinchless hinge seat, eight position rotation, sliders and viscoelastic foam inserts for increased comfort
- Long lasting, high intensity LED headlight offers optimal pathway illumination
- Easy access tiller-mounted charger port
- Optional 18"x18" high-back seat
- Optional 22"x18" high-back seat

Color Choices



Candy
Apple Red



Viper
Blue



The *Ultimate in Style & Performance*®

4-WHEEL SCOOTER

VICTORY[®] 10 LX WITH CTS SUSPENSION

Features

- 400 lbs. weight capacity
- Feather-touch disassembly permits simple frame separation with only one hand
- CTS Suspension (Comfort-Trac front and rear suspension)
- Maximum speed up to 5.6 MPH
- Up to 16.16 miles per charge range
- Delta tiller with wraparound handles
- Black vinyl reclining high-back seat with sliders, height adjustable headrest, flip-up, height and width adjustable armrest
- Full LED lighting package that includes curb and rear lights and directional signals
- Easy access tiller-mounted charger port
- Optional 22" x 18" high-back seat

Color Choices



Candy
Apple Red



Viper
Blue





FULL-SIZE LUXURY

*Outstanding Outdoor
Performance with
Luxurious Features*

4-WHEEL SCOOTER
(ALSO AVAILABLE IN HD VERSION)



Features



- Infinitely-adjustable tiller angle with conveniently located handle
- Full lighting package including directional signals and angle-adjustable, lower headlight
- 13" pneumatic wheels
- Front and rear suspension
- Delta tiller with wraparound handle
- Maximum speed of 7.25 mph/ 6.25 mph with HD version
- Standard jumbo front basket and rearview mirror
- 400 lbs. weight capacity
- Black vinyl reclining high-back seat with sliders, swivel, adjustable armrest (height and width) and height-adjustable headrest

Color Choices



Candy
Apple Red



Onyx
Black

4-WHEEL SCOOTER



Features



- 400 lbs. weight capacity
- Infinitely-adjustable tiller angle with conveniently located handle
- Maximum speed up to 8.5 MPH
- Up to 29.51 miles per charge range (G27), 22.13 miles per charge range (G24)
- Front and rear suspension
- Delta tiller with wraparound handle
- Black vinyl reclining high-back seat with sliders, swivel, adjustable armrest (height and width) and height-adjustable headrest
- 13" low-profile, pneumatic tires for excellent outdoor performance

Color Choices



Candy
Apple Red



Silver



MAXIMA

HEAVY-DUTY

Heavy-Duty Performance at an Exceptional Value

3 & 4-WHEEL SCOOTER

Features

- 500 lbs. weight capacity
- Ultra heavy-duty drivetrain with 110 amp controller
- Maximum speed up to 5.25 MPH
- Easily adjustable tiller
- Up to 18 miles per charge range (3W, NF-22), 9.6 miles per charge range (Power elevating seat, U-1), 16 miles per charge range (4W)
- Durable and stylish rear polyurethane bumper
- User-friendly console



Maxima 3-wheel shown in Blue



Maxima 4-wheel

Color Choices



Candy Apple Red



Viper Blue

A wide selection of Pride Mobility Scooter accessories are available to personalize your Pride Mobility Scooter to better meet your needs and make it uniquely yours.



Cane Clip Attachment



Cup Holder



Front Basket



Rear Basket



Oxygen Holder



RAM® X-Grip® Cell Phone Holder



Wishbone Crutch Holder



Lapbelt



Saddlebag



Walker Holder



Rearview Mirror



Weather Breaker



Weather Cover

DISCLAIMER: Not all accessories are available for all products. Please speak with your Provider for details and availability.

INNOVATIVE TECHNOLOGY

CTS SUSPENSION



Front Suspension



Rear Suspension

Comfort-Trac Suspension (CTS) is the next generation in mobility scooter ride comfort. Pride's **industry first** technology incorporates advanced suspension components for improved absorption over rough surfaces, and a smoother ride experience over various terrains.

FEATHER-TOUCH DISASSEMBLY



Scooter/Travel Mobility disassembly



Go Chair® disassembly

Most Go-Go® Travel Mobility and Pride® Mid-Size Scooters feature **feather-touch** — one-hand disassembly for hassle-free transport and storage.

PRIDE® MOBILITY SCOOTERS

LIVE YOUR BEST®

SPECIFICATIONS

	Front Tires	Rear Tires	Maximum Speed ²	Ground Clearance ^{4, 6, 13}	3-W Turning Radius ⁴	4-W Turning Radius ⁴	3-W Length/Width ¹²	4-W Length/Width ¹²	Seat-to-ground Height Range	(3/4-wheel) Weight ^{3, 4, 14}
TRAVEL MOBILITY										
Go-Chair®	3" x 6" SOLID	3" x 9" SOLID	UP TO 3.7 MPH	1.6" AT ANTI-TIP BRACKET 1.9" AT FRAME	N/A	25.4" WITH FOOT PLATFORM	N/A	33" WITH FOOT PLATFORM/ 22.6"	22"-23"	97 LBS.
Go-Go® Ultra X	2" x 8" SOLID (3W) 2" x 7" SOLID (4W)	2.5" x 8" SOLID	UP TO 4 MPH	2.25" AT MID-FRAME/ 1.25" AT MOTOR	33"	44"	37"/19.5"	39.5"/19.5"	18.13"-20.13"	76 LBS./ 80 LBS.
Go-Go Elite Traveller®	2" x 8" SOLID (3W) 2" x 7" SOLID (4W)	2.5" x 8" SOLID	UP TO 4 MPH	2"	33"	44"	37"/19.5"	39.5"/19.5"	18.13"-20.13"	76 LBS./ 80 LBS.
Go-Go® Sport¹⁵	3" x 8" SOLID (3W) 3" x 9" SOLID (4W)	3" x 9" SOLID	UP TO 4.7 MPH	3" AT MID-FRAME 1.5" AT MOTOR	36"	50.5"	42"/21.25"	43"/21.25"	20.5"-22"	92.8 LBS./ 99.2 LBS.
Go-Go® LX w/ CTS Suspension¹⁵	2" x 8" SOLID (3W) 2" x 7" SOLID (4W)	2.5" x 8" SOLID	UP TO 4.4 MPH	1.5" AT MOTOR	33"	44.75"	39.75"/20.5"	41"/20.5"	20.75"-22.25"	85.4 LBS./ 87.4 LBS.
Go-Go® Folding Scooter	1.6" x 7" SOLID	2" x 8" SOLID	UP TO 4 MPH	2.3"	N/A	32"	N/A	37"/19"	20.5"	N/A/ 51.4 LBS.
VALUE SERIES										
Victory® 10	3.5" x 10" SOLID	3.5" x 10" SOLID	UP TO 5.25 MPH	3" AT MID-FRAME	42"	54"	45.25"/22.25"	47"/22.25"	20.5"-22"	133.5 LBS./ 137.5 LBS.
Victory® Twin¹⁶	Dual 3" x 9" SOLID	3.5" x 10" SOLID	UP TO 5.25 MPH	3"	N/A	42"	N/A	46"/26"	23"-24.5"	N/A/ 150 LBS.
MID-SIZE LUXURY										
Victory® LX w/ CTS Suspension¹⁶	3.5" x 10" SOLID (PNEUMATIC OPTIONAL)	3.5" x 10" SOLID (PNEUMATIC OPTIONAL)	UP TO 5.6 MPH	2.2" AT MOTOR	N/A	59.9"	N/A	47"/22.2"	24"-26"	N/A/ 179 LBS.
Celebrity® LR¹⁶	4" X 10" PNEUMATIC	4" X 10" PNEUMATIC	UP TO 6.5 MPH	2"	44.5"	56.5"	47.5"/24"	48"/24"	22.5"-24"	160 LBS./ 173 LBS.
Celebrity® DX Sport¹⁶	4" x 10" PNEUMATIC	4" x 13" PNEUMATIC	UP TO 9.5 MPH	4.25"	N/A	66.25"	N/A	47.5"/25.25"	22.75"-24.75"	N/A/ 200 LBS.
HEAVY-DUTY										
Maxima¹⁶	3" X 10" SOLID	4" X 10" SOLID	UP TO 5.25 MPH	3W: 4.25" 4W: 3.25"	50.38"	54"	46"/25.5"	46.9"/25.5"	22.5"-24.5"	171 LBS./ 176 LBS.
Pursuit® 16	4" X 13" PNEUMATIC	4" X 13" PNEUMATIC	UP TO 7.5 MPH HD: UP TO 6.25 MPH	4.5"	N/A	67.75" LEFT/ 66.75" RIGHT	N/A	53.25"/25.25"	24.75"-26.25"	N/A/ 225 LBS.
FULL-SIZE LUXURY										
Pursuit® XL¹⁶	13" PNEUMATIC	13" PNEUMATIC	UP TO 8.5 MPH	5.25"	N/A	82.5"	N/A	56"/27.25"	26.25"-27.75"	259.9 LBS.

1. Batteries included.
 2. Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. This specification can be subject to a variance of (+ or -) 10%.
 3. Battery weight not included.
 4. Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3%.
 5. Battery pack is 18 lbs. with 12AH batteries.
 6. At midframe.
 7. Deep cycle (AGM or gel cell type recommended).
 8. 18 AH battery pack weighs 28 lbs.

9. Tested in accordance with ANSI/RESNA, WC Vol 2, Section 4 and ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
 10. Option available through most providers. Not available through Pride Mobility®.
 11. Battery weight may vary based on manufacturer.
 12. Width of base (not including seat).
 13. Does not account for anti-tip clearance.
 14. Includes standard seat weight.
 15. Includes front light.
 16. Includes full lighting package.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed & range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%.

Get Out There!

Seat Type (Standard)	Seat-to deck Height Range	Battery Weight ¹¹	Battery Requirements ⁷	Per Charge Range (Max.) ^{3,9}	Battery Charger	Maximum Weight Capacity	Weight of Heaviest Piece ⁴	Warranty
COMPACT, FOLDABLE (18" X 17") 5-POINT SWIVEL WITH LATCH RELEASE	N/A	31 LBS. BATTERY PACK	(2)18AH	UP TO 8.74 MILES	OFF-BOARD, 2A	300 LBS.	36 LBS. (REAR)	5-YEAR LIMITED ON FRAME; 2-YEAR LIMITED ON ELECTRONICS; 2-YEAR LIMITED ON DRIVETRAIN
COMPACT, FOLDABLE (17" X 17")	15.75"-17.25"	9 LBS. EA. ⁵	(2)12AH ¹	3W: UP TO 6.9 MILES 4W: UP TO 7.2 MILES	OFF-BOARD, 2A	260 LBS.	29.5 LBS. (3-WHEEL) (FRONT SECTION) 34 LBS. (4-WHEEL) (FRONT SECTION)	LIFETIME LIMITED ON FRAME; 2-YEAR LIMITED ON ELECTRONICS; 2-YEAR LIMITED ON DRIVETRAIN
COMPACT, FOLDABLE (17" X 17")	15.75"-17.25"	14 LBS. EA. ⁸	(2)18AH ¹	3W: UP TO 10.4 MILES ⁵ 4W: UP TO 9.7 MILES ¹⁶	OFF-BOARD, 2A	300 LBS.	29.5 LBS. (3-WHEEL) (FRONT SECTION) 34 LBS. (4-WHEEL) (FRONT SECTION)	LIFETIME LIMITED ON FRAME; 2-YEAR LIMITED ON ELECTRONICS; 2-YEAR LIMITED ON DRIVETRAIN
COMPACT, FOLDABLE (17" X 17")	17"-18"	14 LBS. EA. ⁸	(2)18AH ¹	3W: UP TO 12.7 MILES 4W: UP TO 12 MILES	OFF-BOARD, 2A	325 LBS.	35 LBS. (3-WHEEL) (FRONT SECTION) 44 LBS. (4-WHEEL) (FRONT SECTION)	LIFETIME LIMITED ON FRAME; 2-YEAR LIMITED ON ELECTRONICS; 2-YEAR LIMITED ON DRIVETRAIN
COMPACT, FOLDABLE (17" X 17")	16.25"-17.25"	14 LBS. EA. ⁸	(2)18AH ¹	3W: UP TO 9.44 MILES 4W: UP TO 10.5 MILES	OFF-BOARD, 2A	300 LBS.	34 LBS. (3-WHEEL) (FRONT SECTION) 36.6 LBS. (4-WHEEL) (FRONT SECTION)	LIFETIME LIMITED ON FRAME; 2-YEAR LIMITED ON ELECTRONICS; 2-YEAR LIMITED ON DRIVETRAIN
COMPACT, FOLDABLE (16" X 13")	13"	9 LBS. EA. (SLA) 5 LBS. EA. (LITHIUM-ION)	(2) 12 AH SLA OR (1) 11.5 AH LITHIUM-ION	UP TO 12.5 MILES (LITHIUM-ION); UP TO 9.3 MILES (SLA)	2-AMP OFF-BOARD (SLA) OR 2.5-AMP OFF-BOARD (LITHIUM-ION)	250 LBS.	51.4 LBS. (UNIT DOES NOT DISASSEMBLE - UNIT WEIGHT WITHOUT BATTERY PACK)	LIFETIME LIMITED ON FRAME; 2-YEAR LIMITED ON ELECTRONICS; 2-YEAR LIMITED ON DRIVETRAIN
FOLDABLE, MOLDED PLASTIC (18" X 17")	17.5"-19"	24.5 LBS. EA/ 32.5 LBS.	(2) U1 ¹ OR (2) 40AH	UP TO 15.5 MILES	OFF-BOARD, 3.5A (U-1) OFF-BOARD, 5A (40AH)	400 LBS.	54.5 LBS. (3-WHEEL) (FRONT SECTION) 81 LBS. (4-WHEEL) (FRONT SECTION)	LIFETIME LIMITED ON FRAME; 2-YEAR LIMITED ON DRIVETRAIN; 2-YEAR LIMITED ON ELECTRONICS
FOLDABLE, MOLDED PLASTIC (18.5" X 17.25")	17.5"-19"	24.5 LBS. EA/ 32.5 LBS. EA.	(2) U1 ¹ OR (2) 40AH	UP TO 15.5 MILES	OFF-BOARD, 3.5A (U-1) OFF-BOARD, 5A (40AH)	350 LBS.	55 LBS. (FRONT SECTION)	LIFETIME LIMITED ON FRAME; 2-YEAR LIMITED ON DRIVETRAIN; 2-YEAR LIMITED ON ELECTRONICS
RECLINING HIGH-BACK SEAT WITH SLIDERS, SWIVEL, ADJUSTABLE ARMREST (HEIGHT AND WIDTH) AND HEIGHT-ADJUSTABLE HEADREST (18" X 18")	17"-19"	24.5 LBS. EA/ 32.5 LBS. EA.	(2) U1 OR (2) 40AH	UP TO 16.16 MILES	OFF-BOARD, 5A	400 LBS.	77.6 LBS. (FRONT SECTION)	LIFETIME LIMITED ON FRAME; 2-YEAR LIMITED ON DRIVETRAIN; 2-YEAR LIMITED ON ELECTRONICS
RECLINING HIGH-BACK SEAT WITH SLIDERS, SWIVEL, ADJUSTABLE ARMREST (HEIGHT AND WIDTH) AND HEIGHT-ADJUSTABLE HEADREST (18" X 18")	16.5"-18.5"	36.5 LBS. EA.	(2) NF-22	3W: UP TO 17.99 MILES 4W: UP TO 18.97 MILES	OFF-BOARD, 5A	350 LBS.	61 LBS. (3-WHEEL) (FRONT SECTION) 73 LBS. (4-WHEEL) (FRONT SECTION)	LIFETIME LIMITED ON FRAME; 2-YEAR LIMITED ON DRIVETRAIN; 2-YEAR LIMITED ON ELECTRONICS
RECLINING HIGH-BACK SEAT WITH SLIDERS, SWIVEL, ADJUSTABLE ARMREST (HEIGHT AND WIDTH) AND HEIGHT-ADJUSTABLE HEADREST (18" X 18")	16.75"-18.75"	24.5 LBS. EA/ 32.5 LBS. EA.	(2) U1 OR (2) 40AH	UP TO 11 MILES	OFF-BOARD, 5A	350 LBS.	83 LBS. (FRONT SECTION)	LIFETIME LIMITED ON FRAME; 2-YEAR LIMITED ON DRIVETRAIN; 2-YEAR LIMITED ON ELECTRONICS
DELUXE CONTOURED MEDIUM-BACK (22" X 18")	3W: 22.3"-22.6" 4W: 24.3"-24.6"	24.5 LBS. EA/ 38 LBS. EA.	(2) U-1 (POWER ELEVATING SEAT) OR (2) NF-22 ¹	3W: UP TO 18 MILES (NF-22) UP TO 9.6 MILES (U-1) 4W: UP TO 16 MILES	OFF-BOARD, 5A	500 LBS.	76 LBS. (3-WHEEL) (REAR SECTION) 76 LBS. (4-WHEEL) (REAR SECTION)	LIFETIME LIMITED ON FRAME; 2-YEAR LIMITED ON DRIVETRAIN; 2-YEAR LIMITED ON ELECTRONICS
RECLINING HIGH-BACK SEAT WITH SLIDERS, SWIVEL, ADJUSTABLE ARMREST (HEIGHT AND WIDTH) AND HEIGHT-ADJUSTABLE HEADREST (18" X 18")	17"-18.75"	38 LBS. EA/ 53 LBS. EA.	(2) NF-22 OR (2) GROUP 24	UP TO 25 MILES (WITH GROUP 24 BATTERIES)	OFF-BOARD, 8A	400 LBS.	95 LBS. (FRONT SECTION)	LIFETIME LIMITED ON FRAME; 2-YEAR LIMITED ON DRIVETRAIN; 2-YEAR LIMITED ON ELECTRONICS
RECLINING HIGH-BACK SEAT WITH SLIDERS, SWIVEL, ADJUSTABLE ARMREST (HEIGHT AND WIDTH) AND HEIGHT-ADJUSTABLE HEADREST (18" X 18")	21"-22.5"	70 LBS. EA. (G31) 75.5 LBS. EA. (G27) 53 LBS. EA. (G24)	(2) GROUP 31 ¹⁰ , 10AH (2) GROUP 27 100AH OR (2) GROUP 24 75 AH	G31: UP TO 31.83 MILES G27: UP TO 29.51 MILES G24: UP TO 22.13 MILES	OFF-BOARD, 8A	400 LBS.	205.5 LBS. (MAIN FRAME)	LIFETIME LIMITED ON FRAME; 2-YEAR LIMITED ON DRIVETRAIN; 2-YEAR LIMITED ON ELECTRONICS



PRIDE® MOBILITY SCOOTERS
LIVE YOUR BEST®

AUTHORIZED PRIDE® PROVIDER

PROVIDER INFO