



	SOLUTION <sup>®</sup> 1	SOLUTION <sup>®</sup>	SPECTRUM <sup>®</sup> GEL	SPECTRUM <sup>®</sup> FOAM	TRU-COMFORT 2 SEAT CUSHION
<b>FEATURES</b>	<ul style="list-style-type: none"> <li>Deep contoured, high-density molded foam base</li> <li>Coccyx and sacral relief</li> <li>Silicone soft gel top layer with airflow design</li> <li>Moisture wicking four-way stretch incontinence cover</li> </ul>	<ul style="list-style-type: none"> <li>Lateral and pre-ischial shelf contoured into high-density molded foam base</li> <li>Twin-cell, aqueous gel insert under at-risk bony areas</li> <li>Moisture wicking four-way stretch cover</li> <li>Inner liner for incontinence and to reduce shear</li> <li>Viscoelastic total top layer</li> </ul>	<ul style="list-style-type: none"> <li>Lateral and pre-ischial shelf contoured into high-density molded foam base</li> <li>Twin-cell, lightweight fluid insert under at-risk bony areas</li> <li>Moisture wicking four-way stretch cover</li> <li>Inner liner for incontinence and to reduce shear</li> <li>Viscoelastic total top layer</li> <li>Solid seat insert included</li> </ul>	<ul style="list-style-type: none"> <li>Lateral and pre-ischial shelf contoured into high-density molded foam base</li> <li>Multi-layer, viscoelastic foam insert under at-risk bony areas</li> <li>Moisture wicking four-way stretch cover</li> <li>Inner liner for incontinence and to reduce shear</li> <li>Viscoelastic total top layer</li> <li>Solid seat insert included</li> </ul>	<ul style="list-style-type: none"> <li>Deep contoured, high-density molded foam seat base</li> <li>Viscoelastic total top layer</li> <li>Ergonomically contoured back cushion</li> <li>Moisture wicking four-way stretch incontinence cover</li> <li>Additional 1" of viscoelastic foam in seat cushion for maximum comfort and pressure management</li> <li>Pelvic well</li> </ul>
<b>BENEFITS</b>	<ul style="list-style-type: none"> <li>Provides lateral support for clients in need of moderate positioning</li> <li>Designed for clients who need mild to moderate skin protection</li> <li>Reduces shear, moisture build-up and heat build-up</li> <li>Promotes better circulation</li> </ul>	<ul style="list-style-type: none"> <li>Provides lateral and anti-thrust support for clients in need of moderate positioning</li> <li>Designed for clients who need mild to moderate skin protection</li> <li>Reduces shear and moisture build-up</li> <li>Promotes better circulation</li> </ul>	<ul style="list-style-type: none"> <li>Provides lateral and anti-thrust support for clients in need of moderate positioning</li> <li>Provides excellent stability</li> <li>Designed for clients who need moderate to high skin protection</li> <li>Reduces shear and moisture build-up</li> <li>Promotes better circulation</li> </ul>	<ul style="list-style-type: none"> <li>Provides lateral and anti-thrust support for clients in need of moderate positioning</li> <li>Provides excellent stability</li> <li>Designed for clients who need moderate to high skin protection</li> <li>Reduces shear and moisture build-up</li> <li>Promotes better circulation</li> </ul>	<ul style="list-style-type: none"> <li>Provides lateral and anti-thrust support for clients in need of mild to moderate positioning</li> <li>Designed for clients who need mild to moderate skin protection</li> <li>Helps promote good posture</li> <li>Reduces moisture build-up</li> <li>Promotes better circulation</li> <li>High contour side for femur positioning</li> <li>Reduces pressure (ischial tuberosities)</li> </ul>
<b>COVER</b>	Moisture wicking four-way stretch incontinence cover	Moisture wicking four-way stretch cover	Moisture wicking four-way stretch cover	Moisture wicking four-way stretch cover	Moisture wicking four-way stretch incontinence cover
<b>SOLID SEAT INSERT</b>	Available option	Available option	Included	Included	Available option
<b>WATER-RESISTANT INNER LINER</b>	Available option	Included	Included	Included	Available option
<b>AVAILABLE SIZES*</b>	14"x14" - 24"x24"	12"x12" - 24"x22"	16"x16" - 22"x22"	16"x16" - 28"x22"	16" X 14" - 20" X 24" 22" X 16" - 22" X 24"
<b>WEIGHT CAPACITY</b>	350 lbs.	350 lbs.	350 lbs.	650 lbs.	450 lbs.
<b>OVERALL HEIGHT</b>	2.5"	3.5"	3"	3"	3.25" - 5"
<b>MSRP</b>	\$285-\$335	\$420-\$475	\$580-\$630	\$580-\$630	\$420-\$525
<b>WARRANTY</b>	2-year limited	2-year limited	2-year limited	2-year limited	2-year limited

\* Custom sizes available

The information contained herein is correct at the time of publication. We reserve the right to alter specifications without prior notice.



### Power Recline

- 85° - 175° of recline
- 300 lbs. weight capacity
- Does not increase seat-to-floor height
- Standard shear reduction and cantilever armrests
- Available in 12 back shroud colors
- Available with TRU-Balance<sup>®</sup> 3 Power Tilt and TRU-Balance<sup>®</sup> 3 Power Adjustable Seat Lift



### Heavy Duty Power Tilt & Recline

- 450 lbs. weight capacity
- 50° of tilt
- 85° - 175° of recline
- Available in 10 back shroud colors



### Power Tilt

- 50° of tilt
- Lower standard seat-to-floor height
- Unmatched adjustment and durability
- Maintains center of gravity shift
- Available in 12 back shroud colors



### Static Seating (for Edge 2.0 only)

- Easily adjustable and versatile design that can grow with the client
- Easy customization
- 300 lb. weight capacity
- Accepts a wide variety of after-market seats and backs
- 20° of limited recline adjustment (-5 / +15) via gas strut within back cover
- Available in 12 back shroud colors



### 10" Power Adjustable Seat Lift with iLevel<sup>®</sup>

- iLevel<sup>®</sup> technology enables the unit to drive up to 3.5 mph while elevated at 10" of lift
- Offers 10" of power adjustable seat lift
- 300 lbs. weight capacity
- Available in 12 back shroud colors



### Manual Recline (for Edge 2.0 only)

- Solid curved back pan and 2.5" raised pivot to reduce shear forces
- Attendant activated recline seat
- Width: 12" - 24"; depth: 12" - 24"
- 300 lb. weight capacity
- 175° of recline
- Can be ordered with TRU-Balance<sup>®</sup> 3 Tilt and TRU-Balance<sup>®</sup> 3 Static Seating
- Available in 12 back shroud colors



END-TO-END SOLUTIONS



TRU-Balance® 3 color options



	TRU-COMFORT	TRU-COMFORT FOR TRU-BALANCE® 3	SYNERGY® BACK
<b>FEATURES</b>	<ul style="list-style-type: none"> <li>Deep contoured, high-density molded foam seat base</li> <li>Viscoelastic total top layer</li> <li>Ergonomically contoured back cushion</li> <li>Moisture wicking four-way stretch incontinence cover</li> <li>Improved design makes electronics more accessible</li> <li>Each back has 1" of height adjustment</li> </ul>	<ul style="list-style-type: none"> <li>Four-way stretch, recovery, fluid-resistant, breathable cover</li> <li>Additional 1" of viscoelastic foam in back cushion for maximum comfort</li> <li>Additional bolster kit for positioning of back cushion and client</li> <li>Compatible with both TRU-Balance® 3 Power Positioning Systems and the Static Seat Frame System</li> </ul>	<ul style="list-style-type: none"> <li>High-density foam base with viscoelastic foam top layer to prevent bottom out</li> <li>Wraparound cover keeps back cushion from shifting position</li> <li>Foam 1" above and below pan to prevent user contact with the metal pan</li> <li>Available full-length side accessory tracks, for ease of adding accessories such as laterals</li> <li>Available in most sizes from 10"W x 13"H up to 32"W x 25"H</li> </ul>
<b>BENEFITS</b>	<ul style="list-style-type: none"> <li>Provides lateral and anti-thrust support for clients in need of mild to moderate positioning</li> <li>Designed for clients who need mild to moderate skin protection</li> <li>Helps promote good posture</li> <li>Reduces moisture build-up</li> <li>Promotes better circulation</li> </ul>	<ul style="list-style-type: none"> <li>Adjustable positioning and pressure management</li> <li>4" of width adjustability in each back shell</li> <li>Independently adjustable from the seat width</li> <li>Can be mounted to the TB3 back spine system for a more integrated look</li> <li>Provides enhanced comfort and support</li> </ul>	<ul style="list-style-type: none"> <li>Deeper contour for stability positioning and comfort</li> <li>Available for use on all Synergy® Seating Systems</li> <li>Contoured design with 1.5" thick back pad, adjustable height and accepts lateral supports</li> </ul>
<b>COVER</b>	Moisture wicking four-way stretch incontinence cover	Four-way stretch cover	Four-way stretch cover
<b>AVAILABLE BACK HEIGHTS</b>	22" - 25"	14" - 23"	9" - 21"
<b>AVAILABLE BACK WIDTHS</b>	16" - 22"	12" - 28"	10" - 32"
<b>WEIGHT CAPACITY</b>	300 lbs.	450 lbs.	650 lbs.
<b>MSRP</b>	\$730-\$790	\$835	\$555-\$870
<b>WARRANTY</b>	2-year limited	2-year limited	2-year limited


\* Width for each back can be adjusted from 1/2" under to 1-1/2" over the specified size.



# END-TO-END SOLUTIONS

## The Most Complete Selection of Solutions for the Broadest Range of Needs

Synergy® Cushions and Backs from Quantum Rehab® offer an excellent selection of pressure management and positioning solutions available in a range of sizes that accommodate pediatric to bariatric patients. Each Synergy cushion and back uses unique components to deliver specialized pressure management and maximize comfort.



Beamsville, Ontario ▪ 888-570-1113 ▪ [www.quantumrehab.com](http://www.quantumrehab.com)

