

VICTORY[®] 10 LX WITH CTS SUSPENSION

4-WHEEL SCOOTER

THE ULTIMATE MID-SIZE LUXURY MOBILITY SCOOTER

The **Victory[®] 10 LX with CTS Suspension** Mobility Scooter features CTS (Comfort-Trac Suspension) Suspension and brings a new level of performance and features to mid-size scooters. Featuring CTS independent front and rear suspension, sleek style, feather-touch disassembly and standard lighting, the 4-wheel Victory 10 LX with CTS Suspension is the ultimate mid-size luxury scooter.



Front Suspension



Rear Suspension

PRIDE[®] MOBILITY SCOOTERS
LIVE YOUR BEST[®]

www.pridemobility.com

SPECIFICATIONS

VICTORY® 10 LX WITH CTS SUSPENSION 4-wheel Model # S710LX

Weight capacity:	400 lbs.
Maximum speed ¹ :	Up to 5.6 mph
Ground clearance ² :	2.2" at motor
Turning radius ² :	59.9"
Overall length ² :	47"
Overall width ² :	22.2"
Seat-to-floor height:	17.7"
Front tires ² :	3.5" x 10" solid (standard); 10" pneumatic (optional)
Rear tires ² :	3.5" x 10" solid (standard); 10" pneumatic (optional)
Suspension:	Front and rear
Range per charge ^{1,4} :	Up to 16.16 miles
Maximum allowable incline:	6°
Total weight without batteries ⁵ :	156.2 lbs.
Heaviest piece when disassembled:	77.6 lbs. (front section)
Standard seating:	Type: CRS14 - High-back with limited recline, seat slide and swivel with height, width and angle-adjustable armrests Weight: 52 lbs. Material: Black Vinyl Dimensions: 18" width, 17.5" depth
Drivetrain:	Rear-wheel drive, sealed transaxle, 24-volt DC motor
Brakes:	Regenerative and electromechanical
Battery requirements ⁵ :	(2) 12-volt, deep cycle Size: U-1 or 40 AH Weight ³ : (U-1) 24.5 lbs. each (40 AH) 32.5 lbs. each
Battery charger:	5 amp, off-board charger

Color:



Candy Apple Red



Viper Blue

¹ Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
² Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
³ AGM or gel cell type required.
⁴ Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
⁵ Battery weight may vary based on manufacturer.

NOTE: All specifications subject to change without notice.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

FEATURES

- CTS Suspension (Comfort-Trac front and rear independent suspension)
- Sleek, modern design
- Exclusive black, non-scuffing tires
- Feather-touch disassembly permits simple frame separation with only one hand
- Modular design for easy serviceability
- Delta tiller with wraparound handles
- Easily accessible charger port in tiller
- Full LED lighting package that includes curb and rear lights and directional signals
- Courtesy light under tiller
- Storage compartment below tiller
- Enhanced flush-mounted horn buttons
- Auto-connecting front-to-rear cable
- Front tiller-mounted basket

OPTIONS

- Forearm crutch holder
- Cane/crutch holder
- Oxygen tank holder
- Cup holder
- Rear basket
- Safety flag
- Walker holder
- Saddlebag
- Rearview mirror
- Cell phone holder

Pride® Lifts & Ramps

Compatible Lifts:

- Backpacker® Plus
- Olympian AutoLift
- Outlander



PRIDE® MOBILITY SCOOTERS
LIVE YOUR BEST®

Beamsville, Ontario • 888-570-1113
www.pridemobility.com