

# ENHANCE

your Active Lifestyle  
with added Comfort  
and Performance



4-Wheel

**VICTORY® 10**



3-Wheel shown with  
optional 22" wide high-back seat

The **Victory® 10** series offers added style, comfort and performance for your active lifestyle. Standard features include LED pathway light, LED battery meter, feather touch disassembly, wraparound delta tiller, 10" solid tires and a front basket. A 400 lb. weight capacity, a per charge range of up to 16 miles and a maximum speed up to 5.3 mph compliment the versatility of the Victory 10. Available in Viper Blue and Candy Apple Red.

PRIDE® **MOBILITY SCOOTERS**  
LIVE YOUR BEST®

| [pridemobility.com](http://pridemobility.com) |





## Specifications 3W:SC610 4W:SC710

Weight Capacity	400 lbs.
Maximum Speed <sup>1</sup>	Up to 5.3 mph
Ground Clearance <sup>3</sup>	3W: 2.2" at motor, 1.26" at anti-tip 3.4" at side deck 4W: 2" at motor, 1.5" at anti-tip 3.1" at side deck
Turning Radius <sup>3</sup>	3W: 42" 4W: 54"
Overall Length <sup>3</sup>	3W: 45.25" 4W: 47"
Overall Width <sup>3</sup>	22.25"
Seat-to-Ground Height	20.5"-22"
Seat-to-Deck Height	17.5"-19"
Front Tires	3.5" x 10" solid
Rear Tires	3.5" x 10" solid
Suspension	N/A
Range Per Charge <sup>1,2</sup>	Up to 16 miles
Total Weight Without Batteries <sup>3,7</sup>	3W: 133 lbs. 4W: 137 lbs.
Heaviest Piece When Disassembled <sup>3</sup>	3W: 54 lbs. (front section) 4W: 61 lbs. (front section)
Standard Seating	Type: Foldable molded plastic Weight: 28 lbs. Material: Black Vinyl Dimensions: 18" width, 17" depth
Drive System	Rear-wheel drive sealed transaxle with 24-volt DC motor
Dual Braking System	Electronic, regenerative and electromechanical
Battery Requirements <sup>3,4,5</sup>	Size: (2) U-1 standard, Weight U-1: 24 lbs. each Size: (2) 40 AH (optional), Weight 40 AH: 32 lbs. each
Battery Charger	3.5-amp, off-board (U-1), 5-amp, off-board (40 AH)
Colors	 

Candy Apple Red Viper Blue

- Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.
- Tested in accordance with ANSI/RESNA, MC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculation based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
- AGM or gel cell type required.
- Battery weight may vary based on battery manufacturer.
- Includes standard seat weight.

## Features

- A standard front basket provides ample storage
- Feather-touch disassembly permits simple frame separation with only one hand for easy transport and storage
- Charger port located on the tiller enables you to conveniently charge your scooter
- Delta tiller provides ergonomic wraparound handles enables easy-to-use left or right hand operation and the ability to rest your wrist on the handle for added comfort
- Easy-to-read LED battery meter
- Exclusive stylish, lightweight, non-scuffing, black, low-profile tires
- Long-lasting, maximum-intensity LED headlight offers optimal pathway illumination
- Lightweight pinchless hinge seat eight position rotation, sliders and viscoelastic foam inserts for increased comfort
- Easy-grip tiller adjustment knob

## Accessories

- Cane clip attachment
- Cane/crutch holder (*single, double*)
- Cell phone holder
- Cup holder
- Forearm crutch holder
- 50", 60", 70" Lap belt
- Oxygen tank holder
- Rear basket
- Rearview mirror
- Saddlebag
- Safety flag
- Walker holder
- Weather cover
- XLR USB charger

## Options

- 40 AH batteries with 5-amp, off board charger

## Optional Seating

- Limited recline, black vinyl, high back with height adjustable headrest, flip-up height and width adjustable armrests
  - 18" x 18" (sliders) • 20" x 20" (no sliders) • 22" x 18" (no sliders)
- Deluxe contour high back black vinyl with height adjustable headrest, flip-up height and width adjustable armrest (no sliders)
  - 18" x 18" • 18" x 20" • 20" x 20" • 22" x 20"



Find a compatible lift for your vehicle online:  
[PRIDEMOBILITY.COM/LIFT-FINDER](http://PRIDEMOBILITY.COM/LIFT-FINDER)

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for scooter only. Overall measurements will vary based on seating and accessory selections.

