



SPECIFICATIONS	Victory® LX with CTS Suspension	Victory® 10	Victory® Twin
Max. weight capacity:	400 lbs.	400 lbs.	350 lbs.
Maximum speed: ¹	Up to 5.6 mph	Up to 5.25 mph	Up to 5.25 mph
Ground clearance: ²	2.2" (at motor)	3" (at midframe)	3" (at midframe)
Turning radius: ²	59.9"	3W: 42" / 4W: 54"	42"
Overall length: ²	47"	3W: 45.25" / 4W: 47"	46"
Overall width: ²	22.2" (at wheels)	22.25" (at wheels)	22.25" (at wheels)
Battery charger:	Off-board, 5 amp	Off-board, 3 amp (U-1 batteries) Off-board, 5 amp (40 AH batteries) ¹	Off-board, 3 amp (U-1 batteries) Off-board, 5 amp (40 AH batteries) ¹
Per-charge range: ⁵	Up to 16.16 miles	Up to 11.1 miles (w/ U-1 batts.)	Up to 11.1 miles (w/ U-1 batts.)
Battery requirement: ³	(2) U-1 or 40 AH	(2) U-1 (included) or 40 AH	(2) U-1 (included) or 40 AH
Battery weight: ⁴	24.5 lbs. each (U-1) 32.9 lbs. each (40 AH)	24.5 lbs. each (U-1) 32.5 lbs. each (40 AH)	24.5 lbs. each (U-1) 32.5 lbs. each (40 AH)
Tire size:	Front: 10" solid standard (pneumatic optional) Rear: 10" solid standard (pneumatic optional)	Front: 10" solid Rear: 10" solid	Dual 9" solid 10" solid
Total weight without batteries: ^{2, 6}	156.2 lbs.	3W: 131.5 lbs./4W: 132.5 lbs.	149 lbs.
Heaviest piece when disassembled: ²	77.6 lbs. (front section)	3-wheel: 54.5 lbs. (front section) 4-wheel: 61 lbs. (front section)	53.5 lbs. (front section)

¹Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).

²Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).

³AGM or gel cell type required.

⁴Battery weight may vary based on manufacturer.

⁵Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance test conducted at maximum weight capacity.

⁶Includes standard seat weight.

NOTE: All specifications subject to change without notice.

VICTORY® SERIES

VICTORY® 10
LX WITH CTS SUSPENSION



VICTORY® 10



VICTORY® Twin



Pride
Mobility Products Co.

— LIVE YOUR BEST —

Pride
Mobility Products Co.

— LIVE YOUR BEST —

VICTORY[®] 10

LX WITH CTS SUSPENSION



FEATURES

- CTS Suspension (Comfort-Trac front and rear independent suspension)
- Sleek, modern design
- Exclusive black, non-scuffing tires
- Feather-touch disassembly permits simple frame separation with only one hand
- Modular design for easy serviceability
- Delta tiller with wraparound handles
- Easily accessible charger port in tiller
- Full LED lighting package that includes curb and rear lights and directional signals
- Courtesy light under tiller
- Storage compartment below tiller
- Enhanced flush-mounted horn buttons
- Auto-connecting front-to-rear cable
- Front tiller-mounted basket

VICTORY[®] 10



FEATURES

- Available in 3 & 4-wheel models
- 4-pole motor standard
- Exclusive one-hand feather-touch disassembly provides the ultimate quick transport system
- Auto-connecting front-to-rear harness
- Long-lasting, high-intensity LED headlight offers optimal pathway illumination
- Delta tiller with wraparound handles
- Push/pull throttle control standard
- Crushproof basket
- LED battery meter
- Easy access, tiller-mounted charger port for standard dual voltage off-board charger
- Easy-grip tiller adjustment knob
- Exclusive stylish, lightweight, solid, non-scuffing, black, low-profile wheels
- Lightweight seat with viscoelastic foam inserts for increased comfort
- Seat includes pinchless hinges and 8-position seat rotation and sliders
- U-1 batteries included in unit

VICTORY[®] Twin



FEATURES

- Dual front wheels for added stability
- 4-pole motor standard
- Exclusive one-hand feather-touch disassembly provides the ultimate quick transport system
- Auto-connecting front-to-rear harness
- Full lighting package, including high-visibility automatic rear brake light system and turn signals
- Delta tiller with wraparound handles
- Push/pull throttle control standard
- Crushproof basket
- LED battery meter
- Easy access, tiller-mounted charger port for standard dual voltage off-board charger
- Easy-grip tiller adjustment knob
- Exclusive stylish, lightweight, solid, non-scuffing, black, low-profile wheels
- Seat includes pinchless hinges and 8-position seat rotation and sliders
- U-1 batteries included in unit