

FOOTNOTES

1. Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
2. Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition.
This specification can be subject to a variance of (+ or -) 10%.
3. Battery weight not included.
4. Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3%.
5. Battery pack is 19.5 lbs. with 12AH batteries.
6. Deep cycle (AGM or gel cell type recommended).
7. 18 AH battery pack weighs 29.5 lbs.
8. Tested in accordance with ANSI/RESNA, WC Vol 2, Section 4 and ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
9. Battery weight may vary based on manufacturer.
10. Width of base (not including seat).
11. Includes standard seat weight.

12. Must be billed to Medicare as either a beneficiary requested upgrade with a properly executed ABN, or as an upgrade at no additional charge to the beneficiary.
13. Light package.
14. Interchangeable shrouds (all shroud panels included).
15. Results derived from theoretical calculation based on battery specifications and drive system performance. Test conducted at maximum weight capacity.



END-TO-END



	i-Go™ Non-Medical Device*	Go Go® Super Portable Non-Medical Device*	Go Go® Carbon FDA Class II Medical Device¹	Go Go® Endurance Li FDA Class II Medical Device¹	Go Go® Ultra X FDA Class II Medical Device¹	Go Go Elite Traveller® FDA Class II Medical Device¹	Go Go Elite Traveller® 2 FDA Class II Medical Device¹	Go Go® Sport FDA Class II Medical Device¹	Revo® 2.0 FDA Class II Medical Device¹	Zero Turn 10 FDA Class II Medical Device¹	Victory® 9 FDA Class II Medical Device¹	Victory® 10 FDA Class II Medical Device¹	Victory® 10.2 FDA Class II Medical Device¹	Victory® LX Sport¹³ FDA Class II Medical Device¹	Pursuit 2 FDA Class II Medical Device¹	PX4 FDA Class II Medical Device¹	Maxima¹³ FDA Class II Medical Device¹	
MODEL NUMBER	SC20	SC15	SC25M	S54XLIT	S39 (3-Wheel) S49 (4-Wheel)	SC40E (3-Wheel)	SC442E	S73 (3-Wheel) S74 (4-Wheel)	S66 (3-Wheel) S67 (4-Wheel)	S710ZT	SC609 (3-Wheel) SC709 (4-Wheel)	S610 (3-Wheel) S710 (4-Wheel)	S6102	S710LXW	SC7132	SC134	SC901 (3-Wheel) SC941 (4-Wheel)	
WEIGHT CAPACITY	275 lbs.	300 lbs.	300 lbs.	325 lbs.	300 lbs.	300 lbs.	300 lbs.	325 lbs.	400 lbs.	400 lbs.	300 lbs.	400 lbs.	400 lbs.	400 lbs.	400 lbs.	500 lbs.	500 lbs.	
SPEED (MAXIMUM)²	Up to 4 mph	Up to 4.5 mph	Up to 3.7 mph	Up to 4.1 mph	3-Wheel: Up to 4.3 mph 4-Wheel: Up to 4.1 mph	Up to 4 mph	Up to 4.2 mph	Up to 4.7 mph	Up to 5.2 mph	Up to 7.3 mph	Up to 5.3 mph	Up to 5.3 mph	Up to 5.2 mph	Up to 8 mph	Up to 9.3 mph	Up to 9.2 mph	Up to 5.8 mph	
GROUND CLEARANCE⁴	4.2" at mid-frame	1.85" at motor 3.7" at mid-deck	1.5" at motor	1.5" at motor 3.25" at side deck	1.25" at motor 3.2" at mid-frame	0.75" at anti-tip 2.25" at mid-frame	1.34" at transaxle	1.5" at motor 3" at mid-frame	1.25" at motor	2" at motor 3" at mid frame	3-Wheel: 1.4" at anti-tip 1.8" at transaxle 2" at mid-frame 4-Wheel: 1.25" at anti-tip 1.75" at motor 2.5" at mid-frame	3-Wheel: 1.26" at anti-tip 2.2" at motor 3" at mid-frame 4-Wheel: 2" at motor 1.5" at anti-tip 3" at mid-frame	2.4" at transaxle 2.1" at anti-tip 3" at side deck	2.2" at motor 3.15" at rear mid-frame	3.78" at transaxle 4.6" at mid-frame	3.5" at motor 4.4" at mid-frame	3-Wheel: 2.5" at motor 3" at mid-frame 4-Wheel: 2.25" at motor	
3-WHEEL TURNING RADIUS⁵	45.9"	N/A	N/A	N/A	36.5"	33"	N/A	36"	40"	N/A	40.25"	42"	42"	N/A	N/A	N/A	49"	
4-WHEEL TURNING RADIUS⁵	N/A	51.62"	53.15	48.25"	52.25"	N/A	37"	50.5"	52"	43"	51.5"	54"	N/A	61.3"	69.3"	82.3"	53"	
3-WHEEL LENGTH/WIDTH⁴,¹⁰	44"/20.3"	N/A	N/A	N/A	37"/19.25"	37"/19.5"	N/A	42"/21.25"	46"/20.75"	N/A	43.5"/22.25"	45.25"/22.25"	45"/22.25"	N/A	N/A	N/A	47.5"/27"	
4-WHEEL LENGTH/WIDTH⁴,¹⁰	N/A	37.2"/21.54"	40.35"/17.52"	42"/20.5"	40.25"/19.25"	N/A	39.5"/19.5"	43"/21.25"	47"/20.75"	48"/24.25"	45.5"/22.5"	47"/22.25"	N/A	47"/25.5"	57.5"/29.7"	56"/27"-27.8"	47.5"/27"	
FRONT TIRES⁴	8" hub motor	7" solid	7" solid	8" solid	8" solid	8" solid	7" solid	3-Wheel: 8" solid 4-Wheel: 9" solid	3-Wheel: 8" solid 4-Wheel: 9" solid	9" solid	3-Wheel: 8" solid 4-Wheel: 9" solid	10" solid	10" solid	10" solid	13" pneumatic	13" pneumatic	3-Wheel: 10.5" solid 4-Wheel: 10" solid	
REAR TIRES⁴	9" perforated	8" solid	8" solid	8" solid	8" solid	8" solid	8" solid	9" solid	9" solid	10.75" solid	9" solid	10" solid	10" solid	10" solid	13" pneumatic	13" pneumatic	10.75" solid	
SUSPENSION	N/A	N/A	N/A	Front & Rear	N/A	N/A	N/A	N/A	Front & Rear	Front & Rear	N/A	N/A	N/A	Front & Rear	Front & Rear	Front & Rear	N/A	
RANGE PER CHARGE²,³	Up to 9.3 miles (220 lbs.)¹⁵ Up to 7 miles (275 lbs.)¹⁵	Up to 10.1 miles (220 lbs., 11.7Ah) Up to 8.8 miles (300 lbs., 11.7Ah)	Up to 10 miles	Up to 6.6 miles (200 lbs., 8Ah)¹⁵ Up to 13.2 miles (200 lbs., 16Ah)¹⁵	3-Wheel: Up to 8 miles (200 lbs.) Up to 6.2 miles (300 lbs.) 4-Wheel: Up to 8.9 miles (200 lbs.) Up to 6.6 miles (300 lbs.)	Up to 7.3 miles (12Ah, 300 lbs.) Up to 10.5 miles (18Ah, 300 lbs.)	Up to 6.6 miles (12Ah, 300 lbs.) Up to 10.4 miles (18Ah, 300 lbs.)	3-Wheel: Up to 12.75 miles (325 lbs.) 4-Wheel: Up to 12 miles (325 lbs.)	3-Wheel: Up to 17.4 miles (200 lbs.) Up to 22.8 miles (400 lbs.) 4-Wheel: Up to 17.8 miles (200 lbs.) Up to 14.8 miles (400 lbs.)	Up to 21.2 miles (200 lbs.) Up to 15.9 miles (400 lbs.)	3-Wheel: Up to 12.5 miles (300 lbs.) 4-Wheel: Up to 16.3 miles (300 lbs.)	Up to 15.5 miles (400 lbs.)	Up to 13 miles (40Ah, 400 lbs.) Up to 11 miles (U-1, 400 lbs.)	Up to 19.7 miles (200 lbs.) UP to 13 miles (400 lbs.)	Single Battery: Up to 17.4 miles (220 lbs.) Up to 13.8 miles (400 lbs.) Dual Battery: Up to 34.8 miles (220 lbs.) Up to 30.9 miles (350 lbs.)	Up to 23.4 miles (220 lbs.) Up to 14 miles at (500 lbs.)	3-Wheel: Up to 19.8 miles (500 lbs.) 4-Wheel: Up to 19.5 miles (500 lbs.)	
TOTAL WEIGHT WITHOUT BATTERIES¹,⁴,¹¹	60.6 lbs.	41.6 lbs.	35 lbs.	93 lbs.	3-Wheel: 76.3 lbs. 4-Wheel: 81.5 lbs.	75.5 lbs.	83 lbs.	3-Wheel: 92.75 lbs. 4-Wheel: 100 lbs.	3-Wheel: 120 lbs. 4-Wheel: 132 lbs.	181.8 lbs.	3-Wheel: 117 lbs. 4-Wheel: 131 lbs.	3-Wheel: 133 lbs. 4-Wheel: 137 lbs.	130.25 lbs.	160 lbs.	220.6 lbs.	253.8 lbs.	3-Wheel: 183.5 lbs. 4-Wheel: 196.5 lbs.	
WEIGHT OF HEAVIEST PIECE⁴	32.8 lbs. (Rear Base with Seat Back)	N/A (unit does not disassemble)	N/A (unit does not disassemble)	40.75 lbs. (Front Section)	3-Wheel: 29.5 lbs. (Front Section) 4-Wheel: 34.5 lbs. (Front Section)	29.5 lbs. (Front Section)	35 lbs. (Front Section)	3-Wheel: 35.75 lbs. (Front Section) 4-Wheel: 44 lbs. (Front Section)	3-Wheel: 48 lbs. (Rear Section) 4-Wheel: 57 lbs. (Front Section)	75 lbs. (Front Section)	3-Wheel: 48.5 lbs. (Front Section) 4-Wheel: 59 lbs. (Front Section)	3-Wheel: 54 lbs. (Front Section) 4-Wheel: 63 lbs. (Front Section)	56.5 lbs. (Front Section)	77 lbs. (Front Section)	163 lbs. (Base)	201.2 lbs. (Base)	N/A (unit does not disassemble)	
STANDARD SEATING	Compact, Foldable with Removable Back (16" x 16")	Foldable, (15" x 15")	Foldable, 14" x 11.5"	Mesh Back (17" x 17")	Compact, Foldable (17" x 17")	Compact, Foldable (17" x 17")	Mesh Back (17" x 17")	Compact, Foldable (17" x 17")	Compact, Foldable (20" x 18")	Compact, Foldable (18" x 17")	Compact, Foldable (18" x 17")	Compact, Foldable (18" x 17")	Compact, Foldable (18" x 17")	Limited Recline High Back with Sliders & Headrest (18" x 18")	Limited Recline High Back, with Sliders (18" x 18")	Limited Recline High Back (18" x 18")	Depth-Adjustable Limited Recline High Back (22" x 18" x 20")	
BATTERY REQUIREMENTS	25.2V, 9.6Ah (Lithium-ion)	11.7Ah (Lithium-ion) or 16Ah (Lithium-ion)	12Ah (Lithium-ion)	8Ah (Lithium-ion) or 16Ah (Lithium-ion)	(2) 12Ah (SLA/AGM)	(2) 12Ah (SLA/AGM) or (2) 18Ah (SLA/AGM)	(2) 12Ah (SLA/AGM) or (2) 18Ah (SLA/AGM)	(2) 18Ah (SLA/AGM)	(2) U-1 (SLA/AGM)	(2) 40Ah (SLA/AGM)	(2) U-1 (SLA/AGM)	(2) U-1 (SLA/AGM) or (2) 40Ah (SLA/AGM)	(2) U-1 (SLA/AGM) or (2) 40Ah (SLA/AGM)	(2) 50Ah (SLA/AGM)	(1) or (2) 24V, 50Ah Lithium-ion	(2) 12V 75Ah (SLA/AGM)	(2) GEL, 58Ah, NF-22 (Nano-silica)	
BATTERY WEIGHT⁹	4.8 lbs.	8 lbs. (11.7Ah) 9.8 lbs. (16Ah)	4 lbs.	8 lbs. (8Ah) 12 lbs. (16Ah)	9 lbs. ea.⁵	9 lbs. ea.⁵ / 14 lbs. ea.⁷	9 lbs. ea.⁵ / 14 lbs. ea.⁷	14 lbs. ea.⁷	25 lbs. ea.	30.25 lbs. ea	24 lbs. ea	24 lbs. ea./ 32 lbs. ea.	24 lbs. ea./ 32 lbs. ea.	30.2 lbs. ea.	23.8 lbs. ea.	48.8 lbs. ea.	36.5 lbs. ea.	
BATTERY CHARGER	Off-board, 3.5A	Off-board, 3.5A	Off-board, 2A	Off-board, 3.5A	Off-board, 2A	Off-board, 2A	Off-board, 2A	Off-board, 2A	Off-board, 3.5A	Off-board, 3.5A	Off-board, 3.5A	Off-board, 3.5A (U-1/ Off-board, 5A (40Ah)	Off-board, 3.5A (U-1/ Off-board, 5A (40Ah)	Off-board, 5A	Off-board, 8A	Off-board, 8A	Off-board, 5A	
SEAT-TO-DECK HEIGHT RANGE	14.6"	16.61"	16.14"	16.25"-17.75"	16.5"-18"	15.75"-17.25"	16.38"-17.83"	17"-18"	16.5" - 18"	17.75"-19.25"	17.25"-18.75"	17.5"-19"	17.5" - 18"	16.5"-18.5"	17.9"-20.9"	16.54"-18.39"	17.13" - 19.13"; 17.36"-23.61" (with power elevating seat)	
SEAT-TO-GROUND HEIGHT RANGE⁴	22.25"	21.42"	20.87"	21.25"-22.75"	21.5"-23"	18.25"-20.25"	21.38"-22.83"	20.5"-22"	22"-23.5"	23.25"-24.75"	22.5"-24"	23.5"-25"	23.5" - 25"	23" - 25"	28.74"	24.17"-25.59"	22.1" - 24.1" 22.32"-28.57" (with power elevating seat)	
STANDARD BASE COLORS	Sugar Plum (matte) Iceberg Blue (matte) Black (matte) White (matte)	Black Mint	Black Carbon Fiber	Garnet Red (Matte) Euro Grey ¹⁴	Black	Red Blue¹⁴	Cardinal Red, Black, Iceberg Blue, White Robin's Egg Blue, Chell-O, Sugar Plum, Tangerine	Red Blue¹⁴	Grey Street True Blue	Candy Apple Red Ocean Blue (matte) Pearl White	Candy Apple Red	Candy Apple Red Viper Blue	Candy Apple Red	Candy Apple Red True Blue	Black Grey (matte) Orange	Satin Aluminum Dark Violet, Red Peacock Blue, Sunflower, White Camo, Black Camo	Garnet Red (matte) Ocean Blue (matte)	
HCPCS CODE	Non-Medical*	Non-Medical*	K0899**	K0899**	K0800 - Group 1 Standard	K0800 - Group 1 Standard	K0800 - Group 1 Standard	K0801 - Group 1 Heavy Duty	K0899**	K0899**	K0800 - Group 1 standard	K0801 - Group 1 Heavy Duty	K0801- Group 1 Heavy Duty	K0899**	K0899**	K0899**	K0802- Group 1 Very Heavy Duty (3-Wheel) K0808- Group 2 Very Heavy Duty (4-Wheel)¹²	
WARRANTY	3-yr limited on Frame 2-yr limited on Electronics 1-yr limited on Drivetrain, Motor, Battery, & Charger	3-yr limited on Frame 2-yr limited on Electronics 1-yr limited on Drivetrain, Motor, Battery, and Charger	Limited Lifetime on Frame 1-yr limited on Electronics, Drivetrain and Battery	Lifetime limited on Frame 2-yr limited on Electronics 2-yr limited on Drivetrain 2-yr limited on Battery	Lifetime limited on Frame 2-yr limited on Electronics 2-yr limited on Drivetrain 6-month limited on Batteries	Lifetime limited on Frame 2-yr limited on Electronics 2-yr limited on Drivetrain 6-month limited on Batteries	Lifetime limited on Frame 2-yr limited on Electronics 2-yr limited on Drivetrain 6-month limited on Batteries	Lifetime limited on Frame 2-yr limited on Electronics 2-yr limited on Drivetrain 6-month limited on Batteries	Lifetime Limited on Frame 2-yr limited on Electronics 2-yr limited on Drivetrain 6-month limited on Batteries	Lifetime Limited on Frame 2-yr limited on Electronics 2-yr limited on Drivetrain 6-month limited on Batteries	Lifetime Limited on Frame 2-yr limited on Electronics 2-yr limited on Drivetrain 6-month limited on Batteries	Lifetime Limited on Frame 2-yr limited on Electronics 2-yr limited on Drivetrain 6-month limited on Batteries	Lifetime limited on Frame 2-yr limited on Electronics 2-yr limited on Drivetrain 6-month limited on Batteries	Lifetime limited on Frame 2-yr limited on Electronics 2-yr limited on Drivetrain 6-month limited on Batteries	Lifetime limited on Frame 2-yr limited on Electronics 2-yr limited on Drivetrain 6-month limited on Battery	Lifetime limited on Frame 2-yr limited on Electronics 2-yr limited on Drivetrain 6-month limited on Battery	Lifetime limited on Frame 2-yr limited on Electronics 2-yr limited on Drivetrain 6-month limited on Battery	Lifetime limited on Frame 2-yr limited on Electronics 2-yr limited on Drivetrain 6-month limited on Batteries

TRAVEL MOBILITY

MID-SIZE LUXURY

FULL-SIZE

HEAVY-DUTY

*Non-Medical Devices: This product is not a medical device and is not intended to assist, treat, diagnose, or alleviate any medical condition or disability.
**K0899: Power mobility device, not coded by DME PDAC.

DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage, and bundling guidelines may apply.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed & range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%.