

**Footnotes**

- All shroud panels included.
- Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. This specification can be subject to a variance of (+ or -) 10%.
- Battery weight not included.
- Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3%.
- Battery pack is 18 lbs. with 12AH batteries.
- Deep cycle (AGM or gel cell type recommended).
- 18 AH battery pack weighs 28 lbs.

- Tested in accordance with ANSI/RESNA, WC Vol 2, Section 4 and ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- Battery weight may vary based on manufacturer.
- Width of base (not including seat).
- Does not account for anti-tip clearance.
- Includes standard seat weight.
- Must be billed to Medicare as either a beneficiary requested upgrade with a properly executed ABN, or as an upgrade at no additional charge to the beneficiary.
- At midframe.
- Light package.

# PRIDE® MOBILITY SCOOTERS

## LIVE YOUR BEST®

# END-TO-END



	Go-Chair	Go-Go® ES 2	Go-Go® Ultra X	Go-Go Elite Traveller®	Go-Go Elite Traveller® Plus	Go-Go® Sport	Go-Go® LX w/ CTS Suspension	Go-Go® Folding Scooter	Victory® 9	Victory® 10	Victory 10® LX w/ CTS Suspension <sup>15</sup>	Victory® Sport <sup>15</sup>	Maxima	Pursuit® <sup>15</sup>	Pursuit® XL <sup>15</sup>	
<b>FRONT TIRES</b>	6" SOLID (CASTER)	2" x 7.5" SOLID	2" x 8" SOLID (3W) 2" x 7" SOLID (4W)	2" x 8" SOLID (3W) 2" x 7" SOLID (4W)	2" x 8" SOLID	3" x 8" SOLID (3W) 3" x 9" SOLID (4W)	2" x 8" SOLID (3W) 2" x 7" SOLID (4W)	1.6" x 7" SOLID	3W: 3" x 8" SOLID 4W: 3" x 9" SOLID	3.5" x 10" SOLID	3.5" x 10" SOLID	3.5" x 10" SOLID	3" x 10" SOLID	3" x 13" SOLID	13" PNEUMATIC	
<b>REAR TIRES</b>	9" SOLID (DRIVE)	3" x 7.5" SOLID	2.5" x 8" SOLID	2.5" x 8" SOLID	3" x 9" SOLID	3" x 9" SOLID	2.5" x 8" SOLID	2" x 8" SOLID	3" x 9" SOLID	3.5" x 10" SOLID	3.5" x 10" SOLID	4" x 10" SOLID	4" x 10" SOLID	13" SOLID	13" PNEUMATIC	
<b>MAXIMUM SPEED<sup>2</sup></b>	UP TO 3.7 MPH	UP TO 4 MPH	UP TO 4 MPH	UP TO 4 MPH	UP TO 4.5 MPH	UP TO 4.7 MPH	UP TO 4.4 MPH	UP TO 4 MPH	UP TO 5.25 MPH	UP TO 5.25 MPH	UP TO 5.6 MPH	UP TO 8 MPH	UP TO 5.25 MPH	UP TO 8 MPH	UP TO 8.5 MPH	
<b>GROUND CLEARANCE<sup>4,11,14</sup></b>	1.6" at anti-tip bracket 1.9" at frame	1.5" at motor	2.25" at mid-frame 1.25" at motor	2"	3" at frame 1.5" at motor	3" at frame 1.5" at motor	1.5" at motor	2.3"	3W: 2" 4W: 2.25"	3" at mid-frame	2.2" at motor	2.25"	4.25" 3W 3.25" 4W	3.56"	5.25"	
<b>3-W TURNING RADIUS<sup>4</sup></b>	N/A	31"	33"	33"	34.5"	36"	33"	32"	40.25"	42"	N/A	N/A	50.38"	N/A	N/A	
<b>4-W TURNING RADIUS<sup>4</sup></b>	25.4" with foot platform	N/A	44"	44"	45.72"	50.5"	44.75"	N/A	51.5"	54"	59.9"	51.75"	54"	67.75" LEFT 66.75" RIGHT	82.5"	
<b>3-W LENGTH/WIDTH<sup>4,10</sup></b>	N/A	36.75"/19.5"	37"/19.5"	37"/19.5"	40"/21.25"	42"/21.25"	39.75"/20.5"	37"/19"	43.5"/22.25"	45.25"/22.25"	N/A	N/A	46"/25.5"	N/A	N/A	
<b>4-W LENGTH/WIDTH<sup>4,10</sup></b>	28.5" without footplate 33"/22.6"	N/A	39.5"/19.5"	39.5"/19.5"	40.25"/21.25"	43"/21.25"	41"/20.5"	N/A	45.5"/22.5"	47"/22.25"	47"/22.2"	47"/26"	46.9"/25.5"	53.25"/25.25"	56"/27.25"	
<b>SEAT-TO-GROUND HEIGHT RANGE</b>	22"-23"	19"-20"	18.13"-20.13"	18.13"-20.13"	20.5"-22"	20.5"-22"	20.75"-22.25"	20.5"	19.25"-20.75"	20.5"-22"	23.5"-25.5"	21.5"-24.25"	22.5"-24.5"	24.75"-26.25"	26.25"-27.75"	
<b>SEAT TYPE (STANDARD)</b>	COMPACT, FOLDABLE (18"X17")	COMPACT, FOLDABLE (17"X15")	COMPACT, FOLDABLE (17"X17")	COMPACT, FOLDABLE (17"X17")	COMPACT, FOLDABLE (17"X17")	COMPACT, FOLDABLE (17"X17")	COMPACT, FOLDABLE (17"X17")	COMPACT, FOLDABLE (16"X13")	FOLDABLE, MOLDED PLASTIC (18"X17")	FOLDABLE, MOLDED PLASTIC (18"X17")	RECLINING H/B W/SLIDERS & HEADREST (18"X18")	RECLINING H/B W/SLIDERS & HEADREST (18"X18")	DELUXE CONTOURED MEDIUM-BACK (22"X18")	RECLINING H/B W/SLIDERS & HEADREST (18"X18")	HIGH-BACK PILLOW TOP W/SLIDERS & HEADREST (19.5"X16")	
<b>(3/4-WHEEL) WEIGHT<sup>3,4,12</sup></b>	67 LBS.	58.4 LBS.	76 LBS./80 LBS.	76 LBS./80 LBS.	86.5 LBS./87 LBS.	92.8 LBS./92.2 LBS.	85.4 LBS./87.4 LBS.	51.4 LBS.	117.5 LBS./131 LBS.	133.5 LBS./137.5 LBS.	179 LBS.	173 LBS.	171 LBS./176 LBS.	225 LBS.	253.9 LBS.	
<b>SEAT-TO-DECK HEIGHT RANGE</b>	16.25"-17.25"	16.25"-17.25"	15.75"-17.25"	15.75"-17.25"	16.25"-17.75"	17"-18"	16.25"-17.25"	13"	17.25"-18.75"	17.5"-19"	16.5"-18.5"	16"-18.5"	3W: 16"-18" 4W: 16"-18"	17"-18.75"	21"-22.5"	
<b>BATTERY WEIGHT<sup>9</sup></b>	14 LBS. EA.	9 LBS. EA.	9 LBS. EA. <sup>6</sup>	9 LBS. EA. <sup>6</sup> / 14 LBS. EA. <sup>7</sup>	14 LBS. EA. <sup>7</sup>	14 LBS. EA. <sup>7</sup>	9 LBS. EA. <sup>6</sup> / 14 LBS. EA. <sup>7</sup>	9 LBS. EA. (SEALED LEAD ACID)/5 LBS. EA. (LITHIUM-ION)	24.5 LBS. EA.	24.5 LBS. EA./ 32.5 LBS. EA.	32.5 LBS. EA.	32.5 LBS. EA.	38 LBS. EA. (NF-22 OR 55 AH) 24.5 LBS. EA. (U-1)	38 LBS. EA. (NF-22 OR 55 AH) 53 LBS. EA. (75 AH)	53 LBS. EA. (75 AH) 60 LBS. EA. (100 AH)	
<b>BATTERY REQUIREMENTS<sup>6</sup></b>	(2)18AH	(2)12AH	(2)12AH	(2)12AH OR (2)18AH	(2)18AH	(2)18AH	(2)12AH OR (2)18AH	(2)12AH OR (1) 11.5 AH (LITHIUM-ION)	(2) U1	(2) U1 OR (2) 40AH	(2) 40AH	(2) 40AH	(2) NF-22 (2) U-1 (POWER ELEVATING SEAT)	(2) NF-22 OR (2) GROUP 24	(2) GROUP 27 OR (2) GROUP 24	
<b>PER CHARGE RANGE (MAX.)<sup>2,8</sup></b>	UP TO 8.74 MILES	UP TO 6.8 MILES	3W: UP TO 6.9 MILES 4W: UP TO 7.2 MILES	3W: UP TO 7.2 MILES (12AH) UP TO 10.4 MILES (18AH) 4W: UP TO 6 MILES (12AH) UP TO 9.7 MILES (18AH)	3W: UP TO 14.4 MILES 4W: UP TO 8.3 MILES	3W: UP TO 12.7 MILES 4W: UP TO 12 MILES	3W: UP TO 5.9 MILES (12AH) UP TO 9.44 MILES (18AH) 4W: UP TO 5.5 MILES (12AH) UP TO 10.5 MILES (18AH)	4W: UP TO 12.5 MILES (LITHIUM-ION) 4W: UP TO 9.3 MILES (SEALED LEAD ACID)	UP TO 13 MILES	UP TO 15.5 MILES	UP TO 16.16 MILES	UP TO 11 MILES	3W: UP TO 18 MILES, U-1 (POWER ELEVATING SEAT) UP TO 9.6 MILES 4W: UP TO 16 MILES	G24: UP TO 19 MILES NF-22: UP TO 14 MILES	G27: UP TO 29.51 MILES G24: UP TO 22.13 MILES	
<b>BATTERY CHARGER</b>	OFF-BOARD, 2A	OFF-BOARD, 1.5A	OFF-BOARD, 2A	OFF-BOARD, 2A	OFF-BOARD, 2A	OFF-BOARD, 2A	OFF-BOARD, 2A	OFF-BOARD, 2A (SEALED LEAD ACID) OFF-BOARD, 2.5A (LITHIUM-ION)	OFF-BOARD, 3.5A (U-1)	OFF-BOARD, 3.5A (U-1) OFF-BOARD, 5A (40 AH)	OFF-BOARD, 5A (40 AH)	OFF-BOARD, 5A	OFF-BOARD, 5A	OFF-BOARD, 5A (NF-22) OFF-BOARD, 8A (GROUP 24)	OFF-BOARD, 8A	
<b>SUSPENSION</b>	LIMITED	N/A	N/A	N/A	N/A	N/A	FRONT & REAR	N/A	N/A	N/A	FRONT & REAR	FRONT & REAR	3-W REAR/ 4-W FRONT & REAR	FRONT & REAR	FRONT & REAR	
<b>MAXIMUM WEIGHT CAPACITY</b>	300 LBS.	250 LBS.	260 LBS.	300 LBS.	300 LBS.	325 LBS.	300 LBS.	250 LBS.	300 LBS.	400 LBS.	400 LBS.	350 LBS.	500 LBS.	400 LBS.	400 LBS.	
<b>STANDARD BODY COLORS</b>	RUBY RED, SAPHIRE BLUE, CITRINE YELLOW, AMBER ORANGE, PEARL WHITE, ROSE QUARTZ	BLACK	RED	RED, BLUE <sup>1</sup>	RED, BLUE <sup>1</sup>	RED, BLUE <sup>1</sup>	RED, BLUE <sup>1</sup>	WHITE	CANDY APPLE RED, VIPER BLUE	CANDY APPLE RED, VIPER BLUE	CANDY APPLE RED, TRUE BLUE	CANDY APPLE RED, VIPER BLUE	CANDY APPLE RED, VIPER BLUE	CANDY APPLE RED, SILVER	CANDY APPLE RED, SILVER	
<b>WEIGHT OF HEAVIEST PIECE<sup>4</sup></b>	36 LBS. (REAR SECTION)	41.4 LBS. (3-WHEEL) (SCOOTER BASE)	29.5 LBS. (3-WHEEL) (FRONT SECTION) 34 LBS. (4-WHEEL) (FRONT SECTION)	29.5 LBS. (3-WHEEL) (FRONT SECTION) 34 LBS. (4-WHEEL) (FRONT SECTION)	32.5 LBS. (3-WHEEL) (REAR SECTION) 34 LBS. (4-WHEEL) (FRONT SECTION)	35 LBS. (3-WHEEL) (FRONT SECTION) 44 LBS. (4-WHEEL) (FRONT SECTION)	34 LBS. (3-WHEEL) (FRONT SECTION) 36.6 LBS. (4-WHEEL) (FRONT SECTION)	51.4 LBS. (UNIT DOES NOT DISASSEMBLE - UNIT WEIGHT WITHOUT BATTERY PACK)	48.5 LBS. (3-WHEEL) (FRONT SECTION) 59 LBS. (4-WHEEL) (FRONT SECTION)	54.5 LBS. (3-WHEEL) (FRONT SECTION) 61 LBS. (4-WHEEL) (FRONT SECTION)	77.6 LBS. (4-WHEEL) (FRONT SECTION)	68 LBS. (4-WHEEL) (FRONT SECTION)	76 LBS. (3-WHEEL) (REAR SECTION) 76 LBS. (4-WHEEL) (REAR SECTION)	95 LBS. (FRONT SECTION)	205.5 LBS. (MAIN FRAME)	
<b>HCPCS</b>	K0899	K0800 - GROUP 1 STANDARD (3W)	K0899	K0800 - GROUP 1 STANDARD	K0806 - GROUP 2 STANDARD <sup>13</sup>	K0801 - GROUP 1 HEAVY DUTY	K0899	K0899	K0800 - GROUP 1 STANDARD	K0801 - GROUP 1 HEAVY DUTY	K0899	K0899	K0802 - GROUP 1 VERY HEAVY DUTY (3W) K0808 - GROUP 2 VERY HEAVY DUTY (4W)	K0899	K0899	
<b>WARRANTY</b>	5-YEAR LIMITED ON FRAME; 13-MONTH LIMITED ON ELECTRONICS; 13-MONTH LIMITED ON DRIVE MOTORS	LIFETIME LIMITED ON FRAME; 2-YEAR LIMITED ON ELECTRONICS; 2-YEAR LIMITED ON DRIVETRAIN	LIFETIME LIMITED ON FRAME; 2-YEAR LIMITED ON ELECTRONICS; 2-YEAR LIMITED ON DRIVETRAIN	LIFETIME LIMITED ON FRAME; 2-YEAR LIMITED ON ELECTRONICS; 2-YEAR LIMITED ON DRIVETRAIN	LIFETIME LIMITED ON FRAME; 2-YEAR LIMITED ON ELECTRONICS; 2-YEAR LIMITED ON DRIVETRAIN	LIFETIME LIMITED ON FRAME; 2-YEAR LIMITED ON ELECTRONICS; 2-YEAR LIMITED ON DRIVETRAIN	LIFETIME LIMITED ON FRAME; 2-YEAR LIMITED ON ELECTRONICS; 2-YEAR LIMITED ON DRIVETRAIN	2-YEAR LIMITED ON FRAME; 18-MONTH LIMITED ON CONTROLLERS; 1-YEAR LIMITED ON ELEC. COMPONENTS	LIFETIME LIMITED ON FRAME; 2-YEAR LIMITED ON ELECTRONICS; 2-YEAR LIMITED ON DRIVETRAIN	LIFETIME LIMITED ON FRAME; 2-YEAR LIMITED ON ELECTRONICS; 2-YEAR LIMITED ON DRIVETRAIN	LIFETIME LIMITED ON FRAME; 2-YEAR LIMITED ON ELECTRONICS; 2-YEAR LIMITED ON DRIVETRAIN	LIFETIME LIMITED ON FRAME; 2-YEAR LIMITED ON ELECTRONICS; 2-YEAR LIMITED ON DRIVETRAIN	LIFETIME LIMITED ON FRAME; 2-YEAR LIMITED ON ELECTRONICS; 2-YEAR LIMITED ON DRIVETRAIN	LIFETIME LIMITED ON FRAME; 2-YEAR LIMITED ON ELECTRONICS; 2-YEAR LIMITED ON DRIVETRAIN	LIFETIME LIMITED ON FRAME; 2-YEAR LIMITED ON ELECTRONICS; 2-YEAR LIMITED ON DRIVETRAIN	LIFETIME LIMITED ON FRAME; 2-YEAR LIMITED ON ELECTRONICS; 2-YEAR LIMITED ON DRIVETRAIN
	<b>TRAVEL MOBILITY</b>								<b>MID-SIZE LUXURY</b>				<b>HEAVY-DUTY</b>	<b>FULL-SIZE LUXURY</b>		

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed & range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement.