

**New Improved Motor
&
Enhanced Durability**




HERCULES[®]

3000

Pride Mobility Products[®] Corporation

The Hercules[®] 3000 Lift is the only lift specifically engineered for Pride products, and like the products it is designed to complement, it is strong, durable and dependable. The Hercules 3000 Lift fits readily into most vehicles and makes transporting your Pride[®] Scooter or Jazzy[®] Power Chair quick and convenient.

HERCULES 3000
Designed to lift most of our Pride Jazzy Power Chairs, Pride Scooters and PMVs.*



Various docking devices are available for most Pride Scooters and Jazzy Power Chairs. Attachments are available that will allow you to leave the seat on your Pride Scooter or Jazzy Power Chair depending on the cargo space of your vehicle.*

* Please contact your Pride Provider to confirm models/seating systems that the Hercules is capable of lifting and also to determine if the Hercules can be installed in your vehicle. (The lifting capacity of the Hercules is 275 lbs.)



FEATURES

- Easily installs in most vehicles
- Conveniently adjustable with 2 quick-release pins
- Complete line of docking devices available for Pride products that weigh up to 275 lbs.
- Easily assembles and disassembles
- Our exclusive new pivoting hook for easy maneuvering of Scooter/Jazzy into vehicle
- Transport lock pin added for boom arm stability during transport
- Improved motor gear box design offers smoother and more efficient operation than ever before
- 1-year limited warranty
- Optional corded remote available for ease of operation

800-800-8586 (U.S.) • 888-570-1113 (Canada) • www.pridemobility.com

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

Pride[®]
Mobility Products Corp.