# JUNIT OF L

## MORE ADVANCED FEATURES & EXCELLENT PERFORMANCE

# Jazzy® 1107

The Jazzy® 1107 brings the most advanced achievements of our engineering efforts together in one compact yet superbly balanced package. Direct connect components eliminate wiring harnesses and provide clean style and unquestionable convenience.

> Jazzy® 1107 shown in Candy Apple Red with gray contoured mediumback seat, foot platform and Pride Flight controller.



800-800-8586 (US) · 888-570-1113 (Canada) · www.pridemobility.com



The Jazzy® 1107 series models include: • Jazzy 1107-2S-C • Jazzy 1107-2S-SS

# Jazzy® 1107

Weight capacity:	300 lbs.
Maximum speed:	Up to 4.5 mph
Ground clearance:	2.5"
Turning radius:	21.875"
Overall length:	36.375"
Overall width:	23.375"
Wheels: Front anti-tips: Rear casters:	Drive: 10" 5" solid 6" solid
Suspension:	Advanced Active-Trac® with independent rear casters
Drive train:	Two-motor, mid-wheel drive
Braking system:	Intelligent braking (electronic, regenerative disc brakes)
Std. Electronics:	60 amp, Pride Flight controller
Specialty controls:	N/A
Battery charger:	Off-board, 5 amp (standard)
Per-charge range:	Up to 15 miles
Battery requirements:	(2) 12 volt, deep cycle
Battery size:	U-1
Transport weight Battery weight: Base weight: Std. seat weight:	24.5 lbs. each 105 lbs. 37 lbs. (medium-back)
Std. Seat dimensions:	18" × 18"
Warranty:	Lifetime limited warranty on frame; 2-year limited warranty on electronics; 18-months limited warranty on drive motors

### FEATURES

- Aluminum frame rails and center locking mechanism provide a very rigid platform for ultimate stability
- Advanced Active-Trac<sup>®</sup> Suspension for a smoother ride over challenging terrain
- Innovative battery boxes and front module eliminate battery and motor harnesses

### OPTIONS

- Contoured high-back seat
- Contoured medium/high-back with solid seat pan
- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Cane/crutch holder
- Walker holder
- Solid tires
- Swing-away joystick mount
- Safety flag

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure.



- Backpacker<sup>™</sup>
- · Power chair lift with straps
- Full Platform





800-800-8586 (US) • 888-57<u>0-1113 (Canada) • www.pridemobility.com</u>

Viper

Blue

Onyx

Black

Candy

Apple Red

Colors: