

# Jazzy®

## 600 SERIES

The Jazzy® 600 provides excellent performance and maximum stability. Its compact design and standard ATX Suspension allows it to excel indoors and out making it a superb all-around power chair choice.



Jazzy® 600  
shown in *Candy Apple Red*


(US) 800-800-8586 ▪ (Canada) 888-570-1113  
Visit us on the web: [www.pridemobility.com](http://www.pridemobility.com)

**Pride®**  
Mobility Products Corp.

## SPECIFICATIONS

<b>Weight capacity</b>	300 lbs.
<b>Maximum speed</b>	Up to 5 mph
<b>Ground clearance</b>	3"
<b>Turning radius</b>	20"
<b>Overall length</b>	36.5" without front riggings
<b>Overall width</b>	23"
<b>Wheels</b>	Drive: 14"
<b>Front casters</b>	5" solid
<b>Rear casters</b>	6" solid
<b>Suspension</b>	Active-Trac® ATX Suspension
<b>Drive train</b>	Two-motor, in-line, Mid-Wheel 6®
<b>Braking system</b>	Intelligent braking (electronic, regenerative disc brakes)
<b>Standard electronics</b>	60 amp, PG VR2
<b>Specialty controls</b>	Yes
<b>Battery charger</b>	On-board, 4 amp (standard)
<b>Per-charge range</b>	Up to 20 miles
<b>Battery requirements</b>	(2) 12 volt, deep cycle
<b>Battery size</b>	NF-22
<b>Battery weight</b>	38 lbs. each
<b>Base weight</b>	127 lbs.
<b>Standard seat weight</b>	37 lbs. (medium-back)
<b>Max. seat dimensions</b>	22" x 20"
<b>Seat-to-floor range</b> (high-back seat)	21" - 24"
<b>Warranty</b>	Lifetime limited warranty on frame; 2-year limited warranty on electronics; 18-month limited warranty on drive motors

## COLOR CHOICES



## MODELS

Jazzy 600 2S-C  
 Jazzy 600 2S-SS  
 Jazzy 600 2SP-C (250 lbs. weight capacity)  
 Jazzy 600 2SP-SS (250 lbs. weight capacity)  
 Jazzy 600 3S-C  
 Jazzy 600 3S-SS  
 Jazzy 600 3SP-C (250 lbs. weight capacity)  
 Jazzy 600 3SP-SS (250 lbs. weight capacity)

## FEATURES

- OMNI-Casters (nylon, spherical-shaped casters) on front and rear to prevent wheel hang-ups
- Side-mounted, easily accessible freewheel levers
- Active-Trac® ATX Suspension incorporates front OMNI-Casters for enhanced performance over more varied terrain
- Easy front access to batteries
- Easily accessible side charger port

## OPTIONS

- Solid seat pan
- Cup holder
- Swing-away joystick
- Swing-away footrests
- TRU-Balance® 2 Power Tilt
- Elevating leg rests
- Oxygen tank holder
- Cane/crutch holder
- Walker holder
- Synergy® Seat

## RECOMMENDED LIFTS

- Silver Star® Backpacker® Plus
- Silver Star® Backpacker® MV
- Silver Star® Full Platform
- Silver Star® Silver Boom 300

**Pride**®  
 Mobility Products Corp.

(US) 800-800-8586 • (Canada) 888-570-1113  
 Visit us on the web: [www.pridemobility.com](http://www.pridemobility.com)

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%.

\* Special order seats or non-standard dimensions will not be fold-flat and depth adjustable.