

# Jazzy®

## 614 HD

The Jazzy® 614 HD features patented Mid-Wheel 6® technology and Active-Trac ATX Suspension for an outstanding level of performance, maneuverability and tight turning radius, all in one heavy-duty package. Independent front casters raise and lower as needed to compensate instantly for changes in terrain, making transitions easy. The Jazzy 614 HD features dual enhanced in-line motors for greater efficiency, torque, range and performance. Independent drive wheel suspension allows each drive wheel to be independently supported by a dynamic spring-strut dampening system, reducing the vibrations caused by driving surfaces. All-around performance indoors and out over most surfaces, inclines and slopes, makes the Jazzy 614 HD an excellent heavy-duty power chair choice.

*Jazzy® 614 HD  
shown in Candy Apple Red  
with high-back seat*



**Pride**®  
Mobility Products Corp.

(US) 800-800-8586 • (CAN) 888-570-1113 • [pridemobility.com](http://pridemobility.com)



## SPECIFICATIONS

Weight capacity	450 lbs.
Maximum speed	Up to 4.3 mph
Ground clearance	2.25"
Turning radius	22.75"
Overall length	35.5"
Overall width	25.125"
Drive wheels	14" solid
Front wheels	6" solid casters
Rear wheels	6" solid casters
Suspension	Active-Trac® ATX
Drivetrain	Two-motor, in-line, Mid-Wheel 6®
Braking system	Regenerative and electro-mechanical
Standard electronics	5A, Dynamic Shark controller
Battery charger	Off-board, 5 amp
Per-charge range	Up to 11 miles
Battery requirements	(2) 12 volt, deep cycle, NF-22
Battery weight	37 lbs. each
Base weight	121.5 lbs.
Standard seat weight	51.5 lbs.
Standard seat	20" x 18-20"
Max. seat dimensions	24" x 24"
Warranty	Lifetime limited on frame; 13-month limited on electronics; 13-month limited on drive motors

## COLOR CHOICES



Candy  
Apple Red



Viper Blue

## OPTIONAL SEATING

- Synergy® Seat
- Solid Seat Pan

## MODELS

Jazzy 614HD 2HD-SS	Code: K0824
Jazzy 614HD 2HD-C	Code: K0825

## FEATURES

- Depth adjustable, black high-back seat
- Universal, fixed, in-line joystick mount
- Flip-up height adjustable foot platform
- Lap belt (50", 70", 80", 90")

## OPTIONS & ACCESSORIES

- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Walker holder
- Cane/crutch holder
- Saddle bag
- Swing-away footrest
- Elevating leg rests
- Heel loops
- Military patches
- Swing-away joystick

## RECOMMENDED LIFTS

- Backpacker® Plus
- Backpacker® AVP
- Outlander Full Platform
- Olympian
- Commander



Use the Pride Lift Finder to find a compatible lift for your mobility product:

[prideliftsandramps.com](http://prideliftsandramps.com)

**Pride**®  
Mobility Products Corp.

(US) 800-800-8586 • (CAN) 888-570-1113 • [pridemobility.com](http://pridemobility.com)



The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.