



The future of mobility is here.



# Jazzy<sup>®</sup> air

## EMPOWERS

The Jazzy Air<sup>®</sup> elevates 10" in 16 seconds, providing quick environmental access to countertops, table tops, cabinetry and shelving. Users can experience face-to-face social engagement while safely driving up to 3.5 mph while elevated.

- Elevate up to 10 inches in 16 seconds
- Enjoy face-to-face social engagement
- Safely drive elevated at walking speed (3.5 mph)



[pridemobility.com](http://pridemobility.com)

# Specifications

<b>Weight Capacity</b>	300 lbs.
<b>Maximum Speed<sup>1</sup></b>	Up to 4.0 mph (non-elevated position) Up to 3.5 mph (elevated position)
<b>Ground Clearance<sup>2</sup></b>	2.42" at footplate
<b>High-Back Seat to Floor Range</b>	23" - 33"
<b>Turning Radius<sup>2</sup></b>	22.5"
<b>Overall Size<sup>2,6</sup></b>	Length: 42.23" with foot rigging; 38" w/out foot rigging Width: 23.375" at the wheels 25.5" with armrests
<b>Drive Controller</b>	40 amp, PG VR2
<b>Drive Wheels</b>	10"
<b>Caster Wheels</b>	6" solid caster
<b>Suspension</b>	Active-Trac <sup>®</sup> with Mid-Wheel 6 <sup>®</sup> Technology
<b>Range<sup>1,4</sup></b>	Up to 15.29 miles w/ U1 batteries Up to 18.4 miles w/ 40 Ah batteries
<b>Component Weights<sup>5</sup></b>	Base: 225 lbs. w/ U1 batteries 245 lbs. w/ 40 Ah batteries U1 Batteries: 25 lbs. 40 Ah Batteries: 32.2 lbs. (optional) 181 lbs. (base & seat)
<b>Drivetrain</b>	Two motor, Mid-Wheel 6 <sup>®</sup> drive
<b>Brakes</b>	"Intelligent Braking" (electronic regenerative, disc park brakes)
<b>Batteries<sup>3</sup></b>	(2) 12 V, U1 (2) 40 Ah
<b>Battery Charger</b>	3.5 amp, off-board (U1); 5 amp, off-board (40 Ah)
<b>Max. Allowable Incline</b>	6°
<b>Max. Angle Capacity</b>	6°
<b>Maximum Obstacle Climbing Ability</b>	Lowered Position: 2.2" Elevated Position: 1.8"
<b>Warranty</b>	5 year limited (frame) 13 month limited (batteries) 13 month limited (electronics) 13 month (drive motors)

- 1) Varies with base model, user weight, terrain type, battery amp hour rating (Ah), battery charge, battery condition, motors, controller type, tire type, and tire condition. This specification is subject to a variance of +10%, -5%.
- 2) Due to manufacturing tolerances and continued product improvement, this specification can be subject to a variance of (+ or -) 3%.
- 3) AGM or Gel-Cell type required. See VI. "Batteries and Charging."
- 4) Tested in accordance with ANSI/RESNA, WC Vol. 2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- 5) Battery weight may vary based on manufacturer.
- 6) Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

NOTE: This product conforms to all applicable ANSI/RESNA, ISO 7176 Series, and EN12184 standards. All specifications subject to change without notice.



Find a compatible lift for your vehicle online:  
**PRIDEMOBILITY.COM/LIFT-FINDER**

# Features

- Achieve 10" of power adjustable seat height in 16 seconds
- Elevate or lower seat height while on the move with the "Air" button or PG VR2 joystick drive controller
- Seat and integrated footplate design provides maximum comfort in both elevated and non-elevated positions
- Angle adjustable footplate with three height adjustments
- Choice of 16" x 16"-18", 18" x 18"-20" or 20" x 20"-22" depth adjustable, fold-flat, high-back reclining sport seat or solid seat pan
- Angle, height and depth adjustable armrest
- Bright LED marker lights
- 60" lap belt

# Accessories

- Rear basket
- Weather cover
- Cup holder
- Flag kit
- Cell phone holder
- Cane holder
- Walker holder
- Oxygen holder
- Saddle bag
- Replacement shroud (ten available colors)
- Alternative size lap belt: 50", 70", 80", 90"
- XLR USB Charger

# Options

- 40 amp AGM batteries
- Swing-away joystick
- Stealth<sup>®</sup> cushions

# Colors



Authorized Pride<sup>®</sup> Provider

(US) 800-800-8586 • (CAN) 888-570-1113 • pridemobility.com