

# Jazzy®

## SELECT ELITE

The Jazzy Select® Elite delivers high performance and maneuverability in a compact, stylish design for your active lifestyle. In-line, front-wheel drive technology gives the Jazzy Select Elite excellent stability, provides a tight turning radius around corners for optimal indoor maneuverability, climbs obstacles and grade changes with ease and handles various terrains indoors and out for an overall solid performance. The Jazzy Select Elite comes standard with Pride's exclusive Jazzy® Armor® for today's rental environment, depth-adjustable high-back seat, and flip-up height adjustable foot platform.



*Jazzy Select® Elite  
shown in Jazzy Red with high-back  
seat and Jazzy® Armor® (standard)*

**Pride®**  
Mobility Products Corp.

(US) 800-800-8586 • (CAN) 888-570-1113 • [pridemobility.com](http://pridemobility.com)



## SPECIFICATIONS

Weight Capacity	300 lbs.
Maximum Speed	Up to 4 mph
Ground Clearance	1.5"
Turning Radius	24.75"
Overall Length	34.75"
Overall Width	22.75"
Drive Wheels	9" solid
Front Wheels	3" solid anti-tips
Rear Wheels	6" solid casters
Suspension	Limited
Drivetrain	Two-motor, in-line, front-wheel drive
Braking System	Regenerative and electro-mechanical
Standard Electronics	40 amp, PG GC 3 controller
Battery Charger	Off-board, 3 amp
Per-Charge Range	Up to 15 miles
Battery Requirements	(2) 12 volt, deep cycle, U-1
Battery Weight	24.5 lbs. each
Base Weight	68.5 lbs.
Standard Seat Weight	43 lbs. (comfort, high-back)
Standard Seat	18" x 18-20"
Max. Seat Dimensions	20" x 20"
Warranty	5-year limited on frame; 13-month limited on electronics; 13-month limited on drive motors

## COLOR CHOICES



Jazzy  
Red



Jazzy  
Blue

## OPTIONAL SEATING

- Synergy® Seat
- Solid seat pan

## MODELS

Jazzy Select Elite-2S-SS	Code: K0822
Jazzy Select Elite-2S-C	Code: K0823

## FEATURES

- Color-through shroud
- Jazzy® Armor® includes shroud and controller guards to protect from daily wear and tear
- Depth adjustable, black high-back seat
- Universal, fixed, in-line joystick mount
- Flip-up height adjustable foot platform
- Black, flat-free, non-marking tires
- 60" lap belt

## OPTIONS & ACCESSORIES

- Cup holder
- Weather cover
- Rear basket
- Oxygen tank holder
- Walker holder
- Cane/crutch holder
- Saddle bag
- Lap belt (50", 70", 80", 90")
- Swing-away footrest
- Elevating leg rests
- Heel loops
- Military patches
- Swing-away joystick

## RECOMMENDED LIFTS

- Backpacker® Plus
- Backpacker® AVP
- Outlander Full Platform
- Olympian
- Commander

**LIFT FINDER**

Use the Pride Lift Finder to find a compatible lift for your mobility product:

[prideliftsandramps.com](http://prideliftsandramps.com)

**Pride**<sup>®</sup>  
Mobility Products Corp.

(US) 800-800-8586 • (CAN) 888-570-1113 • [pridemobility.com](http://pridemobility.com)



The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.