

REHAB POWER CHAIR SOLUTIONS



RIVAL®
shown in *Orange Crush*
with TRU-Balance® 3
Power Positioning

RIVAL®

The Rival® rear-wheel drive power chair features innovative design that provides exceptional outdoor performance and handling while delivering outstanding tight space maneuverability.

The Rival accepts a complete range of rehab seating and electronics options.

QUANTUM®

Specifications

DRIVE WHEELS	Standard: 14" Solid Optional: 14" Pneumatic	
CASTER WHEELS	Front: 9" Pneumatic	
ANTI-TIP WHEELS	N/A	
MAXIMUM SPEED ⁴	Up to 6 mph (Optional: 8 mph)	
GROUND CLEARANCE	3.25" (frame base)	
TURNING RADIUS ⁵	35.75"	
OVERALL LENGTH	36" without foot riggings	
BASE WIDTH	25"	
SEATING SIZES	TRU-Balance® 2	TRU-Balance® 3
Synergy®/Static Seat	Width: 10-22" Depth: 10-20"	Width: 12-22" Depth: 12-22"
Power Tilt	Width: 14-22" Depth: 14-20"	Width: 12-22" ³ Depth: 12-22" ³
Lift & Tilt	N/A	Width: 12-22" ³ Depth: 12-22" ³
SEAT-TO-FLOOR HEIGHTS	TRU-Balance® 2	TRU-Balance® 3
Synergy®/Static Seat ⁸	16.5-18"	16.5-18.1"
Power Tilt ⁸	16.5-18"	18-19.6"
Lift & Tilt ⁸	N/A	18.5-20.1"
MANUAL TILT/RECLINE ¹⁰	Yes	
BATTERY SIZE ⁶	Group 24	
BATTERY WEIGHT ¹⁴	52 lbs.	
AVAILABLE ELECTRONICS ¹²	75A Q-Logic 2 NE 75A Q-Logic 2 NE+ 100A Q-Logic 2 EX	
BATTERY CHARGER	8A Off-board	
MOTOR PACKAGES	4-Pole 6 mph ¹ or 8 mph ²	
WEIGHT CAPACITY ⁷	300 lbs.	
BASE WEIGHT	129 lbs.	
MEDICARE CODE COVERAGE	Group 3: Standard (K0848/K0849), Single Power (K0856), Multiple Power (K0861)	

- 1) Standard
- 2) Optional
- 3) Requires HD Suspension
- 4) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
- 5) Without front riggings
- 6) Battery capacity may change when adding power positioning systems.
- 7) Weight capacity is dictated by the selection of seating.
- 8) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
- 9) 250 lbs. weight capacity with TRU-Balance® Systems
- 10) The Synergy® Manual Tilt weight capacity is 75 lbs. and maximum dimensions are 14" W x 14" D. Maximum dimensions of the Manual Recline vary by power base and seating system. Call for details.
- 11) Requires Group 34 batteries
- 12) All Quantum Rehab® units can be ordered with specialty controls
- 13) Required with high-speed motors
- 14) Battery weight may vary +/- 3 lbs. based on manufacturer.

Features

- Compact size for outstanding maneuverability in tight spaces
- Maximum speed up to 8 mph²
- High-performance, high-efficiency motor package
- Accepts a wide range of seating and positioning systems
- Sliding battery tray for easy access
- Two suspension packages

Options

- TRU-Balance® 3 Power Tilt
- TRU-Balance® 3 Power Recline
- TRU-Balance® 3 (10") Lift
- Manual and Power Elevating Leg Rests
- TRU-Comfort 2 Seating
- Q-Logic 2 Controller
- Full Lighting Package
- Centermount Adjustable Foot Platform
- Accu-Trac Advanced Tracking Technology
- Power Articulating Foot Platform

Colors



QUANTUM®