

Z-CHAIR

The first power chair that combines super portability, maneuverability and style at a low, low cost.



Specifications

Drive wheels	8" (solid)
Front casters	5" (solid)
Rear anti-tips	3" (solid)
Maximum speed	3.5 mph*
Ground clearance	1.75"
Turning radius	26.75" with foot platform
Length	29" without foot platform
Width	19"
Seat type	Foldable compact lightweight seat
Battery requirements/size	(2) 12V, 12 amp
Range	Up to 8 miles*
Charger	Off-board, 2 amp
Electronics	40 amp, PG GC2
Weight capacity	250 lbs.
Base weight	57 lbs. (front and rear assembled without batteries)
Standard seat weight	25 lbs.
Battery pack weight	20 lbs.
Suspension	Limited
Heaviest piece (rear)	34.5 lbs.
Available colors	Red or Blue
Warranty	5-year limited warranty on frame; 13-month limited warranty on electronics and drive motors

Portability

- Easy-to-remove battery pack and compact, lightweight frame
- Conveniently disassembles into four manageable pieces
- Fits in most vehicle trunks

Performance

- Top speed of 3.5 mph and enough power to get you where you want to go
- Twin motors for traction and precision handling

Maneuverability

- Compact design for a tight turning radius



PRIDE[®]
PERFORMANCE MOBILITY
— LIVE YOUR BEST[®] —

(US) 800-800-8586 ▪ (Canada) 888-570-1113
Visit us on the web: www.pridemobility.com

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

Pride[®]
Mobility Products Corp.