



Kids ROCKTM 2



The Kids ROCKTM 2 chair from Kids UPTM offers highly adjustable and fully customizable seating and positioning to meet the needs of a growing child. Its easy to operate design makes it the perfect choice for children and parents.

The ReactionTM Dynamic Seating System allows a child the freedom to move into flexion or extension in the chair while maintaining proper pelvic positioning and therapeutic support.

SPECIFICATIONS

| | |
|---------------------------------------|--------------------------|
| Weight Limit | 110 lbs. |
| Seat Depth Range | 10" - 16" |
| Seat Width Range | 8" - 14" |
| Back Height Range | 12" - 22" |
| Lower Leg Length Range | 6" - 16" |
| Hip Angle Adjustment | 83°, 90°, 97°, 104° |
| Tilt-in-Space Range | -5° to +45° |
| Push Handle Height Range | 32" - 48" |
| Seat-to-Back Active Range | 35° |
| Knee Active Range | 36° |
| Dimensions (with back folded down) | (W)25" x (H)27" x (L)37" |

(US) 877-454-3787 • FAX: 406-388-7995
(Canada) 888-570-1113 • FAX: 905-682-7256
Visit us on the web: www.quantumrehab.com

The information contained herein is correct at the time of publication.
We reserve the right to alter specifications without prior notice.



QUANTUM REHAB[®]