

**Kids
UP**™

Now
Billable as
Code E1161*

Kids ROCK™ 3



The Kids ROCK™ 3 chair from Kids UP™ offers highly adjustable and fully customizable seating and positioning to meet the needs of a growing client. Its easy to operate design makes it the perfect choice for clients and caregivers.

The Reaction™ Dynamic Seating System allows a child the freedom to move into flexion or extension in the chair while maintaining proper pelvic positioning and therapeutic support.

SPECIFICATIONS

Weight Limit	215 lbs.
Seat Depth Range	14" - 20"
Seat Width Range	12" - 17"
Back Height Range	18" - 26"
Lower Leg Length Range	12" - 20"
Hip Angle Adjustment	83°, 90°, 97°, 104°
Tilt-in-Space Range	-5° to +45°
Push Handle Height Range	30" - 48"
Seat-to-Back Active Range	35°
Knee Active Range	36°
Dimensions (with back folded down)	(W)27" x (H)27" x (L)37"

(US) 877-454-3787 • FAX: 406-388-7995
(Canada) 888-570-1113 • FAX: 905-682-7256
Visit us on the web: www.quantumrehab.com

The information contained herein is correct at the time of publication.
We reserve the right to alter specifications without prior notice.



QUANTUM REHAB®

* MANUAL ADULT SIZE WHEELCHAIR,
INCLUDES TILT-IN-SPACE