

Manual Chairs

L-400

The L-400 by Pride Mobility Products is designed to meet the needs of active users. It's extremely lightweight, highly durable, and provides enhanced fit and function. The L-400 offers all this in a very cost effective package.

Specifications

Overall width:	66 cm (26")
Overall height:	97 cm (38")
Overall length:	84 cm (33")
Folded width:	30 cm (11")
Total weight:	13 kg (29 lbs.)
Weight capacity:	113 kg (250 lbs.)
Seat dimensions:	46 cm x 40 cm (18" x 16")
Seat-to-floor height:	47 cm - 50 cm (19"-20")
Seat back height:	38 cm - 48 cm (15"-19")



Features

- Adjustable tension seat upholstery increases comfort and positioning
- Flip-up, desk-length armrests improve accessibility
- Multiple-position rear wheel and caster mounts optimize centre-of-gravity and seat-to-floor height adjustments
- Quick-release, removable rear wheels
- Full-profile 61 cm x 3.5 cm rear wheels for enhanced outdoor capabilities
- Long-handle wheel locks ensure comfort and security
- Through-bolted, industry-compatible leg rests with aluminum foot plates for durability
- Warranty: 1 year limited

The information contained herein is correct at the time of publication; due to Pride's commitment to constant improvement and development, we reserve the right to alter specifications without prior notice*.

Pride
Mobility Products
New Zealand Ltd.

38 Lansford Crescent
Avondale, Auckland NZ 1007
www.pridemobility.com

Tel: 09 828 0337
Fax: 09 828 0339

XRHT

The XRHT folding manual wheelchair delivers excellent features and solid all-around performance at an exceptional price. Convenience touches like a detachable armrest and detachable swing-away footrest and a weight capacity of 113 kg (250 lbs.) help make the XRHT a superb manual wheelchair choice.



Specifications

Overall width:	66 cm (26")
Overall height:	92 cm (36")
Overall length:	105 cm (41")
Folded width:	32 cm (12.5")
Total weight:	18 kg (40 lbs.)
Weight capacity:	113 kg (250 lbs.)
Seat dimensions:	46 cm x 41 cm (18" x 16")
Seat-to-floor height:	50 cm (19.5")
Seat back height:	46 cm (18")

Features

- Detachable armrest and detachable swing-away footrest
- PU (polyurethane)
- 61 cm (24") rear magwheel with black solid tyre
- Maximum weight capacity of 113 kg (250 lbs.)
- 20 cm (8") front casters
- Aluminum caster fork
- Black durable plastic footplate
- Black coated frame
- Black plastic side panel
- Warranty: 1 year limited

QRHT

The QRHT folding manual wheelchair delivers excellent features and solid all-around performance. Features include detachable armrests and detachable swing-away footrest, quick release rear spoked wheels and a weight capacity of 113 kg (250 lbs.).



Specifications

Overall width:	68 cm (26.5")
Overall height:	89 cm (35")
Overall length:	104 cm (41")
Folded width:	32 cm (12.5")
Total weight:	19 kg (42 lbs.)
Weight capacity:	113 kg (250 lbs.)
Seat dimensions:	46 cm x 41 cm (18" x 16")
Seat-to-floor height:	47 cm (18.5")
Seat back height:	46 cm (18")

- Detachable armrest and detachable swing-away footrest
- PU (polyurethane)
- 61 cm (24") rear, quick release spoked magwheel with black solid tyre
- Maximum weight capacity of 113 kg (250 lbs.)
- 20 cm (8") front casters
- Aluminum caster fork
- Black durable plastic footplate
- Black coated frame
- Black plastic side panel
- Warranty: 1 year limited

BRHT

The BRHT transit chair is designed to deliver solid performance and maximum portability. Features include 12" rear wheels, 8" front castors, detachable arm and leg rests and a folding backrest that makes the BRHT perfect for storage during travel.



Specifications

Overall width:	58 cm (23")
Overall height:	94 cm (37")
Overall length:	88 cm (34.5")
Folded width:	22 cm (8.5")
Total weight:	17 kg (37 lbs.)
Weight capacity:	113 kg (250 lbs.)
Seat dimensions:	46 cm x 41 cm (18" x 16")
Seat-to-floor height:	50 cm (19.5")
Seat back height:	46 cm (18")

- Transit chair
- Detachable armrest and detachable swing-away footrest
- PU (polyurethane)
- 30 cm (12") rear wheels
- Maximum weight capacity of 113 kg (250 lbs.)
- 20 cm (8") front casters
- Aluminum caster fork
- Black durable plastic footplate
- Black coated frame
- Black plastic side panel
- Folding backrest for maximum portability
- Warranty: 1 year limited