



Options & Extras

Milford Person Lift - Sling Options



Standard Slings
X-Small to X-Large



Amputee Slings
Small to Large



Fleece Slings
for sensitive skin



Net Slings for wet
environments

Milford Person Lift - Mounting Options



Cill Bracket & Low Dash Adapter
For allowing use over large door
cills or extended dashboards



Dropped Extension Arms
Available in 2" - 4" - 6" - 10"
and 12" left hand and right
hand drop extensions



Wall Bracket Spacers
Available in 3" - 6" - 9" - 12"
extension lengths




Spreader Bar Covers
Soft spreader bar covers
to protect the user against
involuntary movements

Milford Person Lift Specifications

Lifting capacity:	330lbs
Lifting arm weight:	18lbs
Maximum lift travel:	22.5"
Power supply:	12V dc
Maximum current draw:	22 Amps
Product materials:	Aluminum and Mild/Stainless Steels
Material finish:	Powder Coated - Gray

Locate a Dealer:

 (800) 800-8586 (U.S.)
(888) 570-1113 (Canada)

 info@pridemobility.com

 www.prideliftsandramps.com



10,000
Sold
Worldwide



Simple and easy transfers from the wheelchair to a **car of your choice**

The versatile transfer solution with multiple applications



In the **car**



Around the **home**



On the **move**



Wheelchair To Car Transfer

Whether you have one vehicle you love or several you can't be without, you can transfer from wheelchair to the **front back, passenger or driver side** at your discretion. The Milford Person Lift allows for simple wheelchair to car transfer **without limiting your choice of vehicles or seating positions.**



1) Connect the sling to the Milford



2) Lift and maneuver into the vehicle



3) Lower down onto the seat



4) Detach the sling and remove the Milford

Transfer into an enormous range of vehicles



✓ Car



✓ Minivan



✓ 4x4



✓ SUV



✓ Front seat
✓ Rear seat



One lift - Multiple Transfers

Milford Wall Bracket - For Transfers **Around The Home**



Transfer **In The Living Room**



Transfer **In The Bathroom**



Transfer **In The Bedroom**

Milford Smartbase - For Transfers **On The Go**

Make **effortless transfers** wherever you go with the Milford Smartbase.



Ideal for taking on vacation, visiting relatives or even going to work, the Milford Smartbase enables **transfers from one seated position to another** wherever you might be.



On Vacation



At Work



Visiting Friends