


## Celebrity® XL 4-wheel

Model number:	SC4450DX
Weight capacity:	400 lbs.
Maximum speed <sup>1</sup> :	Up to 5.8 mph
Ground clearance <sup>2</sup> :	4.5"
Turning radius:	62"
Overall length <sup>2</sup> :	51.25"
Overall width <sup>2</sup> :	25.5"
Seat-to-deck <sup>2</sup> :	16.25" - 17.5"
Seat-to-ground:	24.5" - 25.25"
<b>Tires</b>	
Front:	4" x 13" pneumatic
Rear:	4" x 13" pneumatic
Anti-tips:	Rear
Suspension:	Front/Rear
Battery charger:	Off-board, 5 amp
Per-charge range <sup>1,5</sup> :	Up to 18 miles
Battery requirements <sup>3</sup> :	(2) 12V deep-cycle
Compatible battery size(s):	40AH or U-1
<b>Weight<sup>4</sup></b>	
Battery (2 required):	40AH: 32.5 lbs. each; U-1: 23.5 lbs. each
Total w/o batteries <sup>2,6</sup> :	186.5 lbs.
Heaviest piece when disassembled <sup>2</sup> :	78 lbs. (front section)

Colors:



<sup>1</sup>Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).

<sup>2</sup>Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).

<sup>3</sup>AGM or Gel-Cell type required.

<sup>4</sup>Battery weight may vary based on manufacturer

<sup>5</sup>Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance test conducted at maximum weight capacity.

<sup>6</sup>Includes standard seat and weight.

*NOTE: All specifications subject to change without notice.*

### Benefits

- Four-wheel design and full suspension for excellent outdoor performance
- Durable high-impact bumpers provide extraordinary resistance to scratches, chips and dings
- Full lighting package including directional signals and angle adjustable, dual bulb lower headlight
- Fully adjustable, deluxe reclining high-back seat with full-range sliding mechanism
- 110-amp, microprocessor-based motor controller for exceptional power, torque and efficiency
- Easily accessible charger port and circuit breaker
- 3-year limited warranty
- Independent angle-adjustable delta tiller

### Options

- Cane/crutch holder
- Oxygen tank holder
- Rear basket
- Safety flag
- Walker holder
- Durable weather cover
- Solid tires

## Silver Star® Lifts & Ramps ★

### QUAD XL SCOOTER LIFT

The Silver Star® Quad XL Scooter Lift is a sturdy, longer scooter lift for use with Pride's larger products.



This lift is compatible with the Celebrity® XL.

# Celebrity® XL

INDEPENDENCE AT A PRICE  
YOU CAN AFFORD.

## CELEBRITY® XL

The Celebrity® XL Outdoor Scooter delivers unprecedented performance, unmatched durability and an impressive range of standard features. The Celebrity® XL is rugged enough to tackle rough terrain while remaining comfortable and easy to use.



**Pride®**  
Mobility Products Corp.

800-800-8586 (US) • 888-570-1113 (Canada) • [www.pridemobility.com](http://www.pridemobility.com)