

VICTORY[®] 9

THE ULTIMATE 9" SCOOTER
WITH LUXURY FEATURES AT
AN ECONOMY PRICE.

VICTORY[®] 9

The **Victory[®] 9** delivers the high-performance operation and feature-rich design of the Victory series of scooters. Feather-touch disassembly, exclusive, low-profile tires, backlit battery gauge, and wraparound delta tiller all combine to make the Victory 9 an exceptional value.



Pride[®]
Mobility Products Corp.
— LIVE YOUR BEST™ —

800-800-8586 (US) • 888-570-1113 (Canada) • www.pridemobility.com

Victory® 9	3-wheel	4-wheel
Model number:	SC609	SC709
Weight capacity:	300 lbs.	300 lbs.
Maximum speed ¹ :	Up to 5.25 mph	Up to 5.25 mph
Ground clearance ² :	2"	2.25"
Overall length ² :	43.5"	45.5"
Overall width ² :	22.25"	22.5"
Turning radius:	40.25"	51.5
Range per charge ^{1,6} :	13 miles ⁵	
Tire size:	Front: 8" solid Rear: 9" solid	Front: 9" solid Rear: 9" solid
Total weight without batteries ^{2,7} :	117.5 lbs.	131 lbs.
Heaviest piece when disassembled ² :	Front: 48.5 lbs.	Front Rear: 45 lbs.
Standard seating:	Type: Foldable molded plastic Weight: 28 lbs. Material: Black vinyl Dimensions: 18" overall width 17" effective depth	
Seat-to-deck range:	17.25" - 18.75"	
Drivetrain:	Rear-wheel drive, sealed transaxle, 24-volt DC motor	
Dual braking system:	Regenerative and electromechanical	
Battery requirements ^{3,4} :	Two 12-volt, deep cycle Size: U-1 Weight: U-1:24.5 lbs. each	
Battery charger:	3-amp, off-board	

Colors:



Viper Blue

Candy Apple Red

¹Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. This specification can be subject to a variance of (+ or -) 10%.

²Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of + or - 3%.

³AGM or Gel-Cell type required.

⁴Battery weight may vary based on manufacturer.

⁵This range was obtained with the 4-wheel model.

⁶Tested in accordance with ANSI/RENSA, WC Vol 2, Section 4, and ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

⁷Includes standard seat weight.

Features


- Exclusive, one-hand feather-touch disassembly provides the ultimate quick transport system
- Pride's exclusive stylish, lightweight, low-profile, non-scuffing wheels
- Auto-connecting front-to-rear harness
- Delta tiller with wraparound handles
- Long-lasting, high-intensity LED headlight offers optimal pathway illumination
- Includes pinchless hinges & eight-position seat rotation
- 300 lbs. weight capacity

Options

- Forearm crutch holder
- Cane/crutch holder
- Oxygen tank holder
- Cup holder
- Rear basket
- Safety flag
- Walker holder
- Weather cover
- Saddlebag

Silver Star Lifts & Ramps ★

- Backpacker® Plus
- Full Platform
- Silver Boom 300
- Backpacker® MV



Full Platform

Backpacker® Plus

These lifts are compatible with the Victory® 9.



The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed & range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%.

800-800-8586 (US) • 888-570-1113 (Canada) • www.pridemobility.com

