

# VICTORY<sup>®</sup> 9 PS

FULLY LOADED WITH AN  
IMPRESSIVE LIST OF FEATURES

## VICTORY<sup>®</sup> 9 PS

The **Victory<sup>®</sup> 9 PS** delivers a long list of exciting features including a power elevating seat, low-profile tires, Feather-touch disassembly, backlit battery gauge and wraparound delta tiller. With all of these luxury touches standard, along with high performance operation and exhilarating style, the Victory 9 PS is a truly exceptional value.


*Standard  
with a  
power  
elevating  
seat*



# Victory® 9 PS 3-wheel

Model number:	SC609PS
Weight capacity:	300 lbs.
Maximum speed <sup>1</sup> :	Up to 5.5 mph
Ground clearance <sup>2</sup> :	2" at midframe
Overall length <sup>2</sup> :	43.5"
Overall width <sup>2</sup> :	23"
Turning radius:	41"
Range per charge <sup>1,5</sup> : <small>(Theoretical value based on test standards)</small>	13 miles
Tire Size:	Front: 8" solid Rear: 9" solid
Total Weight Without Batteries <sup>1,2</sup> :	131 lbs.
Heaviest piece when disassembled <sup>2</sup> :	Front: 59.5 lbs.
Standard seating:	Type: Foldable molded plastic Weight: 28 lbs. Material: Black Vinyl Dimensions: 18" overall width 17" overall depth
Seat-to-Deck Range <sup>2</sup> :	17.25" - 23.25"
Seat-to-Ground Range <sup>2</sup> :	23.5" - 29.5"
Drive System:	Rear-wheel drive, sealed transaxle, 24-volt DC motor
Dual Braking System:	Regenerative, and electromechanical
Battery Requirements <sup>3</sup> :	Two 12-volt, deep cycle Size: U-1 Weight: U-1:24.5 lbs. each <sup>4</sup>
Battery Charger:	3-amp, off-board

Colors:



<sup>1</sup>Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).

<sup>2</sup>Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).

<sup>3</sup>AGM or Gel-Cell type required.

<sup>4</sup>Battery weight may vary based on manufacturer

<sup>5</sup>Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance test conducted at maximum weight capacity.

*NOTE: All specifications subject to change without notice.*

## Features

- Standard power-elevating seat
- Exclusive, one-hand Feather-touch disassembly provides the ultimate quick transport system
- Pride's exclusive stylish, light-weight, low profile, non-scuffing wheels
- Auto-connecting front-to-rear harness
- Delta tiller with wraparound handles
- Long-lasting, high-intensity LED headlight offers optimal pathway illumination
- Includes pinchless hinges and eight-position seat rotation
- 300 lbs. weight capacity

## Options

- Forearm crutch holder
- Cane/crutch holder
- Oxygen tank holder
- Cup holder
- Rear basket
- Safety flag
- Walker holder
- Weather cover
- Saddlebag

## Silver Star Lifts & Ramps

- Backpacker® Plus
- Full Platform
- Backpacker® MV

