



# GO-GO

## TRAVEL SCOOTER



The Go-Go® Travel Scooter takes easy transportability to a whole new level! It quickly and conveniently disassembles, so feel free to take it anywhere, even without a lift!

- Easy-to-remove battery pack
- Flat-free tires
- Easily adjustable tiller
- 4-position seat
- Weight capacity of 250 lbs.
- Maximum speed up to 4 mph
- Tremendous overall value

**Super  
portability...  
at a low,  
low price!**



**800-800-8586 (U.S.)  
888-570-1113 (Canada)  
www.pridemobility.com**

## Specifications

Model numbers	SC40 (3-wheel); SC44 (4-wheel)
Weight capacity	250 lbs.
Maximum speed	Up to 4 mph
Ground clearance	1.75"
Turning radius	3-wheel: 32" 4-wheel: 43"
Overall length	3-wheel: 37" 4-wheel: 39"
Overall width	19"

## Tires

Front	2" x 7.5" solid
Rear	2" x 7.5" solid
Anti-Tips	Rear
Standard seating	Foldable, compact, vinyl-covered molded foam
Drive train	24-volt DC motor, sealed mini-transaxle
Dual braking system	Electronic, regenerative and electromechanical

Battery charger	Off-board, 2 amp
Per-charge range	Up to 10 miles
Battery requirements	(2) 12V deep-cycle
Compatible batteries	12 AH

## Weight

Batteries (2 required)	9 lbs. each <sup>1</sup>
Total weight	3-wheel: 87.5 lbs. 4-wheel: 96.5 lbs.
Heaviest piece	3-wheel: 30 lbs. (rear section) 4-wheel: 31 lbs. (front section)

## Standard Colors



Red



Blue

**1.** Battery pack weighs 20 lbs.

The Go-Go easily disassembles into 4 lightweight pieces for convenient transport and storage.



**The Go-Go Travel Scooter is not classified as an E-1230 (POV/Scooter) for Medicare reimbursement purposes.**

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed & range vary with user weight, terrain type, battery charge, battery condition and tire pressure.

**800-800-8586 (U.S.) • 888-570-1113 (Canada) • [www.pridemobility.com](http://www.pridemobility.com)**

**Pride**  
Mobility Products Corp.

## Features

- Disassembles into 4 lightweight pieces for convenient transport and storage
- 4-position seat
- Easily adjustable tiller
- Easily accessible freewheel release lever
- Easy-to-remove battery pack
- Flat-free tires
- Easy lock-up system for quick and simple disassembly
- Scratch resistant shroud



Small enough to fit in your closet