




### Celebrity® XL 4-wheel


### Legend® XL 4-wheel

Model number:	SC4450DX	SC3450
Weight capacity:	400 lbs.	400 lbs.
Maximum speed:	Up to 5.8 mph	Up to 5.8 mph
Ground clearance:	4.5"	4"
Turning radius:	62"	65"
Overall length:	51.25"	52"
Overall width:	25.5"	24"
Seat-to-deck:	16.25" - 17.5"	16.5" - 19.5"
Seat-to-ground:	24.5" - 25.25"	23" - 26"
<b>Tires</b>		
Front:	4" x 13" pneumatic	4" x 13" pneumatic
Rear:	4" x 13" pneumatic	4" x 13" pneumatic
Anti-tips:	Rear	Rear
Suspension:	Front/Rear	Front/Rear
Battery charger:	O -board, 5 amp	O -board, 5 amp
Per-charge range:	Up to 18 miles	Up to 25 miles
Battery requirements:	(2) 12V deep-cycle	(2) 12V deep cycle
Compatible battery size(s):	40AH or U-1	NF-22 or U-1
<b>Weight<sup>1</sup></b>		
Battery (2 required):	40AH: 32.5 lbs. each; U-1: 23.5 lbs. each	NF-22: 38 lbs. each; U-1: 24.5 lbs. each
Total w/o batteries:	186.5 lbs.	207 lbs.
Heaviest piece when disassembled:	78 lbs. (front section)	82 lbs. (rear section)

Colors:



Viper Blue      Candy Apple Red



Viper Blue      Candy Apple Red      Champagne

<sup>1</sup>Battery weight may vary based on manufacturer

### Benefits

- Four-wheel design and full suspension for excellent outdoor performance
- Durable high-impact bumpers provide extraordinary resistance to scratches, chips and dings
- Full lighting package including directional signals and angle adjustable, dual bulb lower headlight
- Fully adjustable, deluxe reclining high-back seat with full-range sliding mechanism
- 110-amp, microprocessor-based motor controller for exceptional power, torque and efficiency
- Easily accessible charger port and circuit breaker
- 3-year limited warranty
- Independent angle-adjustable delta tiller (Celebrity® XL only)
- Accepts NF-22 batteries for maximum range per charge (Legend® XL)

### Options

- Cane/crutch holder
- Oxygen tank holder
- Rear basket
- Safety flag
- Walker holder
- Durable weather cover
- Solid tires

Pride Silver Star Lifts & Ramps 

### QUAD XL SCOOTER LIFT

The Pride Silver Star® Quad XL Scooter

Lift is a sturdy, longer scooter lift for use with Pride's larger products.



This lift is compatible with the Celebrity® XL and Legend® XL.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed & range vary with user weight, terrain type, battery charge, battery condition and tire pressure.

800-800-8586 (US) • 888-570-1113 (Canada) • www.pridemobility.com

**Pride®**  
Mobility Products Corp.

# Celebrity® XL & Legend® XL

INDEPENDENCE AT A PRICE  
YOU CAN AFFORD.

## CELEBRITY® XL & LEGEND® XL

The Celebrity® XL and Legend® XL Outdoor Scooters deliver unprecedented performance, unmatched durability and an impressive range of standard features. Both scooters are rugged enough to tackle rough terrain while remaining comfortable and easy to use.



**Pride**<sup>®</sup>  
Mobility Products Corp.

800-800-8586 (US) • 888-570-1113 (Canada) • [www.pridemobility.com](http://www.pridemobility.com)

CELXL\_LEGXL\_09\_08