

Rally®

3 and 4-wheel scooters

The Rally® is sure to cause a stir with its futuristic styling, technologically advanced performance and overall great value—a combination irresistible and unmatched in the small scooter class. With the Rally, you'll stake your claim on a lifestyle of fun and freedom.



SC155
Rally® 4-wheel scooter
shown in Candy Apple
Red



SC151
Rally® 3-wheel scooter
shown in Viper Blue

***Stylish Design & Premium Performance
...at an Excellent Value***

Pride®
Mobility Products Corp.

1-800-800-8586 (U.S.)
1-888-570-1113 (Canada)
www.pridemobility.com

Benefits

- Fully adjustable tiller for supreme comfort
- Controls at your fingertips for ease of use
- Standard headlight for increased safety
- 4-wheel model now available for stable outdoor performance
- Restyled sleek tiller
- Now with a 300 lbs. weight capacity
- Seat offers height adjustment and swivel with flip-up armrests for comfortable positioning
- Up to a 25-mile range per battery charge; on-board charger and easy-access external circuit breaker
- Speed range adjustable up to 4.25 mph
- External freewheel release lever
- Front basket for convenience
- Aluminum alloy mag wheels for great style

Options

Forearm crutch holder • Taillight • Cane/crutch holder • Oxygen holder • Cup holder • Walker holder • Rear view mirror • Rear basket • Safety flag • Two battery options (32 AH; 17 AH) Dust cover • External battery charger



RALLY-NG-1001A


1-800-800-8586 (U.S.)
1-888-570-1113 (Canada)
www.pridemobility.com

Specifications*

Model Numbers:	(3-wheel) SC151 (4-wheel) SC155
Available Colors:	Candy Apple Red, Viper Blue
Overall Length:	(3-wheel) 41" (4-wheel) 44"
Overall Width:	21"
Total Weight Without Batteries:	(3-wheel) 89 lbs. (4-wheel) 116 lbs.
Heaviest Piece When Disassembled:	(3-wheel) 40.5 lbs. (4-wheel) 45 lbs.
Turning Radius:	(3-wheel) 39" (4-wheel) 54"
Maximum Speed:	4.25 mph
Range Per Charge:	up to 20-25 miles
Ground Clearance:	3"
Weight Capacity:	300 lbs.
Standard Seating:	Foldable molded plastic
Drive System:	Sealed transaxle
Dual Braking System:	Electronic, regenerative and electromechanical
Tires (front):	2" x 8"
Tires (back):	2.5" x 9"
Battery Requirements:	
Type:	(2) 12V deep cycle (SLA or gel cell)
Size:	32AH or optional 17AH
Battery Charger:	On-board 3A
Warranty:	3-year limited

*Information subject to change.

Colors



Blue



Red