



sonic™

3-wheel scooters



Pride's® Sonic™ Scooter offers true innovation in the realm of ultra-light, transportable scooters. It easily disassembles for storage in even the smallest car trunks and, in use, is perfect for navigating through crowded stores and maneuvering in tight hallways. With the Sonic, you'll have a scooter that, like you, is ready for fun at a moment's notice.

**The ultimate
portable scooter.
Take the Sonic™
wherever
you go!**



The Sonic™ disassembles into 4 pieces.



**800-800-8586 (U.S.)
888-570-1113 (Canada)
www.pridemobility.com**

Pride®
Mobility Products Corp.

Specifications

Model number(s)	SC 50
Weight capacity	250 lbs.
Maximum speed	Up to 4.25 mph
Ground clearance	1.75"
Turning radius	34"
Overall length	37"
Overall width	21"

Tires

Front	2" x 7.5" solid
Rear	2.5" x 7.5" solid
Anti-tips	Rear
Standard seating	Foldable molded plastic
Power elevating seat	Option: NO
Suspension	NA
Drive train	Sealed transaxle, 24 volt DC motor
Dual braking system	Electronic, regenerative and electromechanical
Battery charger	On-board, 2 amp
Battery Requirements	(2) 12V deep-cycle (SLA or gel cell)
Battery size(s)	12AH

Weight

Battery	9 lbs. each
Total w/o batteries	93 lbs.
Heaviest piece when disassembled	35 lbs. (rear section)

Colors



Blue



Red



Black



Yellow

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed & range vary with user weight, terrain type, battery charge, battery condition and tire pressure.

**Take the Sonic™
wherever you go!**

Benefits

- Disassembles into 4 easy-to-transport pieces
- Heaviest piece weighs only 35 lbs.
- Adjustable tiller and fingertip controls for comfort and ease of use
- Easy to maneuver around the home
- Solid, flat-free tires for a worry-free ride
- Up to a 10-mile range per battery charge; removable battery pack
- Maximum weight capacity of 250 lbs.
- Speed range adjustable up to 4.25 mph

Options

- Forearm crutch holder
- Lap belt
- Cane/crutch holder
- Wishbone crutch holder
- Oxygen tank holder
- Cup holder
- Rear basket
- Safety flag
- Walker holder
- Heavy-duty dust cover
- Saddle bag
- External battery charger
- Quad cane holder
- Front basket
- Clip, cane attachment

