



- Feather-touch disassembly permits simple frame separation with only one hand
- Maximum speed up to 8 MPH
- Up to 11 miles per charge range
- Delta tiller with wraparound handles
- Black vinyl reclining high-back seat with sliders, height adjustable headrest, flip-up, height and width adjustable armrest
- Long-lasting, high-intensity LED headlight offers optimal pathway illumination with directional signals
- Easy access tiller-mounted charger port
- Front and rear suspension
- 10" drive wheels
- 350 lbs. weight capacity



Delta Tiller
Wraparound handle design with easy to use forward and reverse controls



Feather Touch
One-Hand Disassembly



Crush proof front basket



Black, non-scuffing tires



Long lasting, high intensity LED headlight



Easy grip tiller adjustment knob



Backlit battery gauge



Easy access tiller-mounted charger port

COLOR CHOICES



Candy Apple Red



Viper Blue

VICTORY[®] SERIES



PRIDE[®] MOBILITY SCOOTERS
LIVE YOUR BEST[®]

Pride[®]
Mobility Products Corp.

(US) 800-800-8586 • (Canada) 888-570-1113

pridemobility.com

THE ULTIMATE IN STYLE
& PERFORMANCE[®]

pridemobility.com

VICTORY[®] 9



- Feather-touch disassembly permits simple frame separation with only one hand
- Maximum speed up to 5.25 MPH
- Up to 13 miles per charge range
- Pride's exclusive black, non-scuffing tires
- Lightweight pinchless hinge seat, eight position rotation, sliders and viscoelastic foam inserts for increased comfort
- Delta tiller with wraparound handles
- Long-lasting, high-intensity LED headlight offers optimal pathway illumination
- Easy access tiller-mounted charger port
- 9" drive wheels
- 300 lbs. weight capacity

COLOR CHOICES



Candy Apple Red Viper Blue

VICTORY[®] 10



- Feather-touch disassembly permits simple frame separation with only one hand
- Maximum speed up to 5.25 MPH
- Up to 15.5 miles per charge range
- Pride's exclusive black, non-scuffing tires
- Lightweight pinchless hinge seat, eight position rotation, sliders and viscoelastic foam inserts for increased comfort
- Delta tiller with wraparound handles
- Lightweight seat with sliders
- Long lasting, high intensity LED headlight offers optimal pathway illumination
- Easy access tiller-mounted charger port
- 10" drive wheels
- 400 lbs. weight capacity

COLOR CHOICES



Candy Apple Red Viper Blue

VICTORY[®] 10 LX WITH CTS SUSPENSION



- CTS Suspension (Comfort-Trac front and rear suspension)
- Maximum speed up to 5.6 MPH
- Up to 16.16 miles per charge range
- Feather-touch disassembly permits simple frame separation with only one hand
- Pride's exclusive black, non-scuffing tires
- Black vinyl reclining high-back seat with sliders, height adjustable headrest, flip-up, height and width adjustable armrest
- Full LED lighting package that includes curb and rear lights and directional signals
- Easy access tiller-mounted charger port
- 10" drive wheels
- 400 lbs. weight capacity

COLOR CHOICES



Candy Apple Red True Blue