

GO-GO

ELITE TRAVELLER® *SPORT*

4-WHEEL SCOOTER

THE ULTIMATE HIGH- PERFORMANCE SCOOTER

The stylish, high-performing **Go-Go Elite Traveller® Sport** features an innovative design with a full complement of unique standard features. A charger port conveniently located in the tiller, front and rear LED lighting, 9-inch wheels at front and rear, an additional 2-inch deck space and an increased weight capacity of 23 Stone 3 lbs. are just a few features that make the Go-Go Elite Traveller Sport an exceptional scooter choice.



Go-Go Elite Traveller® Sport 4-wheel



Wraparound Delta Tiller



Rear LED Light



Larger 9" Wheels



/Pride.MobilityUK



/PrideMobilityUK

PRIDE® MOBILITY SCOOTERS
LIVE YOUR BEST®

(Tel) 01869 324 600 • (Fax) 01869 323 070
32 Wedgwood Road, Bicester, Oxfordshire OX26 4UL
pridemobility.co.uk

SPECIFICATIONS

GO-GO ELITE TRAVELLER® SPORT 4W: S74

Weight capacity:	23 Stone 3 lbs. (147 kg)
Maximum speed ¹ :	Up to 4 mph (6.4 km/h)
Ground clearance ² :	3.8 cm (1.5") at motor
Turning radius ² :	128.3 cm (50.5")
Overall length ² :	109.2 cm (43")
Overall width ² :	54 cm (21.25")
Front tyres:	7.62 cm x 22.9 cm (3" x 9") solid
Rear tyres:	7.62 cm x 22.9 cm (3" x 9") solid
Suspension:	Suspension seat post (no frame suspension)
Total weight without batteries ⁴ :	44 kg (100 lbs.)
Heaviest piece when disassembled:	20 kg (44 lbs.) (front section)
Drive system:	Rear-wheel drive, sealed transaxle, 24-volt DC motor
Dual braking system:	Regenerative and electromechanical
Standard seating:	Type: Compact foldable molded plastic Weight: 10 kg (23.5 lbs.) Material: Black vinyl Dimensions: 43 cm (17") overall width 43 cm (17") effective depth
Battery requirements ³ :	Two 12-volt, deep cycle (standard) Size: 2 x 18 AH (standard) Weight ⁵ : 18 AH: 13.4 kg (29.5 lbs.) battery pack
Battery charger:	Off-board, 2-amp
Standard colours:	



Red



Blue

2 sets of exciting colours with every scooter!

FEATURES

- Front and rear LED lighting
- 5-piece disassembly
- Additional 2 inches of deck space for added leg room
- Exclusive black, low-profile, non-scuffing tyres with new tri-spoke mags
- Delta tiller with wraparound handles
- Auto-connecting front to rear cable
- Includes 2 sets of easily changeable coloured panels - Red & Blue
- Front frame-mounted seat post offers maximum stability
- Modular design for easy serviceability
- Microprocessor-based controller offers optimal power management and added safety features
- Front tiller-mounted basket
- Charger port conveniently located in the tiller
- Suspension seat post

¹ Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tyre condition. These specifications can be subject to a variance of (+/- 10%).
² Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
³ AGM or gel cell type required.
⁴ Includes standard seat.
⁵ Battery weight may vary based on manufacturer.

NOTE: All specifications subject to change without notice.

The information contained herein is correct at the time of publication; due to the manufacturer's commitment to constant improvement and development, they reserve the right to alter specifications without prior notice. Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressure, ambient temperature, and the gradient and surface of the road or pavement.

NOTE: This product conforms to all applicable ANSI/RESNA testing requirements and ISO 7176 series EN12184 standards. All specifications subject to change without notice.



(Tel) 01869 324 600 • (Fax) 01869 323 070
 32 Wedgwood Road, Bicester, Oxfordshire OX26 4UL

pridemobility.co.uk

