

SPARKY

The Sparky is designed to deliver full rehab capability in a compact, powerful paediatric package. Engineered to operate in tight child-sized spaces and to keep up with today's active child, the Sparky is an exceptional fit from the playground to the classroom and beyond.



Pride[®]
Mobility Products Ltd.
— LIVE YOUR BEST —

SPECIFICATIONS

Weight capacity	8 Stone 13 lbs. (56.7) kg
Maximum speed¹	Up to 4 mph (6.43 km/h)
Ground clearance	9.5 cm (3.75")
Turning radius³	74.3 cm (29.25")
Overall length³	95.9 cm (37.75") with front riggings
Overall width³	61 cm (24")
Drive wheels	30.48 cm (12") (pneumatic)
Anti-tips	7.6 cm (3") solid (rear)
Castors	20.32 cm (8") (pneumatic) (front)
Suspension	Rear
Drivetrain	Two-motor, rear-wheel drive
Braking system	Intelligent braking (electronic, regenerative disc brakes)
Standard electronics	70 amp PG VR2 controller
Standard seat	Sling
Specialty controls	Yes (optional attendant controller)
Battery charger	Off-board, 5 amp
Battery requirements	(2) 12 volt deep-cycle
Battery size	40 AH
Battery weight²	10.5 kg (23 lbs.) each
Base weight	59.2 kg (130.5 lbs.) (including seat)
Standard seat weight	N/A (seat is integrated with base)
Max. seating dimensions	35.5 cm x 35.5 cm (14" x 14")
Warranty	Lifetime warranty on frame; 1-year limited warranty on electronics; 18-month warranty on PG Drive control system; 2-year limited warranty on drivetrain

The information contained herein is correct at the time of publication; due to the manufacturer's commitment to constant improvement and development, they reserve the right to alter specifications without prior notice.

Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressure, ambient temperature, and the gradient and surface of the road or pavement.

- 1) Speed varies with user weight, terrain type, battery charge, battery condition and tyre pressure.
- 2) Battery weight may vary based on manufacturer.
- 3) Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3%.



PRIDE
PERFORMANCE MOBILITY
— LIVE YOUR BEST® —

FEATURES

- Fully adjustable seat - 25 cm x 25 cm (10" x 10") through 35.5 cm x 35.5 cm (14" x 14")
- PG VR2 controller
- Power tilt
- Weight capacity up to 8 Stone 13 lbs. (56.7 kg)
- Maximum speed up to 4 mph (6.43 km/h)
- Swing-away footrests standard
- Manual recline (attendant only)
- Headrest (standard)

OPTIONS

- Manual elevating leg rests
- Attendant control

Sparky

Shown in tilted position with sling seat, swing-away leg rests and PG VR2 controller



COLOUR



White

32 Wedgwood Road, Bicester, Oxfordshire OX26 4UL
(Tel) 01869 324 600 • (Fax) 01869 323 070 • www.pridemobility.co.uk