PRODUCT SPECIFICATIONS SHEET

Model: Select HD

Manufacturer: Pride Mobility Products Corp. Internet Address: www.pridemobility.com

Document Information: INFINFB2598/Rev A/April 2011

| Suspension: | Limited |
|---------------------------------|-----------------------------------------------------------------------------------------------------------------------|
| Drive Wheels: | 14 in. (35.56 cm), solid |
| Caster Wheels: | 8 in. (20.32 cm), solid, rear-mounted |
| Anti-tip Wheels: | 6 in. (15.24 cm), solid, front-mounted |
| Maximum Speed: ¹ | Up to 4.0 mph (6.43 km/h) |
| Range: ^{1, 4} | Up to 13.45 miles (21.65 km) |
| Brakes: | "Intelligent Braking" electronic regenerative, disc park brake |
| Ground Clearance: ² | 2.5 in. (6.35 cm) at bottom of base |
| Maximum Safe Slope: | 6° (10.5%) |
| Turning Radius: ² | 21 in (53.3 cm) without front riggings 27 in. (68.58 cm) with front riggings |
| Overall Size: ^{2,6} | Length: 37.75 in. (95.88 cm) without front riggings 42.25 in. (107.31 cm) with front riggings Width: 25 in. (63.5 cm) |
| Drivetrain: | Two motor, Mid-wheel drive |
| Batteries: ³ | Two 12-volt, NF-22 batteries |
| Battery Charger: | 5-amp, off-board |
| Weight Capacity: | 450 lbs. (204 kg) |
| Component Weights: ⁵ | Base: 96 lbs. (43.54 kg) Standard Seat (High Back): 43 lbs. (19.5 kg) Batteries: 38 lbs. (17.24 kg) each |

Varies with base model, user weight, terrain type, battery amp hour rating (Ah), battery charge, battery condition, motors, controller type, tire type, and tire condition. This specification is subject to a variance of +10%, -5%.

NOTE: This product conforms to all applicable ANSI/RESNA testing requirements and ISO 7176 Series EN12184 standards. All specifications subject to change without notice.

Due to manufacturing tolerances and continued product improvement, this specification can be subject to a variance of (+ or -) 3%.

³ AGM or Gel-Cell type required. See VI. "Batteries and Charging."

Tested in accordance with ANSI/RESNA, WC Vol 2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

⁵ Battery weight may vary based on manufacturer.

⁶ Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

