

# VivaLIFT!®

ATLAS PLUS<sup>1\*\*</sup> & TRANQUIL<sup>1</sup>

NOW AVAILABLE WITH  
**CRYPTON®**  
 for a limited time only

## 4-MOTOR RECLINERS



Atlas Plus\*\*  
 shown in Crypton® Cool Grey fabric

FDA Class II  
 Medical Device<sup>1</sup>



Atlas Plus\*\*  
 shown in Cool Grey

FDA Class II  
 Medical Device<sup>1</sup>



Tranquil  
 shown in Espresso

FDA Class II  
 Medical Device<sup>1</sup>



### VIVALIFT!® INDEPENDENT FUNCTIONS...

- The power headrest lifts your head, neck and shoulders for maximum support
- The upward lift offers convenient stand assist for getting in and out of the recliner
- The power lumbar enables you to personalize your comfort
- The footrest extension allows you to relax and stretch out completely
- The backrest moves up and down independently to help you find the perfect position
- Infinite position recline enables the backrest and footrest to move independently from one another so that a variety of positions can be reached

### CRYPTON® FEATURES AND BENEFITS...

- 100% water resistant with a moisture barrier backing
- Stain, mildew, bacteria and odor resistant
- Easy to clean and maintain



Scan here for  
 more information  
 on **VivaLift!®**



**PRIDE®**  
 MOBILITY  
 LIVE YOUR BEST®

# ATLAS PLUS<sup>1\*\*</sup> & TRANQUIL<sup>1</sup>

**VivaLIFT!**<sup>®</sup>  
POWER RECLINER COLLECTION

COLLECTION	ATLAS PLUS**	TRANQUIL	TRANQUIL	TRANQUIL	TRANQUIL
MODEL	PLR-2985M	PLR-935PW	PLR-935S	PLR-935M	PLR-935LT
POSITION	Infinite Lay Flat Relaxed Reclined Upward Lift 5'4"-5'9"	Infinite Lay Flat Relaxed Reclined Upward Lift 5'3" and below	Infinite Lay Flat Relaxed Reclined Upward Lift 5'3" and below	Infinite Lay Flat Relaxed Reclined Upward Lift 5'4"-5'9"	Infinite Lay Flat Relaxed Reclined Upward Lift 5'10" and above
SIZE	Medium	Petite Wide	Small	Medium	Large/Tall
WEIGHT CAPACITY	400 lbs.	400 lbs.	400 lbs.	400 lbs.	400 lbs.
SEAT-TO-FLOOR	19.5"	18"	18"	19.5"	21.5"
SEAT DEPTH	21"	18.5"	18.5"	20.5"	22.5"
SEAT WIDTH	22"	24"	20"	22"	23.5"
TOP OF BACK TO SEAT	26"	24.5"	25"	27"	30"
FULLY PADDED CHAISE	Yes	Yes	Yes	Yes	Yes
FOOTREST EXTENSION	Standard	Standard	Standard	Standard	Standard
SEATING TYPE	Total Comfort	Total Comfort	Total Comfort	Total Comfort	Total Comfort
CHAIR HEIGHT	42"	38"	40"	42"	47"
DISTANCE FROM WALL	16"	12.5"	13"	16"	16"
OVERALL WIDTH	39"	36"	33.75"	37.5"	38"
CHAIR WEIGHT	164 lbs.	137 lbs.	137 lbs.	145 lbs.	155 lbs.
WARRANTY	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated	7-Yr. Pro-rated

## STANDARD FEATURES

- The USB remote adjusts all functions independently and allows you to charge your smart device.
- Lithium battery backup ensures your recliner is operational when the power goes out.


1. VivaLift!® Power Recliners are FDA Class II Medical Devices designed to aid individuals with mobility impairments.

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

 = Indicates shortest available lead time

Recommended HCPCS: E0627-Seat lift mechanism, electric, any type DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage and bundling guidelines may apply. \*\*This product is available exclusively to Authorized Pride Brick & Mortar Providers. No internet sales are permitted. Available in store only. Call your local Pride provider for further information. Other conditions apply.

## OPTIONAL FABRIC: CRYPTON <sup>LIMITED TIME ONLY!</sup>



ESPRESSO

COOL GREY

(US) 800-800-8586 / prideliftchairs.com / pridemobility.com

**PRIDE**<sup>®</sup>  
MOBILITY  
LIVE YOUR BEST<sup>®</sup>