

## Pride Lift Chairs

*Experience the Luxury*




# 805 Wall Hugger Lift Chair

## A 2-Position Chaise Lounger That Saves Space Without Sacrificing Comfort

The 805 2-position lift chair is a great option for individuals who desire the comfort and convenience of a lift chair without sacrificing living space. The 805's stylish look and its unique design allows it to be placed just four inches from the wall, making it the ultimate lift chair for anyone with limited room.

### FABRIC SELECTION



Almond



Latte

Note: Measurements may vary +/- 2.5 cm (1").

The information contained herein is correct at the time of publication; due to Pride's commitment to constant improvement and development, we reserve the right to alter specifications without prior notice.

### SPECIFICATIONS

#### 805

<b>Back Width:</b>	76 cm (30")
<b>Seat-to-Floor:</b>	48 cm (19")
<b>Seat Depth:</b>	51 cm (20")
<b>Seat Width:</b>	51 cm (20")
<b>Top of Back to Seat:</b>	69 cm (27")
<b>Weight Capacity:</b>	136 kg (300 lbs.)

*Experience the Luxury*



**Sawa Medical Supplies, s.a.r.l.**

Tel: 01/804.881 03/639.650

Raouche, Beirut Lebanon

Website: [www.sawamed.com](http://www.sawamed.com)

EEmail: [info@sawamed.com](mailto:info@sawamed.com)

**sawa**  **سَوَا**  
medical supplies s.a.r.l. **شركة  
للمعدات الطبية ش.م.ل.**