

Pride Lift Chairs

Experience the Luxury



LC-107 Lift Chair

An Infinite-Position Lift Chair at an Exceptional Value

The LC-107 is an infinite-position lift chair which allows for a variety of comfort, from watching TV to sleeping. The LC-107 is available at a low cost, making it an exceptional value.



FABRIC SELECTION



Cocoa



Sandal

Note: Measurements may vary +/- 2.5 cm (1").

The information contained herein is correct at the time of publication; due to Pride's commitment to constant improvement and development, we reserve the right to alter specifications without prior notice.

SPECIFICATIONS

LC-107

Back Width:	74 cm (29")
Seat-to-Floor:	51 cm (20")
Seat Depth:	53 cm (21")
Seat Width:	51 cm (20")
Top of Back to Seat:	70 cm (27.5")
Weight Capacity:	147 kg (325 lbs.)

Experience the Luxury



Sawa Medical Supplies, s.a.r.l.

Tel: 01/804.881 03/639.650

Raouche, Beirut Lebanon

Website: www.sawamed.com

EEmail: info@sawamed.com

sawa  **سَوَا**
medical supplies s.a.r.l. **شركة
للمعدات الطبية ش.م.ل.**