

Model	Victory® LX Sport ¹
Product Number	4W: S710LXW
Weight Capacity	400 lbs.
Maximum Speed ²	Up to 8.6 mph (200 lbs.) Up to 8 mph (400 lbs.)
Ground Clearance ⁴	2.2" at motor, 3.15" at mid-frame
Turning Radius ⁴	61.25"
Overall Length ⁴	47"
Overall Width ⁴	25.5"
Seat-to-Ground Height	23" - 25"
Seat-to-Deck Height	16.5" -18.5"
Front Tires	3.5" x 10" solid
Rear Tires	3.5" x 10" solid
Suspension	Front and rear
Range Per Charge ^{2,3}	Up to 19.7 miles at 200 lbs. Up to 13 miles at 400 lbs.
Total Weight ⁶	Without batteries: 160 lbs. With batteries: 220.4 lbs.
Heaviest Piece When Disassembled ⁴	77 lbs. (front section)
Standard Seating	Type: Limited reclining high-back Weight: 43 lbs. Material: Black Vinyl Dimensions: 18" width, 18" depth
Drive System	Rear-wheel drive, sealed transaxle with 24-volt DC motor
Dual Braking System	Regenerative and electromechanical
Battery Requirements ^{5,7}	Size: (2) 50 AH Weight: 30.2 lbs. each
Battery Charger	5 amp, off-board
Warranty	Frame: lifetime limited Drivetrain: 2-year limited Electronics: 2-year limited Battery: 6-months

HCPCS Code: K0899

POWER MOBILITY DEVICE, NOT CODED BY DME PDAC OR DOES NOT MEET CRITERIA

- 1 Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
- 2 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
- 3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
- 4 AGM or gel cell type required.
- 5 Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- 6 Battery weight may vary based on manufacturer.
- 7 Includes standard seat weight.



FDA Class II Medical Device¹

PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]

(US) 800-800-8586 • (Canada) 888-570-1113

FOLLOW US      
pridemobility.com

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

101524

VICTORY[®] LX¹ SPORT

ENJOY A SMOOTH AND SPORTY RIDE



FDA Class II Medical Device¹

LIVE YOUR BEST[®]

FOLLOW US      

pridemobility.com

SMOOTH, COMFORTABLE RIDE

18" x 18" Vinyl reclining high-back seat with sliders, height adjustable headrest and flip-up, height and width-adjustable armrests

Puddle light under tiller

Standard front basket

Easy-to-grip tiller adjustment knob

Full LED lighting package includes curb light, rear lights and directional signals

Storage compartment under the tiller

CTS Suspension includes adjustable shocks for greater comfort



Tiller console



CTS Suspension

Maintain your active lifestyle with the **Victory® LX Sport**. Equipped with **Comfort-Trac Suspension (CTS)**, the Victory LX Sport easily handles more varied terrain while providing a smooth and comfortable ride. Enjoy convenient storage and a courtesy light under the tiller. A long-lasting LED headlight ensures optimal pathway illumination. Experience great features and a sporty style with the Victory LX Sport.

- CTS Suspension includes adjustable shocks for greater comfort
- Exclusive one-hand feather touch disassembly
- Charger port located on the tiller lets you conveniently charge your scooter
- Delta tiller with ergonomic wraparound handles lets you operate the scooter with one hand and rest your wrist
- Easy-to-read LED voltmeter
- Rearview mirror

Colors



True Blue

Candy Apple Red

Optional Seating

Limited recline, black vinyl, high back with height adjustable headrest, flip-up height and width adjustable armrests:

- 20"x 20"

Black vinyl, high back without sliders:

- 22"x18"

Accessories

- Cane/crutch holder (single, double)
- Cell phone holder
- Cup holder
- Forearm crutch holder
- Oxygen tank holder
- Quad cane holder
- Rear basket
- Additional rearview mirror
- Saddlebag
- Safety flag
- Walker holder
- Weather cover
- Wishbone crutch holder
- XLR USB charger
- Lap belts: 50", 60" and 70"

Military Patches

