

C6



A comfortable Infinite Position Lift Chair that's stylish and durable.

The C6's innovative frame design allows all upholstery components to be removed making it easy to match any décor. All steel frame construction gives the C6 unmatched durability for lasting comfort and style.

AUSTRALIA'S #1
POWER LIFT RECLINER

Pride
POWER
LIFT
RECLINERS

Medium to Large

SPECIFICATIONS	C6
	Infinite Position Chaise Lounger with fully removable upholstery components 163 cm - 183 cm (5'4" - 6'0")
Weight Capacity	136 kg (300 lbs.)
Back Width	71 cm (28")
Seat-to-Floor	48 cm (19")
Seat Depth	53 cm (21")
Seat Width	56 cm (22")
Top of Back to Seat	66 cm - 71 cm (26-28")
Back Style	Blown Fiber Sewn Pillow-Back
Fully Padded Chaise	Yes ^{1,2}
Standard Arm Covers	No
Standard Head Covers	No

FEATURES

- Steel frame for added strength and durability
- Infinite-position enables the backrest and footrest to move independently from one another so that a variety of positions can be reached
- Standard large dual pockets
- Easy-to-use hand control with USB charging port
- Removable back
- Energy-efficient transformer, features self-diagnostic electronics
- Quiet and smooth lift system
- Heavy-duty lift actuator and scissor mechanism
- Integrated battery backup in case of power failure³

AVAILABLE FABRICS



Oatmeal



Chocolate
(Micro-Suede)



Arctic Blue
(Micro-Suede)



Royal Blue
(Levante Vinyl)

Note: Measurements may vary +/- 2.5 cm (1").
Specifications subject to change without prior notice.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

1. Removable fully padded chaise.
2. All upholstery components are fully removable.
3. Batteries are not included.

20-24 Apollo Drive
Hallam, Victoria 3803
ACN 088 609 661 ABN 66 088 609 661

Tel: 03 8770 9600
Fax: 03 9703 2720
www.pridemobility.com.au



Pride
Mobility Products
Australia Pty. Ltd.
— LIVE YOUR BEST —