



Celebrity X

3 AND 4-WHEEL SCOOTERS

The all-new Celebrity X Scooter, available in 3 or 4-wheel format, offers elegant simplicity, unmatched reliability and impressive performance. The look is fresh and exciting, coupled with ease-of-use standard features. Plus, the overall scooter design of the Celebrity X provides sleek, next generation styling, while delivering advanced features usually reserved for higher-priced scooters. The Celebrity X is a perfect blend of style and performance technology, making it the best overall luxury scooter value in its class.



Easily disassembles into 7 components



Includes bright, long-lasting LED lights



Easy Drive tiller with wraparound handles



Specifications

Weight capacity	159 kg (350 lbs.)
Maximum speed	Up to 9.5 km/h (6 mph)
Ground clearance	9 cm (3.5")
Turning radius	3-wheel: 104 cm (41") 4-wheel: 145 cm (57")
Overall length	3-wheel: 117 cm (46") 4-wheel: 119 cm (46.75")
Overall width	3-wheel: 61 cm (24") 4-wheel: 61 cm (24")

Tyres

Front	9 cm x 25 cm (3.5" x 10")
Rear	10 cm x 25 cm (4" x 10")
Anti-tips	Rear
Standard seating	Deluxe medium back seat
Suspension	3-wheel: rear 4-wheel: front and rear
Drive train	Sealed transaxle, 24-volt DC motor
Dual braking system	Electronic, regenerative and electromechanical
Battery charger	Off-board, 5 amp
Battery requirements	(2) 32AH (Standard) (2) 40AH (Optional)

Weight

Battery	(2) 11 kg (24.5 lbs.) each
Total w/o batteries	3-wheel: 64 kg (141.5 lbs.) 4-wheel: 79 kg (149.5 lbs.)
Heaviest piece when disassembled	3-wheel: 25 kg (56 lbs.) 4-wheel: 30 kg (65.5 lbs.)
Warranty	2 year limited

Colours



Red

"The information contained herein is correct at the time of publication; due to Pride's commitment to constant improvement and development, we reserve the right to alter specifications without prior notice".

"Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressures, ambient temperature, and the gradient and surface of the road or pavement".

Benefits

- State-of-the-art rear suspension
- Top speed of up to 9.5 km/h (6 mph)
- Maximum weight capacity of 159 kg (350 lbs.)
- Ultra-simple lock-up for quick disassembly
- Easily accessible tie down points (for transport of unoccupied scooter)
- Front basket and rearview mirror standard
- Tight turning radius:
3-wheel: 104 cm (41")
4-wheel: 145 cm (57")

Options

- Forearm crutch holder
- Cane/crutch holder
- Power Elevating seat
- Cup holder
- Oxygen tank holder
- Rear basket
- Safety flag
- Walker holder
- Heavy-duty dust cover
- Saddlebag
- Lap belt

