

# *The Smoothest Ride* in Travel Mobility



4-wheel



**GO GO**®  
*LX* WITH CTS SUSPENSION



Comfort-Trac Suspension

The **Go-Go® LX with Comfort-Trac Suspension (CTS)** offers a whole new level of comfort and performance. This industry-first CTS technology incorporates advanced suspension for improved absorption over various terrains. The Go-Go® LX with CTS offers a 136 kg (300 lb.) weight capacity, a per charge range up to 15 km (9.4 miles) and a maximum speed of 7 km/h (4.4 mph). Standard features include feather-touch disassembly, a wraparound delta tiller with charger port, black, non-scuffing tires, LED headlight and front basket.

## Specifications 4W: S54LX

Weight Capacity	136 kg (300 lbs.)
Maximum Speed <sup>1</sup>	Up to 7 km/h (4.4 mph)
Ground Clearance <sup>2</sup>	3.6 cm (1.5") at motor
Turning Radius <sup>2</sup>	116.2 cm (44.75")
Overall Length <sup>2</sup>	106 cm (41.75")
Overall Width <sup>2</sup>	52.1 cm (20.5")
Front Tyres	5 x 17.8 cm (2" x 7") solid
Rear Tyres	6.35 x 20.3 cm (2.5" x 8") solid
Range Per Charge <sup>1,4</sup>	Up to 15 km (9.4 miles)
Total Weight Without Batteries <sup>5</sup>	41 kg (90 lbs.)
Heaviest Piece When Disassembled	18 kg (39 lbs.) (front section)
Standard Seating	Type: Foldable, compact, vinyl-covered, molded plastic Weight: 10.5 kg (23 lbs.) Material: Black Vinyl Dimensions: Width: 43.2 cm (17") Depth: 43.2 cm (17")
Motor Output	24-volt DC motor, sealed mini-transaxle
Brakes	Regenerative and electromechanical
Battery Requirements <sup>3,5</sup>	(2) 12V Deep-cycle Size: (2) 18 AH (standard) Weight: 14 kg (30 lbs.) per battery pack
Battery Charger	Off-board, 2 amp
Colours	  Interchangeable shrouds (both included) Red Blue

## Features

- Charger port in tiller
- Comfort-Trac front and rear suspension for added comfort
- Convenient dual voltage charger can charge battery pack on-board or off-board
- Feather-touch disassembly permits simple frame separation with only one hand and easily disassembles into 5 lightweight pieces for convenient transport and storage
- Front frame-mounted seat post offers maximum stability
- Long-lasting, high-intensity, front LED headlight offers optimal pathway illumination
- Microprocessor-based controller offers optimal power management and added safety features
- Pride's exclusive black, non-scuffing tyres
- Standard front tiller-mounted basket
- Standard 18 AH battery pack

## Accessories

- Pride<sup>®</sup> rear bag (includes crutch/walking stick holder)
- Walker holder (must have Pride<sup>®</sup> rear bag to attach)
- Oxygen tank holder (must have Pride<sup>®</sup> rear bag to attach)
- Saddlebag
- Mobile phone holder
- Cup holder
- Flag
- Rearview mirror

<sup>1</sup> Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tyre condition. These specifications can be subject to a variance of (+/- 10%).

<sup>2</sup> Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).

<sup>3</sup> AGM or gel cell type required.

<sup>4</sup> Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

<sup>5</sup> Battery weight may vary based on manufacturer.

**NOTE: All specifications subject to change without notice.**

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tyre pressure.

Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressure, ambient temperature, and the gradient and surface of the road or pavement.