

Economy Class Price, First Class Features



GO-GO[®]
ULTRA X



The **Go-Go[®] Ultra X** is loaded with features such as one-hand disassembly and a convenient drop-in battery box for worry-free travel. With a 118 kg (260 lb.) weight capacity, a maximum speed up to 6.4 km/h (4 mph), a per-charge range up to 11 km (6.9 miles) on the 3-wheel and 11.5 km (7.2 miles) on the 4-wheel, Go-Go[®] Ultra X is an exceptional value.

PRIDE[®] MOBILITY SCOOTERS
LIVE YOUR BEST[®]

pridemobility.com.au



Specifications 3W: 40X 4W: 44X

Weight Capacity	118 kg (260 lbs.)
Maximum Speed ¹	Up to 6.4 km/h (4 mph)
Ground Clearance ²	5.7 cm (2.25") at mid-frame 3.2 cm (1.25") at motor
Turning Radius	3W: 83.8 cm (33") 4W: 112 cm (44.125")
Overall Length ²	3W: 94 cm (37") 4W: 100.3 cm (39.5")
Overall Width ²	49.5 cm (19.5")
Front Tyres	3W: 5 x 20.3 cm (2" x 8") solid 4W: 5 x 17.8 cm (2" x 7") solid
Rear Tyres	6.35 x 20.3 cm (2.5" x 8") solid
Range Per Charge ^{1,4}	3W: Up to 11 km (6.9 miles) (12 AH) 4W: Up to 11.5 km (7.2 miles) (12 AH)
Total Weight Without Batteries	3W: 34.5 kg (76 lbs.) 4W: 36.3 kg (80 lbs.)
Heaviest Piece When Disassembled	3W: 13.4 kg (29.5 lbs.) (front section) 4W: 15 kg (33 lbs.) (front section)
Standard Seating	Type: Foldable, compact, vinyl-covered, molded plastic Weight: 10.4 kg (23 lbs.) Material: Black Vinyl Dimensions: Width: 43 cm (17") Depth: 43 cm (17")
Motor Output	24-volt DC motor, sealed mini-transaxle
Brakes	Regenerative and electromechanical
Battery Requirements ^{3,5}	(2) 12V Deep-cycle Size: (2) 12 AH, 4 kg (9 lbs.) each Weight: 12 AH: 8.8 kg (19.5 lbs.) per battery pack
Battery Charger	Off-board, 2 amp
Colour	



Red

Features

- Feather-touch disassembly permits simple frame separation with only one hand and easily disassembles into five super lightweight pieces for convenient transport and storage
- Black, non-scuffing tires exclusively from Pride[®]
- Front frame-mounted seat post offers maximum stability
- Convenient off-board dual voltage charger can charge battery pack on-board or off-board



Feather-touch disassembly

Accessories

- Pride[®] rear bag (includes crutch/walking stick holder)
- Walker holder (must have Pride[®] rear bag to attach)
- Oxygen tank holder (must have Pride[®] rear bag to attach)
- Saddlebag
- Mobile phone holder
- Cup holder
- Flag

¹ Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tyre condition. These specifications can be subject to a variance of (+/- 10%).

² Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).

³ AGM or gel cell type required.

⁴ Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

⁵ Battery weight may vary based on manufacturer.

NOTE: All specifications subject to change without notice.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tyre pressure.

Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressure, ambient temperature, and the gradient and surface of the road or pavement.