

Jazzy[®] **1650**

The Jazzy[®] 1650 offers heavy-duty construction and performance components including upgraded motors, a reinforced frame, enhanced springs on the front anti-tip wheels, and the use of two 70-amp hour batteries. This allows the 1650 to accommodate a maximum weight capacity of 295 kg (650 lbs.) In addition, the Jazzy 1650 delivers style, high performance, and superb manoeuvrability.



Jazzy[®] 1650

Shown in Black Onyx with heavy-duty contoured seat

21 Healey Rd., Dandenong, Victoria 3175

ACN 088 609 661 ABN 66 088 609 661

(Tel) 03 9706 4611 • (Fax) 03 9706 4622 • www.pridemobility.com

Pride[®]

Mobility Products
Australia Pty. Ltd.

— LIVE YOUR BEST —

SPECIFICATIONS

| | |
|-------------------------------------|---|
| Weight capacity | Up to 295 kg (650 lbs.) |
| Maximum speed¹ | Up to 8 km/h (5 mph) |
| Ground clearance | 11 cm (4.25") |
| Turning radius | 67 cm (26.25") |
| Overall length | 119 cm (47") |
| Overall width | 77 cm (30.375") |
| Drive wheels | 41 cm (16") |
| Anti-tips | 23 cm (9") solid |
| Casters | 20 cm (8") solid |
| Suspension | Limited |
| Drive train | Two-motor, 4-pole, in-line, mid-wheel drive |
| Braking system | Intelligent braking (electronic, regenerative disc brakes) |
| Standard electronics | 100 amp, Pilot Plus controller |
| Standard seat | Heavy-duty contoured seat |
| Specialty controls | No |
| Battery charger | Off-board, 8 amp (standard) |
| Per-charge range¹ | Up to 24 km (15 miles) |
| Battery requirements | (2) 12 volt deep-cycle |
| Battery size (largest) | 70 AH |
| Battery weight² | 24 kg (53.5 lbs.) each |
| Base weight | 80 kg (176 lbs.) |
| Standard seat weight | 34 kg (75 lbs.) (heavy-duty contoured seat) |
| Std. seating dimensions | N/A |
| Warranty | 2-year limited |


COLOURS



Candy
Apple Red



Viper Blue



Onyx Black



Champagne

1) Range and speed vary with user weight, terrain type, battery charge, battery condition and tyre pressure.

2) Battery weight may vary based on manufacturer.

21 Healey Rd., Dandenong, Victoria 3175

ACN 088 609 661 ABN 66 088 609 661

(Tel) 03 9706 4611 • (Fax) 03 9706 4622 • www.pridemobility.com

The information contained herein is correct at the time of publication; due to the manufacturer's commitment to constant improvement and development, they reserve the right to alter specifications without prior notice. Theoretical range is calculated under test conditions in accordance with European standards. Actual range in normal use will depend on many factors, including the condition of the vehicle and its batteries, the weight of the driver, correct tyre pressure, ambient temperature, and the gradient and surface of the road or pavement.

FEATURES

- Heavy-duty, in-line dual motors; patented mid-wheel drive design
- Reinforced frame and heavy-duty
- Weight capacity of 295 kg (650 lbs.)
- Accepts larger seating systems for maximum comfort
- 8 amp, off-board battery charger
- 41 cm (16") drive wheels and 23 cm (9") solid caster wheels for maximum stability
- Top speed up to 8 km/h (5 mph)
- Easy front-access dual freewheel levers
- Lifetime limited warranty on frame
- Standard heavy-duty contoured seating for maximum comfort

OPTIONS

- 70° heavy-duty leg rests
- Walker holder
- Lap belt
- Cane/crutch holder
- Rear basket
- Safety flag
- Oxygen tank holder
- Heavy-duty Synergy® Seat
- (4) 55 AH batteries 17 kg (38 lbs. each) for extended range (20 miles)



PRIDE[®]
PERFORMANCE MOBILITY
— LIVE YOUR BEST[®] —